



A World of Healthcare Solutions



Apollo

USER MANUAL


Rev 1 (16/05/2017)

Apollo Saturn Rotate™

Profiling Chair Bed

Product Code - APH034





- 1 HEAD PANEL
- 2 BED RAILS (ON BOTH SIDES)
- 3 HAND CONTROL x 2
- 4 BED RAIL BUTTON
- 5 FOOT PANEL
- 6 WHEELS WITH BRAKE
- 7 HAND RAILS (ON BOTH SIDES)

INDEX

1. INTRODUCTION	Pg.4
2. GENERAL INFORMATION	Pg.5
2.1 Groups of people mentioned in this manual	Pg.5
2.2 Safety warnings	Pg.6
2.2.1 Safety symbols used	Pg.6
2.2.2 Safety warnings for the owner	Pg.6
2.2.3 Safety instructions for the user	Pg.7
2.3 Product description	Pg.8
2.3.1 Proper use	Pg.8
2.3.2 Special features	Pg.8
2.3.3 Materials used	Pg.8
2.3.4 Frame	Pg.9
3. ASSEMBLY AND OPERATIONS	Pg.10
3.1 Assembly	Pg.11
3.2 Operations	Pg.15
3.2.1 Requirements of the installation site	Pg.15
3.2.2 Operations	Pg.15
4. FUNCTIONING	Pg.16
4.1 Special safety warnings	Pg.16
4.2 Hand control	Pg.16
4.2.1 Hand control use	Pg.17
4.3 Swivel wheels	Pg.18
4.4 Lifting Pole	Pg.18
4.5 Bed rails	Pg.19
4.6 Hand rails	Pg.20
4.7 Mattress	Pg.20
5. CLEANING AND DISINFECTING	Pg.21
5.1 General information	Pg.21
5.2 Clearing and disinfection program	Pg.21
5.3 Training of users	Pg.22
5.4 Detergents and disinfectants	Pg.22
5.4.1 How to use disinfectants	Pg.23
6. MAINTENANCE	Pg.23
6.1 By the user	Pg.24
6.2 By the owner	Pg.25
6.3 Spare parts	Pg.28
6.4 Manufacturer's address	Pg.29
6.5 Electrical components replacement	Pg.29
7. ELIMINATING DEFECTS	Pg.30
8. TECHNICAL DATA	Pg.31
8.1 Size and weight	Pg.31
8.2 Electrical data	Pg.31
8.3 Environmental conditions	Pg.31
8.4 Norms and applied guidelines	Pg.32
9. WARRANTY	Pg.32
10. DECLARATION OF CONFORMITY	Pg.33

1. INTRODUCTION

Dear customer,

Apollo Healthcare Technologies wishes to thank you for purchasing one of our nursing beds and for the trust shown in our products.

The electrical safety and functionality of each bed have been tested by the manufacturer in order to avoid any faults.

This manual contains important information regarding the safe setup and use of the bed. It is important that anybody involved should have read and understood its contents. Failure to do so could result in damage to the bed and may invalidate the warranty. Please be aware that elements of the bed are heavy. Those involved should observe the safe lifting and handling techniques recommended under Health & Safety Regulations

This manual provides useful information to the buyer and to the end user on all the functions which will allow a safe and practical daily use of this bed.

Therefore we recommend you to use this manual as a handy guide and to keep it near your bed.

2. GENERAL INFORMATION



Before using the bed:

- Read through this whole manual in order to avoid causing damages to the bed due to improper use.

Before using the nursing bed, the user must verify its perfect status and complete efficiency. This bed conforms to all the requirements of Directive 93/42/EEC for medical products and it is classified as a Class I medical device.

This manual contains safety instructions which must be fully respected. All the people involved with using the bed must know the complete contents of this manual and observe the safety warnings listed herein.

Information for the owner:

- Every electrical appliance may cause hazards if used incorrectly. Therefore it is necessary to instruct the user on how to use the nursing bed correctly.
- Inform the user of where this manual will be kept.

2.1 Groups of people mentioned in this manual

The manual may refer to the following groups of people somehow involved with using the bed:

- **Owner**
The owner (for example: suppliers of medical products, dealers, health care authorities, etc.) can be any physical person or juridical body which uses the nursing bed or under whose name the bed is being used. The owner is responsible for providing the necessary information to the user.
- **User**
the user is the person authorized to use the nursing bed and to carry out any operations on the bed on the basis of his/her training, experience or education, and who has been trained to use the bed. Furthermore, the user must be able to recognize and avoid potential hazards and evaluate the condition of the patient.
- **Patient**
In this manual the term patient refers to the person who needs to be taken care of, who might be sick or disabled and therefore in need of using this bed.
- **Specialized personnel**
The term specialized personnel refers to the owner's co-operators who are authorized, on the basis of their training or education, to deliver and assemble.

2.2 Safety warnings

The nursing bed is a high technology product, therefore make sure to read the complete manual when delivered. Make sure to use the nursing bed only if it is in good condition. The goal of the safety warnings is to avoid injuries to people.

2.2.1 Safety symbols used

The following symbols appear in this manual:

Warnings against injuries



This symbol indicates the presence of dangerous power voltage which could jeopardize the user's life.




This symbol indicates the presence of hazards which might jeopardize the user's life and cause health problems.

Warnings against damages to property



This symbol indicates a potential damage to things or properties. Damages to the bed mechanism, to its materials or to the environment might occur.

Other warnings



This symbol indicates a general useful suggestion. Following this suggestion, it will be easier to use the bed and to understand the instructions.

The safety symbol used does not replace the safety text it refers to. It is always necessary to read the relevant safety warning in full and to follow it carefully!

2.2.2 Safety warnings for the owner

- Before the bed is used, inform the user about safely operating the bed by strictly following this manual which must be provided together with the bed. Inform each user of any possible dangers due to any improper use. Pay special attention when using the electric functions and the bed rails.
- Nursing beds are classified as active Class I medical devices. You must fulfil your obligations in accordance with the current regulations about installation, operation and use of medical devices, to ensure a safe and durable use of this medical product without any risks for patients, users and any third parties. In case of a prolonged use of the bed, it is necessary to carry out an inspection, after a certain period of time (we recommend at least every year), to check its functionality and any visible damages (see Chapter 6).
- Make sure that only qualified personnel use this bed.

- Make sure that any possible substitute personnel are also sufficiently trained on how to use the bed correctly.
- In case any additional electric accessories are used such as air compressors for overlay mattresses, make sure they are properly installed and work correctly. In particular pay attention to the following:
 - Make sure all cables, flexible tubes etc are fixed correctly.
 - Do not use a multiple socket under the bed (possible fire hazard in case of contact with liquids). If in doubt, please contact the manufacturer of the accessories.
- Make sure all staff complies with the safety warnings.

2.2.3 Safety instructions for the user

- Make sure to be instructed on the safe use of the bed by the owner.
- Before each use check the good condition and the absence of defects in the nursing bed.
- Make sure there are no obstructions such as furniture or inclinations of the ceiling when making any adjustments.
- In case any additional electric components are also used such as hoists, reading lamps or air compressors for overlay mattresses, make sure their cable does not get caught in the mobile parts of the bed or get damaged.
- Multiple sockets must not be used for connecting any external electrical components.
- If you suspect damage or malfunctioning, immediately disconnect the power plug from the wall, clearly mark it as "Damaged" and deactivate the nursing bed. Immediately inform the owner of the bed about this situation. In chapter 6.1 you will find a checklist to assess the condition of the bed.

Safety warnings for the user



- The power cord must be installed so that when operating the bed it is not carried away, crushed or damaged by any moving parts.
 - Before transport totally remove the plug from the socket.
 - When the bed is not used, make sure there is no possibility for the hand control to fall and get damaged (it must remain attached to the hook). Make sure the cable cannot be damaged by any moving parts of the nursing bed.
 - Place the bed in the lowest position when you leave the patient unattended on the bed. In so doing, the risk of injuries to the patient when getting on and off the bed is minimized.
 - The adjustments can be made only by trained personnel or in the presence of experienced staff.
 - Check the power cord regularly for possible mechanical damages (abrasions, uncovered wires, bends, dents, etc...) especially in the following cases:
 - after each mechanical overload (for example: crushing of the cable with the bed itself or with a trolley, after strong tractions and bending of the bed with the power cord plugged in)
 - after each change of position/after every movement before inserting the plug
 - during the normal operation by the user
- Regularly check that the power cord assembly is properly tightened.

2.3 Product description

2.3.1 Proper use

- The nursing bed, hereinafter simply called “bed”, is designed for disabled or ill people and is suitable also for home care.
- The bed was not specifically designed for hospital use.
- This bed can be used for treatment with a doctor's prescription and it can be used for diagnosis, care or observation of the patient.
- This bed does not have any special options suitable for an equipotential connection. Make sure to take this aspect into account before connecting with any other (medical) equipment already linked to the power supply.
- This bed must be used only by adequately trained people.
- This bed is designed for repeated re-use. Please follow the instructions listed below:
 - Cleaning and disinfecting (see chapter 5)
 - Maintenance / Verification of tests (see chapter 6.2)

This bed can be used only according to the use conditions described in this manual. Any other use is not to be considered as proper.

2.3.2 Features

- Electric height adjustment of the bed surface from 27 to 68 cm.
- Electric backrest (with Auto-Regression) section adjustment from 0° to 72°.
- Electric leg section adjustment from 0° to 33°.
- Electric leg section adjustment from 0° to 20°.
- Four swivel wheels, individually lockable.
- Hospital surface 210 x 88 cm.
- Overall dimensions about 232.4 x 105.2cm.
- Drop down 3 bar bed rails on both sides.

The Saturn bed has a number of electrically operated profiling functions – back rest, leg rest and bed height are adjustable together with Fowler position knee-break.

Similarly powered functions raise the user from a lying position and rotate to a chair position to assist transfer from the bed and re-entry from a standing position.

Removable hand-rails provide additional support during transfer.

The functions are powered by linear actuators which are sealed and rely on maintenance-free permanent lubrication. They are ingress protected to IP54 standard.

2.3.3 Materials used

The bed is mostly made up by steel tubing profiles, covered by a polyester powder coating or by a metallic zinc or chrome coating. The head/foot panels are made of wood or similar material, whose surface is waterproof. The bed rails are made of aluminium, coated to mimic the appearance of wood.

The surfaces are not harmful in contact with the skin.

2.3.4 Frame

- The bed is supplied disassembled to ease transportation. It is made up of two wooden panels (head panel and foot panel), a mattress base with a motors, a pantograph mechanism with an electric actuator, six bed rail bars and 2 bed arms. The bed has four swivel wheels, all with a parking brake.

- **Mattress base**

The mattress base is supplied complete with backrest with auto-regression, a fixed central element, a leg support and knee raiser and fixed foot section. The backrest and leg support sections can be adjusted to profile the mattress platform. Once fully profiled, the mattress platform can be rotated around through 90° and tilted to aid ingress and egress from the bed. The whole bed can be adjusted in height horizontally or obliquely compared to the floor. All the adjustments are performed by electric motors by using a hand control.

- **Bed Rails**

To protect patients from possible accidental falls from the bed, the bed is equipped with rails (meet clause 201.9.8.3.3.3 of IEC 60601-2-52) on both sides, which can act as a barrier or be lowered when not in use. The blocking mechanism is locked and it can be released only by pressing the button on the side.

- **Electric system**

The electric operating system of this bed is safe and easy to use. It consists of:

- Electronic controllers. A minimum protection voltage of 24 volts is produced which is not dangerous for the patient and the user. The electric motors and the hand controls are connected to the central controller with connectors, which work with a minimum protection voltage of 24 volts.
- Separate electric motors for the back and the leg supports and to rotate and tilt the bed.
- Two electric motors to adjust the height of the mattress base (one in each bed end).
- 2 hand control units with fixed hook.

3. ASSEMBLY AND OPERATION

- The bed is supplied disassembled.
- The bed should be assembled on location by suitably competent persons.
- The installation must be carried out by a minimum of two persons.
- Ensure that the local electricity mains voltage corresponds to that marked on the main controller label before connecting to the supply.
- Ensure that cables from actuators are plugged into the main controller correctly.
- The fuse in the mains 'safety connector' plug should not exceed 5amps.
- The bed should be located on a level surface & not sited on loose floor coverings.
- The cable from the mains electricity supply must be kept away from the lifting mechanism & casters to avoid danger of shearing or crushing.
- The brake on each caster must be applied whilst the bed is in use.
- Do not position objects under the bed.
- Always ensure that limbs or body parts of both user and carer are not protruding from the bed, are not between the side-rails or in the way of the lifting mechanism when using the powered functions.
- When raising/lowering the bed, be aware of any fixtures such as window sills, shelf units and electrical power sockets which may obstruct the raising of the bed or cause damage to the fixtures. Re-position the bed if necessary.
- If moving the bed to a different location, always ensure that all bed functions are at their lowest position before disassembling and remember to remove the plugs from the mains electricity power supply.
- Never use the side-rails to push or pull the bed. Always use the head and foot boards.
- With or without occupant, the bed should only be moved at slow speed and not pushed over a threshold.
- Once the bed is re-sited, remember to re-apply the brakes.

Packaging units:

Carton 1: Bed frame / mattress platform complete with back rest, knee break, rotating and tilt function actuators, 2 control boxes and handsets, side rail end caps & mattress retainers.

Carton 2: Headboard Bed end with a Hi-Lo actuator.

Carton 3: Footboard Bed end with a Hi-Lo actuator.

Carton 4: Six wood effect aluminium bars to form the bed side rails.

Carton 5: Two hand rails.

3.1 Assembly

Mounting instructions:

- The assembly of the bed must be carried out by a minimum of two persons;

Open cartons 2/5 and 3/5 both labelled Saturn Bed End and remove both beds ends and swing out the two arms 90° towards the inside of the bed. In the image this image the arm on the right has not been swung out but the arm on the left has been swung out to the correct position. Once both arms on both bed ends have been swung out they are ready to be fitted to the bed frame.



Open carton 1/5 labelled APH034 Saturn Bed Frame and remove the entire bed frame from the packaging and unscrew the four hand screws at each corner as shown below in the image on the left.



You can then insert the arms of the bed ends into the either end of the frame. It does not matter which bed end is inserted into each end of the frame. When all bed end arms have been fully inserted into the bed frame, as shown in the right hand image above, you can hand tighten the hand screws until they lock the arms in place.

The tilt function actuator is only connected the bed frame at one end so the other end requires connecting. This is done by removing the spring loaded pin, swinging the actuator up to the bracket (indicated by the red circle in the image to the right) and inserting the pin before locking it in place with the spring loaded loop. The electrical cable is already connected to the control box.




The rotating function actuator is also only attached to the bed frame at one end so the other end needs connecting to the bed frame. This is again done by removing the pin and swinging the actuator round to the bracket (shown in the red oval in the image to the right) and then inserting the pin before locking it in place with the loop. The electrical cable connecting it to the control box is already connected.



Depending on which direction you want the chair/bed to rotate and then rise to aid the transfer between sitting and standing, the rotate and tilt actuators can be moved to the opposite side of the frame which will make the chair/bed rotate the other way. The chair/bed will rotate through 90° toward the opposite side of the frame where the actuators are mounted.


To move the rotate actuator remove the locking pins at either end of the actuator and move the whole actuator to the opposite side of the bed frame.



The image on the left shows the rotating actuator in the standard position on the left of the bed frame connected at both ends (indicated by the red circles)

To move it to the opposite side of the bed frame you need to remove both pins and move the actuator so that it is attached to the bed frame on the brackets indicated in the green circles

Make sure that the actuator remains the same way around, do not rotate the actuator when moving it.




To move the tilt actuator to the opposite side of the bed you will need to remove the actuator and the tilting mechanism as a whole.

The image on the left shows the tilt actuator and tilting mechanism in the standard position (on the same side of the frame as the rotate actuator). You need to remove the locking pins indicated by the red circles and also unscrew the bolt indicated by the blue circle which will enable you to remove the actuator and mechanism together.

The image on the right shows the mounting positions on the bed frame for the tilt actuator and tilting mechanism to be installed on the opposite side.

The green circles indicate where the pins should be inserted to fix the actuator to the frame and the yellow circle indicates where the bolt should be used to secure the tilting mechanism to the bed frame.



Please take note of the way the bolt and washers were installed before removing them so that they can be re-installed in the same manner.

Do not rotate the actuator and mechanism when moving to the opposite side of the bed. Ensure that the actuator is towards the foot of the bed as it were before moving from the other side of the frame.

As well as moving the actuator to the opposite side of the bed frame you will also need to move the tilt locking pins to the opposite side. There is a tilt locking pin at both the head and foot of the bed frame.



The image on the left shows the tilt locking pin at the head of the bed in its standard position, shown by the red circle. Note that in the standard position it is on the opposite side of the bed frame to the actuator.

The tilt locking pin must be removed from the standard position (red circle) and moved to the opposite side (green circle).

Repeat this process with the tilt locking pin at the foot end of the bed frame.



The tilt locking pins must always be on the opposite side of the frame to the rotate and tilt actuators.



Ensure all pins used to mount actuators are fully inserted and locked into position with the spring loaded loop.

Both the head and foot boards come already equipped with a Hi-Lo actuator that can electrically raise and lower each end of the bed, either together at the same time or opposing each other to engage the Trendelenburg tilt function.

The two main control boxes come already attached to the bed frame with all the cables attached both to themselves and to the correct actuator apart from Hi-Lo actuator cables. You will need to carefully cut the ties (indicated in the image to the right by the red circle) that have been used to secure the cables to the bed frame for transit purposes in order to connect them to the Hi-Lo actuators at both the head and foot ends.




Once the cables are free from the ties they can be inserted into the actuator in the bed end. To do this you will first need to remove the securing clip as shown in the image below left.

You can then push the connector at the end of the cable into the socket in the actuator. Once the connector is fully inserted you can replace the securing clip fully (as shown in the image below right) to prevent the connector from accidentally coming out of the actuator.



Make sure not to damage any cables. Verify that the cables do not get tangled and that they do not get caught in any in moving parts.



In case of replacement of any individual components, always make sure that the electric plugs are then properly inserted into the actuator and controller. Only in this way the bed will function properly.

3.2 Operations

Every time before using the bed, the user must make sure that:

- the wheels of the bed are locked (see also chapter 4.3).
- the power supply is compatible with the bed (230 Volt, 50 Hertz).
- the power cable is plugged in and placed in such a way that it cannot be damaged.
- the power cable, the motor cables and the hand control cable are not damaged by any moving parts of the nursing bed.
- there are no obstacles in the way such as bedsides or chairs.
- all control functions have been accurately inspected (see also chapter 6.1) and are working properly.
- only after having carried out the above operations it will be possible to start using the bed.

3.2.1 Requirements of the installation site

- Make sure there is enough space to allow the whole range of adjustments of the bed. There should not be any furniture or other objects close to the bed such as window cills.
- If the bed is to be used on parquet floor, make sure that the use of wheels will not cause any discolouring due to the parquet waterproof treatment. It is possible to use the bed on tiles, carpets, laminated floors and linoleum.
- A 230 volt power outlet, correctly installed, must be located as near as possible to the bed.
- If there is the need to use other accessories (such as compressor systems etc...), make sure they are correctly fixed and that they work properly. In particular, make sure all different cables are safely connected. In case of questions or doubts please contact the manufacturer of such accessories.



During the installation of the bed you must take the following issues into serious consideration to minimize potential fire hazards due to external influences. Inform the user and the personnel about these issues!

- If possible, use only non-flammable mattresses and linens.
- Don't smoke while in bed as there is no guarantee against risks caused by smoking tools (such as lighters, matches) depending on the mattresses and linens used.
- Use any electric accessories (electric blanket) and appliances (lights, radios) only if they are in a good condition.
- Make sure these appliances are used only for their intended purpose and that they cannot accidentally end up on and under the sheets (risk of overheating).
- Absolutely avoid using extension cables or connections to multiple sockets under the bed (fire hazard due to the penetration of liquids).
- Extension cables and/or power multiple sockets should never be used at all.

3.2.2 Operations

- After the assembly of the bed, make sure to carry out an inspection in conformity with chapter 6.2.

After carrying out all the instructions listed in chapters 3.1 - 3.2.2 the bed can start functioning.

4. FUNCTIONING


4.1 Special safety warnings



- When making the adjustments make sure that no parts of the body of any patients, users and other people, especially children, are positioned under any adjustable elements or under the bed frame to avoid being caught and injured.
- Once you have made the necessary adjustments and positioned the bed rails correctly, make sure to lock the backrest and leg adjustments possibilities with the hand control, to protect the patient from unwanted movements (risk of crushing of his/her limbs).
- In case of use of accessories on the electric bed, make sure such accessories do not crush or cut the patient with sharp edges during the legs and backrest adjustments. If this condition is not guaranteed, the user must not proceed with the backrest and legs adjustments.
- Make sure the power cable and the hand control cable do not get caught anywhere or get damaged.



- Before each operation with the bed, make sure the power cable is not stretched, crushed or damaged in any way.
- When making the adjustments, make sure there are no obstructions such as furniture or a ceiling inclination, to avoid possible damages.



As with all electrical appliances, despite the compliance with prescribed values, it is not possible to totally exclude the possibility of interferences from other electrical devices which might be nearby (for eg. the “buzz” of a radio). In such a rare case, increase the distance of the bed from the interfering device, do not use the same socket or temporarily disconnect such device.

4.2 Hand control

The electric functions of the bed are activated by the patient and by the user with the use of two separate hand controls.

The first hand control operates only the Hi-Lo actuators located in each end of the bed in order to raise and lower the whole bed and to activate the Trendelenburg and anti-Trendelenburg tilt functions.

The second hand control operates the profiling back and foot rests and the rotating and tilting functions when in the chair position.

- The electric motors work if the corresponding buttons are pressed.
- It is possible to adjust different sections in both directions.
- The hand control can be attached to the bed where needed with a hook.
- The spiral cable allows freedom of movement.

4.2.1 Hand control use



Height adjustment of the mattress base

Press the corresponding buttons to raise or lower the height of the mattress base. There is an automatic horizontal positioning in the highest or lowest position.

Angle adjustment of head / foot section (Trendelenburg Tilt)

Press the corresponding buttons to lower the foot or head section independently. Before lowering the section you need to release the brake of the two wheels on the head or foot panel side, to avoid possible damages to the floor.



Backrest adjustment

Press the corresponding buttons to change the angle of the backrest.

Leg support adjustment

Press the corresponding buttons to change the angle of the knee break/leg support.

Chair Rotate function

Press the corresponding buttons to rotate the chair/bed round once fully profiled.

Chair Raise function

Press the corresponding buttons to raise and lower the chair/bed once fully rotated.



The chair/bed will only rotate when the bed has been fully profiled into the chair position. Pressing the rotate button before the bed has been fully profiled will automatically complete the profiling process before rotating.

The chair/bed will only raise once it has fully rotated round the full 90° to the fame. Pressing the raise function on the handset will not work until the chair/bed is fully rotated.

4.3 Swivel wheels

The nursing bed has four swivelling double castors which can be blocked individually.



- Bring the mattress base into its lowest position.
- Before each operation of the bed, make sure the power cable is not stretched, crushed or damaged in any way.

How to apply the brake: Press the rocker lever down with your foot until it clicks.



How to release the brake: Press the rocker lever up with your foot until it clicks.

4.4 Lifting Pole

In both corners of the mattress base, towards the head panel, there is a round insert with a cavity into which the lifting pole must be placed. The lifting pole must be inserted in the side from where the patient gets in or out of the bed to help him/her out. Only use the manufacturers lifting pole design for this bed and for specific Bariatric use.

Assembly

- Place the lifting pole into its housing. The metal pin must be placed into the cavity. Make sure to double check the oscillation field of the lifting pole which must be limited.

Disassembly

- Remove the lifting pole from its housing by pulling it upwards.

Height adjustment of the handle

The height of the lifting pole handle can be adjusted with the strap supplied.

4.5 Bed rails

Fitting

- The easiest way is to raise the bed to its highest position.
- Insert each bed rail into the end caps provided and insert the end cap guides into the runners on the head and foot board.
- Make sure to lift one rail at a time ensuring both ends are level before moving on to the next rail.
- You may have to press in the bottom button in order to fit the bottom bed rail.
- Once all rails have been fitted and are above the bottom button lift up the top bed rail until it automatically fits into the highest position on both ends above the top button. You should not be able to move the top rail upwards or downwards.
- Check that the other bed rails are properly inserted by pressing from the top of each bed rail. You should not be able to move them downwards.



Lowering

- Lift the top bed rail slightly.
- Press the upper button and lower bed rail slowly to its lower position.
- The rails must be in the lowered position before rotating the mattress platform.



Special safety rules for the bed rails.

The bed rails are there to protect patients from possible accidental falls from the bed. However, they cannot prevent the patient from voluntarily leaving the bed. To make the best use of their intended protection function, please observe the following precautions:

- Use only undamaged bed rails and in good conditions positioned at the right distance to allow the rails to snap into place properly.
- Make sure to use only the Apollo original bed rails supplied with the bed.
- When using bed rails you must always verify if they are suitable for each individual patient. In particular, it is necessary to pay attention to the distance between the bars in relation to the patients' build. If they are very skinny and thin it could be necessary to use some padding or pads around the bed rails (optional). Only in this way will it be safe to use the bed without the risk of the patient to be caught or to slide out of the bed rails.
- The mechanical elements which attach the bed and the rails should be checked regularly to make sure there are no signs of damage. This must be done not only before applying the bed rails, but also during use within short intervals (at least before any new use).
- When applying the bed rails to the electric bed always make sure to block the backrest and leg supports with the hand control to protect the patient from any accidental electric movements (risk of involuntarily crushing the patient's limbs).

Failure to comply with the above safety rules will lead to the possible risks of the patient getting caught in the bed rails, falling out of the bed if the bed rails are too wide due to possible damages, in case they are not the correct bed rails for the electric bed or if they have not been mounted and inserted correctly.

4.6 Hand rails

The Saturn bed comes complete with two hand rails that must be fitted on either side of the mattress platform before using the rotating feature of the bed. These are provided to further aid the bed user when transferring from a lying position to a sitting position and to a standing position and back again.



- **The hand rails must be fitted before using the rotating feature of the bed.**
- **The bed rails must be in the lowered position on both sides of the bed in order to fit the hand rails. This is a safety feature to ensure that the side rails are lowered before rotating the bed in order to prevent injury to the user and possible damage to the bed.**

Fitting

On either side of the mattress platform, on the fixed section in the middle there are two holes into which the hand rails can be inserted.

The image on the left shows the holes on the left-hand side (indicated by the arrows) and the hand rail fitted on the right-hand side.

Ensure that the longer part of the hand rail is pointing towards the foot of the bed.



Once the hand rails have been inserted into the holes they can be secured in place by using the small hand screws provided. Only one hand screw is required on each hand rail.

4.7 Mattress

The Saturn bed is supplied with either a static foam mattress or a dynamic air alternating mattress. Both mattresses are covered in a polyurethane coated, high frequency welded cover with straps on the underside to secure it to the mattress platform.

The small separate foot section of the mattress should be secured to the foot of the bed onto the mattress platform that is fixed and does not move when the mattress platform is profiled.

The larger section of the mattress should be positioned on the bed with the Apollo logo at the head end and secured to the mattress platform using the straps provided.



Ensure the straps are only secured to the moving parts of the mattress platform and not secured around any other parts of the frame. If the straps are secured to any non-moving parts of the frame, when the mattress platform is moved the straps will restrict the movement and could severely damage the bed.



Only use the Apollo Healthcare Technologies supplied bed rails, hand rails and mattress. Never under any circumstances should any third party accessories be used and your warranty will be void if you do so.

5. CLEANING AND DISINFECTING

5.1 General information

Clearing is fundamental as well as a proper chemical disinfection. In general, regular cleaning of the bed when used by a single patient is sufficient in terms of hygiene. It is necessary to disinfect the bed frame only in case of a visible contamination with contagious or potentially contagious materials (ex. blood, urine, pus etc.) or in case of a contagious disease upon request by the physician.

Every time the bed is used by a new patient, it is mandatory to first clean and then disinfect the bed by scrubbing thoroughly!



Before you start clearing:

- **Remove the plug from the socket and make sure it does not come into contact with water or disinfectant.**
- **Make sure all plugs are properly connected.**
- **None of the electric components must be externally damaged to avoid the penetration of water or disinfectant liquids which might lead to malfunctioning or damages to the electric components.**
- **Before using the bed again, by drying or blowing in the power plug, make sure there is no residual moisture left in the electrical contacts.**
- **The electric components must not be subjected to a water jet, a high pressure cleaner or anything similar! Cleaning should be done with a damp cloth!**
- **If you suspect that water or moisture have entered into the electric components, immediately remove the power plug from the socket outlet and do not plug it in again. Mark clearly the bed with a "Faulty" sign and make sure it cannot be used. Inform immediately the owner of the bed about the problem.**

In case of non-compliance with these safety rules, it is not possible to exclude the risk of serious damages to the bed with related consequences!

5.2 Cleaning and disinfection program

- Remove the sheets from the bed to be washed.
- Clean all surfaces including the mattress platform, the plastic inserts and cables with a mild non aggressive detergent, including the hand control.
- The bed should be disinfected by scrubbing using a disinfectant recommended for hospital use especially in case of visible contamination with contagious or potentially infectious material. The same as applicable to beds in intensive care or infectious diseases units.
- It is necessary to disinfect the wheels only in case of visible contamination with infectious or potentially infectious material.

Note: A continuous disinfection is necessary only in case of patients with multi-resistant pathogens within a hospital

5.3 Training of users and specialized personnel

To ensure a proper cleaning and disinfecting procedure, we recommend you to appropriately train all users and specialized personnel involved. The following points must be respected:

- Make sure the bed does not get soiled or contaminated during transport to the patient's house.
- In case of disassembly of the bed, it is necessary to immediately clean and disinfect it by scrubbing. The personnel must be informed about the correct cleaning and disinfecting procedure and carry out his/her necessary treatment properly (according to the steps of the procedure notified by the owner). The disinfectant used must be suitable for the treated surface.
- For this operation the specialized personnel shall be provided with waterproof disposable gloves and aprons.
- For the treatment only clean cloths should be used which must then be taken to be washed.
- After the treatment the personnel should disinfect their hands before leaving the place to go somewhere else. The qualified personnel must be equipped with a disinfectant with dispenser.

Cleaning of the bed directly on site has the advantage that the bed and/or none of its parts can get "soiled" when then transported in the same vehicle with other clean beds. It prevents potentially infectious germs located in the bed frame from spreading around. These recommendations will help to prevent nosocomial infections. If the beds are not immediately reused, they should be stored away with a cover to protect from dust and accidental contamination.

5.4 Detergents and disinfectants

Make sure to follow the recommendations listed below to maintain the bed functionality as long as possible:

- All surfaces must be intact and any damages must be immediately eliminated.
- We recommend using a method of disinfection (wet) by scrubbing. The detergent should not be aggressive (not harmful for the skin and for surfaces) and it must be ecological. In general you can use a household detergent.
- To clean and disinfect the painted metal parts use a damp cloth with a common non aggressive mild household detergent.
- The disinfectant for scrubbing the bed must be chosen from among those recommended (in the indicated concentration).
- Despite the mechanical strength of the bed it is necessary to remove any scratches and dents in the paint with special repair kits in order to prevent moisture penetration.



Never use abrasive clearing products and stainless steel cleaners containing solvents or clearing pads, to avoid damaging the bed surface.

Note – Aldehyde-based disinfectants normally cover a wide action range and are ecological. However, they are potentially allergenic and irritating. Glucoprotamina preparations do not present this problem and are very effective even though in most cases

they are quite expensive. Disinfectants based on chlorine splitting compounds can have a corrosive effect on metals, plastics, rubber or other materials in case of a prolonged contact or in high concentrations. In addition, these products can be irritating for mucous membranes and do not have a great environmental resistance.

5.4.1 How to use disinfectants

- Pay attention to the proper dosage! We recommend using automatic dispensers.
- Prepare the solution always with cold water to prevent the build-up of vapours which might be irritating for mucous membranes.
- Wear gloves to avoid any direct contact with the skin.
- Do not store the disinfectant solutions in open bottles or containers. Make sure to always close them!
- Use bottles with pump dispensers to wet the cleaning cloth.
- Ventilate the room after the disinfection.
- Carry out the disinfection by scrubbing and not by spraying the disinfectant in a large quantitative with the risk of being inhaled.
- This is why rubbing plays an important role.
- Do not use alcohol for large surfaces.

6. MAINTENANCE

Regulations

The owners of the nursing beds, in accordance with the applicable regulations in force, are requested to guarantee the safety conditions of the medical products during their use. To achieve this, it is therefore necessary to regularly carry out inspections and maintenance interventions.

Notes for the owner

This bed was manufactured to function properly for a long period of time. When used correctly, it can be operated for around 2 to 8 years depending on the conditions and frequency of its use.



After repeated operations of transport, assembly, disassembly, improper and long-term use it is not possible to exclude that any damages, faults, wear and tear problems might occur. These problems may cause safety risks if not recognized and solved in time!

Therefore it is necessary both for the owner and the user to follow the safety norms for medical devices with the implementation of regular inspections to verify the correct functioning of the bed.

Inform the user that it is mandatory to carry out these controls.

6.1 By the user

In addition to regular checks by the specialized staff, even the normal users (such as nurses, assistants, etc.) must perform both a visual and a functional inspection before each new use at regular intervals.

Recommendation: check all mechanical and electrical components once a month. In addition check the power cord and the hand control cable after each mechanical stress and after each change of position.

See the checklist below.

Checklist for the user

Control		ok	not ok	Problem description
Visual inspection of electrical components				
Hand control cables	Possible damages, cable connection			
Power cords	Possible damages, cable connection			
Hand controls	Possible damages			
Visual inspection of mechanical components				
Lifting pole and its housings	Possible damages, deformations			
Bed frame	Possible damages, deformations			
Wooden slats	Possible damages, wooden chips			
Wooden head and foot panels	Possible damages, wooden chips			
Mattress base	Possible damages, deformations			
Wood effect bed rail bars	Possible damages, wooden chips			
Functional control of electrical components				
Hand controls	Test to check if it works properly			
Functional control of mechanical components				
Casters	Test the brakes, moving forwards			
Emergency descent of the backrest	Test according to the instructions			
Bed rails	Locking and unlocking mechanisms			
Knee raiser	Mechanism and attachment			

Accessories	Fixing, possible damages			
Inspector's signature:	Control result:			Date:



- If you think that there might be damage or a functional problem, the bed must not be used at all and it must be disconnected from the power supply until the damaged components have been totally replaced or repaired!
- To replace or repair damaged parts, please contact the owner of the bed.

6.2 By the owner

The owner of this nursing bed is required to carry out regular tests in case of any new installation, any maintenance service and during the normal use operations to ensure the safety of this bed.

These tests must be verified during the regular maintenance services depending on use conditions.

- The functional tests, the evaluations, the documentation of the inspection results must be carried out only by experienced personnel who know the proper way and conditions necessary to implement a correct testing.
- The evaluation and the documentation of the results can be performed only by skilled electricians with additional medical and technical knowledge.
- The controls must be carried out according to the sequence below:
 1. Visual inspections
 2. Electrical measurements
 3. Functional tests

Frequency of the inspections

As a benchmark, we recommend an examination every six months can be extended in accordance with the terms of use on its own responsibility respecting a share of proven error of 2% to a maximum of two years.



- If you think that there might be a damage or a functional problem, the bed must not be used at all and it must be immediately disconnected from the power supply until the damaged components have been totally replaced or repaired!

The pages below can be used as an inspection protocol.

Sheet 1 of 2

Client/ Medical Centre/ Surgery:				
Address:				
Task carried out: <input type="checkbox"/> Testing of the bed <input type="checkbox"/> Testing before the first use of the bed				
<input type="checkbox"/> _____ <input type="checkbox"/> Testing after repair / maintenance				
Device type: <input type="checkbox"/> Nursing Bed		Protection class: <input type="checkbox"/> I		
Bed type:		Inventory number:		
Place:		Serial number:		
Manufacturer: Apollo Healthcare Technologies Ltd.		Parts used: none		
Equipment used (Type/ number of inventory): 1.				
Classification according MPG: class I 2.				
I. Visual inspection		ok	not ok	Problem description
Visual inspection of electrical components				
Identification labels and stickers	In place and clearly readable			
Housing of the control device	Properly positioned, damages			
Housing and lifting tubes of the motors	Properly positioned, damages			
Hand control	Damages			
Power cable of the motors, hand control cable	Properly positioned, damages			
Connector and clamp of the controller	Properly positioned			
Visual inspection of mechanical components				
Identification labels and stickers	In place and clearly readable			
Lifting Pole	Damages, deformations			
Bed frame	Damages, deformations			
Wooden slats	Damages, wooden chips			
Wheels	Damages			
Mattress base	Damages, deformations			
Wooden panels	Damages, wooden chips			
Welding	Welding problems			
Bed rail bars	Damages, wooden chips			

Knobs	Correct fixing			
Worn out parts, such as the hinged points	Damages			
II. Electrical measurement in accordance with VDE 0751-1:2001-10 (power from the device, direct)				
Proceed as follows:		ok	not ok	Problem description
1. Insert the power cord of the bed into the testing outlet of the measuring device.				
2. Connect the probe of the measuring device to a free conductor element of the frame (such as a screw or similar).				
3. Activate the motors with the hand control for the duration of the measurements				
4. Start the measurement procedure with the appropriate device.				
	Limit Value	Value of the first measurement	Current actual value	
Result: bed SK II (type B)	0,1 ma	ma	ma	

Protocol system for testing electrical medical devices

Sheet 2 to 2

III. Functional control		ok	not ok	Problem description
Functional control of electrical components				
End stroke of the motors	Automatic disconnection			
Controller and motors	Unusual noises			
Controller and motors	Test according to use instructions			
Strain relief cord	Connecting cable properly fixed			
Functional control of the mechanical components				
Hinged and pivot points	Working smoothly			
casters	Brakes, moving forwards			
Bed rails	Locking and unlocking mechanisms			
Knee raiser	Hinge point			
Accessories (ex. Trapeze bar, handle)	Mounting, damages, suitability			
Inspection result:				
All values within the allowed range: <input type="checkbox"/> <input type="checkbox"/> yes <input type="checkbox"/> no			Next inspection date:	

The test was passed: <input type="checkbox"/> yes <input type="checkbox"/> no		
If the test was not passed: <input type="checkbox"/> <ul style="list-style-type: none"> <input type="checkbox"/> defective, do not use the bed! => Maintenance <input type="checkbox"/> defective, do not use the bed!! => Elimination <input type="checkbox"/> The bed does not meet the safety requirements 		
Mark test applied: <input type="checkbox"/> yes <input type="checkbox"/> no		
This inspection protocol included:		
Notes:		
Tested on:	Examiner:	Signature:
Evaluated on:	Owner/staff:	Signature:

6.3 Spare parts

Replacement parts are available From Apollo Healthcare Technologies Ltd. by indicating the item number, order number and serial number. The necessary information can be obtained from the identification label placed on the bed.

Model	Bed model
Type	Article number
Series	Serial number

6.4 Manufacturer's address

To obtain a proper functioning and warranty replacements it is necessary to use only original spare parts Apollo Healthcare Technologies Ltd!

To order spare parts, for any requests or assistance and for further information please contact:

Apollo Healthcare Technologies Ltd
 Holme Street
 Liversedge
 West Yorkshire
 WF15 6JF
Tel: 01924 614567 Fax: 01924 607480
Email: info@apollo-ht.co.uk

6.5 Electrical components replacement



Dangerous for your life, life threat, due to electric shock!

- Before starting any operations remove the plug from the wall outlet!
- Any electrical equipment operations must be carried out only by qualified service and repair staff members, by the producer or by qualified and licensed electricians, in accordance with VDE and relevant safety standards!
- Disassemble the operating system only in the horizontal position to avoid the danger of being crushed by any elements falling from the mattress base. The drive components are maintenance free and cannot be opened. In case of malfunction, the component in question must be completely replaced!



7. ELIMINATING DEFECTS

Should any anomalies occur during use, it is necessary to contact the owner's qualified personnel trained for the maintenance and repair of the bed.

The user shall in no case try to eliminate any malfunctions that might occur in the electrical equipment!



Dangerous for your life, life threat, due to electric shock!

Any interventions to the electrical equipment must be carried out only by qualified service and repair staff members, by the producer or by qualified and licensed electricians, in accordance with VDE and relevant safety standards!

Possible solutions

The following table is useful to eliminate possible problems and faults:

Problem	Possible causes	Solution
The hand control and / or the operating system does not work	<ul style="list-style-type: none"> Power cord not connected Plug not under voltage Hand control plug not properly connected Faulty parts 	<ul style="list-style-type: none"> Insert the power cord Check the socket / the fuse box Check all connections Inform the owner and ask for repairs
Pressing the keys, the movements are activated only for a short time	<ul style="list-style-type: none"> Too much weight on the bed The bed hits against an obstacle 	<ul style="list-style-type: none"> Reduce the load Remove the obstacle
Movement in one direction only	<ul style="list-style-type: none"> Defective hand control or operating console 	<ul style="list-style-type: none"> Inform the owner and ask for repairs
Function in contrast with the instruction on the hand control	<ul style="list-style-type: none"> Reversed motor connectors 	<ul style="list-style-type: none"> Inform the Manufacture and ask for repairs
The movements stop immediately after a longer adjustment time	<ul style="list-style-type: none"> Activation of the circuit breaker in the transformer 	<ul style="list-style-type: none"> Put in a rest position and let the system cool down for about 20 to 30 minutes

8. TECHNICAL DATA

8.1 Size and weight

Bed fitted with rails:

	Apollo – Saturn Rotate
Mattress base	88 x 210 cm
Overall dimensions	105.2 x 232.4 cm
Weight capacity	180 kg

8.2 Electrical data

Operating device (Main Control Units):

Type	Linak CB6S Open Bus Control Unit
Input voltage	AC 230 V, ± 10 %, 50 Hz
Classification	Class I

Hand control:

Type	Linak
------	-------

Electric motors for height adjustment of the mattress base backrest, footrest, rotate & tilt:

Type	Linak Linear Actuators
Max speed at no load / full load (mm/s)	38 / 23.2

8.3 Environmental conditions

Noise during adjustments: max. 48 dB(A)

The following environmental conditions must be respected:

Storage:

	Minimum	Maximum
Storage temperature:	+ 5° C	+ 50° C
Relative humidity:	50 %	70 %

Functioning:

	Minimum	Maximum
Environmental temperature :	+ 10° C	+ 40° C
Relative humidity:	20 %	90 % (without condensation)
Air pressure:	700 hPa	1060 hPa

8.4 Norms and applied guidelines

- Directive CE 93/ 42 Medical Devices Directive
- VDE 0751-1:2001-10 Safety verification for electrical devices
- EN 12182:1999-11 Technical aids for disabled persons
- DIN 32977-1:1992-07 Appropriate design for disabled persons
- EN 14971: 2001 Application of risk management to medical devices
- DIN EN 1970 : 2005-10 Adjustable beds for disabled persons
- EN 60601-1: 1996-03 Safety of electrical medical equipment
- EN 60601-1-2:2002-10 Electromagnetic compatibility
- DIN EN 60601-2-38/A1: 2001 Nursing beds electrically controlled (if required for home care)
- Classification as class I medical product (according to MPG § 13)

9. WARRANTY

The beds have a 5 year warranty period for any production faults from their purchase date. Keep the purchase receipt. In particular, the guarantee of the beds will NOT cover:

- Damages due to improper use by the user;
- Damaged or ruined parts due to wear and tear or due to excessive exposure to sunlight;
- Damages caused by loss or spillage of liquids;
- Damages from collisions or impacts against walls, furniture etc.;
- Tears and accidental damages, cigarette or open flames burns etc.

Warranty Terms

- The warranty is not transferable;
- Repairs must be carried out by an authorized dealer. Any repairs or alterations made by someone other than the dealer will void the warranty;
- The warranty of any replaced components will have the same expiry date of the original warranty;

The warranty will be void in the following cases:

- Damages caused by lack of maintenance or use of non-approved accessories;
- Damages caused by incorrect maintenance;
- Modifications or repairs carried out by someone other than the dealer;
- Modifications or repairs with parts not original.

If the device has problems, please contact your local authorized dealer and give all the information necessary to identify the defect and the type of repair. The manufacturer reserves the right to change weights, measures or other technical data contained in this manual without notice. The information contained in this manual is approximate.



10. DECLARATION OF CONFORMITY

The Company

Apollo Healthcare Technologies Ltd
Holme Street
Liversedge
WF15 6JF

Declares that

The articles which it sells, called and related accessories, are classified as Class I non-invasive Medical Devices according to Medical Directive 93/42/EEC and to national law no. 46/97 and comply with the European Directive on CE marking; in addition

- they are not measuring instruments;
- they are not intended for clinical investigations;
- they are sold in non-sterile packaging;

The above mentioned devices meet with the essential requirements and provisions of Directive 93/42/EEC.

Apollo Healthcare Technologies Ltd
Holme Street
Liversedge
WF15 6JF
Tel: 01924 614567
Fax: 01924 607480
Email: info@apollo-ht.co.uk

