Active Hands Looped Exercise Aids are padded to reduce chafing.

Your exercise aids have been sent to you fastened in the position they will be in when you wear them.

To put on the aids please follow these instructions whilst referring to the corresponding diagrams:

1. Undo the Velcro fixing on the wrist strap. Loosen the strap but do not pull the restraint through the plastic loop (Diagram 1).
2. Slip the aid over your hand so that the webbing loop lies in the palm of your hand and the plastic loop lies on the opposite side of the wrist to the thumb (Diagram 2).
3. With the aid of the ring, tighten the wrist strap and fasten across the back of your wrist with the Velcro.
4. Slip the loops over the handles of your exercise machine to use (Diagram 3).
5. Do not try to use your hands to hold on – use the loops instead of your hands.

## Precautions

If you are unsure whether this gripping aid is suitable for you, ask your Doctor, Physiotherapist or Occupational Therapist before use.

Always follow the instructions on this leaflet when putting on and using the gripping aid.

Active Hands gripping aids have been designed to assist with the gripping of objects which you may have been unable to use for some time. Therefore, care should be taken not to overstrain muscles (particularly in the wrists) when using heavy or awkward objects.

Do not use the gripping aid for prolonged periods, particularly if you suffer from problems with circulation.

If the gripping aid becomes damaged, do not use.

When using the gripping aid to perform a new task it is advisable to have another person present.

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