

Abduction Wedges

Posey Abduction Wedges are lightweight, polystyrene filled wedges for full leg and above the knee abduction. They feature a unique construction of an outer breathable vinyl mesh filled with polystyrene beads. This allows air to flow through the wedge, thereby minimizing the discomfort of conventional foam filled products, which do not breathe and can build up heat. And 4" (10 cm) wide, soft foam fabric straps with hook closures keep the wedge secured in place between the knees.



Small Abduction Knee Wedge

Allows abduction and alignment of the thighs while in a seated position. Supplied with one soft foam fabric strap with hook fasteners.

REF 6302S

Full Leg Abduction Wedge

Supplied with two soft foam fabric straps with hook fasteners.

REF 6302L

Monitoring

⚠️WARNING Be sure to follow your facility's policies and guidelines for frequency of patient monitoring. Posey recommends that this product be removed at least every two (2) hours to check for skin integrity, proper circulation and range of motion.

Inspection

⚠️WARNING Inspect before each use: check for broken stitches or parts; torn, cut or frayed material; or buckles that are cracked or broken and do not hold securely. DO NOT use soiled or damaged products.

Application Instructions

1. Position the patient as desired.
2. Place the wedge on the patient ensuring that there is no discomfort.
3. Position the Abduction Wedge in place with the soft foam strap(s) and secure with the hook closures sewn into the strap.
4. Reposition the patient as necessary.

Washing Instructions

Machine wash in warm temperature. Air dry preferred.



Disposal

⚠️WARNING Properly dispose of the product per facility's policy for BIOHAZARDOUS materials.

Storage and Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.



6302S



6302L

DESCRIPTION	POSEY ITEM #	DIMENSIONS		
		LENGTH	THICKNESS	WIDTH
Knee Abduction Wedge	6302S	8"	5"	7" tapering to 5"
		20 cm	13 cm	18 cm tapering to 13 cm
Full Leg Abduction Wedge	6302L	22"	5"	15" tapering to 7"
		56 cm	13 cm	38 cm tapering to 18 cm

Materials used meet the requirements of FMVSS-302 and CA #117 Section E for Flammability.