

## Accessories

A unique feature of the Netti wheelchairs in the vast amount of accessories available.

### Head supports

We have developed three kinds of head supports to meet the various needs of the users. These fit all Netti chairs. The head supports are attached to a bracket on the back of the chair and can be adjusted vertically, horizontally and in angle.



**Head support A (with side support)**  
This head support with side support, both in standard and Mini version, is good for users with reduced stability head control.



**Head support B (Small)**  
This head support is good for users who need some support in a reclined/tilted position.



**Head support C (Large)**  
This support give room for many head positions.



**Vertical stems.** The vertical stems come in three different lengths 200, 300 and 500 mm.



**Horizontal stems.** The horizontal stems come in three different lengths 150, 250 and 380 mm.



**Removable cover for improved hygiene.**



**Head cushion Butterfly.** This is a comfortable soft cushion that can be used as a soft support between the head support and the neck. The cushion is supplied with a strong rubber band for easy attachment.



**Comfort head cushion with kospflex.** Can be mounted on the existing head rest A and C.



**Head support fixation bracket**  
Offers extra stabilization of head support for spastics and other cases where the head support is exposed to hard usage.



**Head support racket – Netti Svipp**  
Allows mounting of head support on the Netti Svipp chair with foldable seat. Easy to mount and remove.

### Belts and straps

Alu Rehab has a number of different belts and straps. The sortiment includes H-Belts and Hip Belt with a choice of padded or not padded plastic lock or a car belt lock. The length of the belt is standard but can easily be cut to the required size after mounting.



**H-belt simple with plastic lock.**



**H-belt upholstered with car belt lock.**



**H-belt attachment bar.** It is fitted onto the push handles and is adjustable in height. It is telescopic and comes in two sizes.



**Abduction belt.**

**Simple hip belt with car belt lock or plastic lock.**



**Upholstered hip belt with car belt lock or plastic lock.**

*Item numbers and prices for the accessories you will find in the export pricelist*

## Seat and back unit

All Netti wheelchairs except for Netti Svipp and Netti 4U have a solid seat plate – comfort plate or sliding plate. The configuration varies between the markets. The comfort plate has a cut hole at the rear covered with a strong mesh fabric. The hole allows the Ischii Tuberositas (seat bones) to sink into the cushion and get stabilized on front. This reduce the tendency of sliding forward in seat.



Seat plate Comfort.



**Adjustment of seat depth**  
Length of thighs might vary to individual differences or to surgical operation of the hip etc. This can easily be accommodated by pulling out the foot support in front on either side of the wheelchair. This ensures a correct movement on the foot.



The depth of the chair can be adjusted on all Netti wheelchairs. On the fixed frame chairs, the depth is always adjusted at the rear. This is ergonomically the best principle for adjusting depth. Avoiding increased load over the front castors and correct proportion to the foot support is achieved. All Netti chairs can be adjusted 10 cm with standard parts.

By changing some components even longer seat depths can be achieved. The minimum seat depth varies from model to model.




Sliding seat plate.



**Back cushions.** All the seat cushions and nearly all back cushions are removable and can be individually fitted to the size and shape of each user.



**Lumbar support and wedge cushion**  
The lumbar support can be used in addition to the standard backrest to increase the curvature of the integrated lumbar support. The wedge cushion is used for increased sideways stability. It can be used on either or both sides. The cushions are placed between the frame and backrest. These are cost effective solutions.



**Velcro straps.** All Netti chairs have Velcro straps as standard.



**Aluminium back**  
This solution replaces the Velcro back on the Netti wheelchair giving a fixed and more solid base for the backrest cushion.



Wedge cushion.



**Back extender.** For tall users or users who need extra high back support a back extender can easily be mounted.

*Item numbers and prices for the accessories you will find in the export pricelist*

## Side support


For users with heavily reduced stability a stronger side support may be necessary. These support are attached to the backrest tubes and have different adjustability.



**Side support Mini**  
This side support rotates around the backrest tube before fixation and it is adjustable in depth and height.



**Side support Correction**  
This is a correctional side support. It is meant for users who need support directly on the body. Best performance is achieved in combination with the Flip backrest cushion as shown in the picture. Good solution for flexible scoliosis. It is adjustable in height, depth in combination with rotating pad.



**Vital Base Integral**  
For some users sliding out of the chair is a problem. Extra pelvis support will prevent this. Vital Base Integral is a separate plate that can be used on all kinds of wheelchairs. The Vital Base Integral pelvis supporter will have a very positive effect on sling seat wheelchairs as it prevents the bad hammock effect. The plate is put into cushions and can easily be fitted to different seat widths.



**Side support Stabil**  
This is a stabilizing side support. It is meant for users with reduced stability in the upper trunk. Used in combination with Smart, Stabil or Super Stabil backrest a very stable seating position is achieved, It is adjustable sideways and in height and has a very safe and user friendly locking system.

## Abductions block



The abduction block is meant to reduce abduction of the knees. For some users the abduction block may also prevent subluxation of the hip.

The abductions block can be delivered in small, medium and large sizes.

## Adductions support



Knee and thigh support reduces adduction.

## Arm support



Arm support – Netti Mini.



Arm support std for Netti 4U CE, CES and alternative for CED.



Arm support with cloth protector for Netti Svipp.



Arm support standard.



Arm support for Netti 4U CED.

Item numbers and prices for the accessories you will find in the export pricelist

## Arm support continue



**Arm support pads**  
A vast number of arm support pads are available to meet individual demands. These can be used on all kinds of Netti chairs.



**Comfort pads**  
In some cases it is beneficial to have the possibility to adjust the width of chair somewhat. This may be when the user needs a wide back rest, but is more narrow over the hips. On all Netti chairs the seat width adjusted by using pads on the skirt guards.

3 sizes are available 2,5, 3,5 and 4,5 cm



**Arm support pad curved**



Comfort pads for Netti Mini.



Comfort pads for Netti 4U CED.



**Hemiplegic arm support**  
The hemiplegic arm support makes the wheelchair suitable for e.g. stroke patients.

## Driving handles



Driving handles can be swung outwards and are delivered in 30 cm or 50 cm length.




**Hemiplegic cushion.**  
This cushion gives freedom to move the affected arm in a required position. It is mounted on a standard arm support.

## Release handles for tilt and recline

In case the user wants to be able to adjust the recline and tilt function themselves, the release handle can be moved from the back of the chair to the arm support. This can be mounted on both right and left side.

## Push bow



Push bow can be delivered to all Netti wheelchair sizes.



Release handle mounted on the arm support.



Release handle on the push bar.

*Item numbers and prices for the accessories you will find in the export pricelist*

## Trays and upholstery



Useful universal trays for all widths of the chairs can be delivered. The standard tray comes without upholstery. Upholstery comes in two different sizes. A reading tray can be mounted to the universal tray.



Universal tray comes in three different sizes.



**Swing able tray**  
It is swing able with the possibility to hang on the side of the chair. It does not conflict with the main wheels in this position. A solid tray with lock.



**Half tray**  
Swing able upholstered half tray. This tray replaces the existing standard arm support.



**Tray with lock**  
This tray has the same constructions as the swing able. This has to taken on and off between use.



**Tray upholstery**  
Can be delivered in two different sizes.



**Reading tray**  
Reading tray and universal tray used together.

## Foot support

Foot support can be delivered in two different models as angle adjustable and universal which is angle adjustable in fixed positions. Both models are swing and removeable and easy to assemble on and off.



Angle adjustable foot support.



Universal foot support.

*Item numbers and prices for the accessories you will find in the export pricelist*

## Foot support continue



Foot plates are delivered in correct widths according to wheelchair sizes. An easily operated lock between the plates ensures a solid foundation for the feet.



Foot board with lock. A very rigid footplate which is easy to mount and take off. It is easily removed by using a locking mechanism on the right side.



Calf supports are depth and height adjustable.



Knee upholstery.



Upholstery for footplate.



Calf support hinged. This is possible to move back under the chair by using a release handle. This makes removing the foot support much easier.



Extra upholstery for the calf support bracket prevents pressure points.



Extension piece for foot support.



Foot boxes can be delivered to each size of the chairs.

## Abduction front



Heel straps in different sizes.



Places the foot support 2.5 cm out from the frame, allowing room for abducted knees.



Amputation support.

## Brakes



Users with reduced mobility and weakened function in the fingers may need longer brake handles for easier function of the brakes.



The brakes can be mounted on the upper or bottom frame tube and there are several alternatives.

*Item numbers and prices for the accessories you will find in the export pricelist*

## Main wheels

The Netti wheelchairs have main wheels with solid tires in different sizes. Quick release axles enable smooth assembly and folding of the chair. Standard width of the tires is 1" but 1 3/8" is another option.



**Main wheel with one hand drive**  
20"-24" main wheels can be delivered with one hand drive.



**Main wheel with companion brake**  
Companion brakes are working directly to the hub of the main wheel with brake handles attached the driving handles. Companion brakes are delivered in complete set from 12" to 26" main wheel.



Quick release on main wheel.



**Spoke protectors**  
Spoke protectors can be delivered in black or transparent.



Frame extender makes it possible to mount the main wheels further back.



Camber adapter. Can be delivered in 2° or 4°.



Spoke protectors can also be delivered with different motives such as Lightening, Blue Bat, Feet, Tiger, Energy and Spider.



Main wheel adapter for Netti Svipp.



Push rims.



Distance block – front castors  
Gives room for the feet under the chair.

Item numbers and prices for the accessories you will find in the export pricelist

## Front castors

Quick release axles enable smooth assembly and folding of the chair.



Different front castors.

## Anti tip

The anti tippers except for Netti 4 U are adjustable in length and height. It can easily be turned around and functions like a tip bar when passing small steps like pavement etc.



Front forks. Different front forks to be used with the front castors.



To adjust the height of the anti tip wheel, remove the three screws holding the adjusting tube.

All tubes are adjustable in height.



Anti tip long/short.

Anti tip long/long.



The anti tip for Netti 4 U can be adjusted in two fixed positions. Short position is for 12", 16" and 24" main wheel in upper position. The long position is for 24" main wheel in lower position.



Anti tip X long/short.

Anti tip X long/long.



As accessory anti tip with pedal can be delivered to Netti 4 U.

## Mounting

Netti I/II/III, Svipp and Mini		
Description	Main wheel	Item nr.
Short/short, Mini	20" - 24"	25416*
Short/long, Mini	26"	25417
Long/short	20" - 24"	25419*
Long/long	26"	25418**
X long/short	20" - 24"	25422***
X long/long	26"	25423***

\* Standard anti tip for standard wheelchairs and for Netti Mini with frame extender

\*\* Netti Mini with frame extender

\*\*\* For wheelchairs with frame extender

Item numbers and prices for the accessories you will find in the export pricelist





### Jap unit

Can only be mounted on chairs with 1 3/8" main wheels.

Jap unit without main wheels.



Jap unit with main wheels.

### Rack for aspirator and respirator



Can be fitted on Netti III in sizes 35cm - 55cm.

### Intravenous bar



Can easily be adjusted from 1 to 4 hooks.

### Service kit



Service kit is a suitcase filled with all the screws, nuts and other small spare parts that are needed for to do service on a Netti wheelchair.

The service kit is delivered in a practical suitcase.

### Netti wheelchair bag



For mounting on all kinds of Netti wheelchairs.



### Electrical functions

Some users may benefit from electrical tilt, recline and elevating foot support. This enable the users to change the positioning and vary the sitting position themselves. Electrical functions can be supplied on Netti I Comfort, Netti II and Netti III Comfort chairs.

### Supportive engines

Most of the models of Netti I, Netti I comfort, Netti II and Netti III comfort can be mounted with little Viking assertive engine, F-16 assertive engine, E-Motion, E-Fix, Via-Mobile and Scala Mobil staircase climber. Alu Rehab does not sell these products but kits to prepare mounting of the engines.

Electrical functions for tilt, recline and elevating foot support.



The Netti III with electrical solutions comes with an integrated battery charger.



TGA with twin wheels.



### Wheelchair Powerpack

Alu Rehab is selling this solution for the Netti range wheelchairs. It is the most technically advanced, easiest to fit, simplest to use power pack ever. The standard single wheel version is capable of transporting up to 115 kg user weight. The standard twin-wheel model has the same user weight limit but much greater traction – really effective on loose gravel and wet grass. The heavy duty version is only available with twin and is capable of transporting up to 165 kg.



Item numbers and prices for the accessories you will find in the export pricelist