

Sito personal lift

Two Lifts in one

As this motto suggests, the new Sito personal lift unites the functions of an elevating lift and an erector lift and is thus an extremely flexible and versatile aid, giving you a very good relationship between price and performance.

Built in simple units, it can be transformed quickly and without tools into an elevating or an erector lift. The base frame remains unchanged. Simply order the lifting or upright arm in addition.

General:

Equipped with electrically adjustable elevating and/or erector function and a manual width adjustment of the chassis

Mechanical emergency lowering

Elevating lift:

Enables a comfortable raising of people, whether from a lying position in bed or a sitting position in the wheelchair

It is also possible to lift someone from the floor

Helps the nursing staff in many situations and protects the back



Erector lift:

Helps the patient stand up

Supports the nursing staff in daily tasks with the patient

Offers the necessary safety through careful movements

Enables the patients to stand up without danger to the back

Makes a standing transfer possible

Makes it better and easy to dress and change people, e.g. for toilet use, nappy change, etc.

Supports standing exercises for daily cardiovascular training and for strengthening of the musculature and bones



Lifting belts:

The Lifting belt comes in two different materials depending on application:

- universal belt made of nylon for everyday use, e.g. for the transfer from the bath to the wheelchair
- lifting belt for the bath and WC area

The shape and materials are determined by the needs of the patients

Easy to put on and off

Offer a pleasant, upright seating position to the patients



Upright belt:

in combination with the erector lift

Fast and simply to adapt to the patients

Offers the necessary support and security





Lifting arm conversion kit

Transforms from the erector lift to the elevating lift



Upright arm conversion kit

Transforms from the elevating lift to the erector lift without tools



	Technical dimensions
Load carrying capacity	max. 160 kg
Height of lift (starting from bracket, without belt)	min. 44 cm, max. 177 cm
Overall height (arms above)	207 cm
Dimensions of chassis	Length: 123 cm Width: 63.5 cm (92 cm, with spread chassis) Height: 11 cm
Emergency lowering	mechanical
Rechargeable battery	2 at 12 V, 7.2 Ah, closed lead gel rechargeable batteries
Hoist motor	24 V, 4 A
Material	Steel, powder coated, RAL 9006
Safety class	IPX 6
Wheels	Standard in front: 75 mm double castors Standard at the back: 75 mm single castors, lockable (alternatively 100 mm)
Total weight	Sito elevating lift: 42,5 kg Sito erector lift: 49 kg
Unit weight of the lifting arm	6,5 kg
Unit weight of the upright arm	5 kg
Unit weight of the shin pads plus footplate	7 kg
Unit weight of the rechargeable battery	6 kg