



## Fealite

### Nasal Pillows System (BMC-PM)

#### Less Weight, More Comfort

Fealite is much lighter than general nasal masks. It is quieter and offers minimal contact to the face, providing a more comfortable experience and making therapy easy to accept.



- **Best fit angle on pillows**

Adjustable pillows provides the most comfortable angle and maintains stability



- **360° rotatable elbow**

Flexible tube position makes minimal effect to sleep



- **Face matching headgear**

Relaxing

- **Anti-choke vent design**

Provides fail-safe operation

- **Diffuse vent hole**

Gentle and quiet exhaust

- **Dual direction adjustable headgear**

Ensure the most comfort and stability

- **Packed with 3 sizes pillows - S, M, L**

Simplify the purchase



- **Light and soft tube retainer**

Provides stable use regardless of tube orientation

- **Breathing Headgear**

Ventilative fabric ensures comfort throughout the night

Manufacturer: BMC Medical Co., Ltd.

NOTE: Specifications subject to change without prior notice.

20140225