

OPERATIONAL - QUALITY CHECK BEFORE SHIPMENT

ENSURE THE UPPER FRAME OPENING AND CLOSING MOVEMENTS EASILY OPERATES

ENSURE BOTH LOCKING PLUNGERS OPERATE SMOOTHLY AND LOCK THE UPPER FRAME INTO OPEN AND CLOSED POSITIONS

ENSURE THE SEAT OPENING AND CLOSING MOVEMENTS EASILY OPERATES SILENTLY THERE SHOULD NOT BE ANY FRICTION BETWEEN THE STUB TUBES AND INSIDE THE TRACKS

ENSURE THE TELESCOPIC GAS STRUT OPENS AND CLOSES THE REAR FRAME SMOOTHLY AND THE BRAKE CABLES ARE CLEAR IN THIS OPERATION BETWEEN THE CARRY HANDLE FRAME AND LOWER FRAME

ENSURE BOTH LOCKING PLUNGERS OPERATE SMOOTHLY WHEN OPENING AND CLOSING THE CARRY HANDLE



ENSURE THE BRAKE HANDLE OPERATES BOTH BRAKE CABLES SMOOTHLY AND CHECK EACH RUBBER TRACK RUNS FREELY OR HELD IN THE STOP POSITION

ENSURE THE 5 POINT SAFETY HARNESS WORKS CORRECTLY

ENSURE THE SEAT BELT FIXING POINT HAS SECURING WASHERS FITTED

ENSURE THE SAFETY LIGHT IS SECURE AND FACING FORWARD

ENSURE BOTH REAR CASTORS SWIVEL AND BRAKE

ENSURE THE SPRINGS ARE IN POSITION AND WORK EASILY WITH THE BRAKE ASSEMBLY

ENSURE BOTH QUICK RELEASE WHEELS ARE SECURE AND ROTATE SMOOTHLY

SURFACE FINISH - QUALITY CHECK BEFORE SHIPMENT

ENSURE ALL FABRIC PARTS FIT CORRECTLY AND DO NOT HAVE ANY TEARS OR LOOSE STITCHING

ENSURE ALL SAFETY BELT FABRIC PARTS FIT CORRECTLY AND DO NOT HAVE ANY TEARS OR LOOSE STITCHING

ENSURE ALL HEAT SHRINK COVERS ARE FITTED CORRECTLY AND DO NOT HAVE ANY TEARS OR HEAT DAMAGE ON THE SURFACE

ENSURE ALL BLACK PAINTED TUBES ARE FINISHED TO A HIGH STANDARD AND DO NOT HAVE ANY SCRATCHES OR SURFACE DEFECTS



ENSURE THE BRAKE HANDLE PART IS FITTED WITH THE NEW COVER AND IT IS GLUED INTO POSITION CORRECTLY.

ENSURE ALL SILVER PAINTED FRAMES AND TRACKS ARE FINISHED TO A HIGH STANDARD AND DO NOT HAVE ANY SCRATCHES OR DEFECTS

ENSURE THE EMBROIDERY IS CLEAR AND IS POSITIONED CORRECTLY ON THE COVER

ENSURE THE SAFETY LIGHT IS RED IN COLOUR

ENSURE BOTH RUBBER TRACKS ARE CLEAN AND WITHOUT ANY DAMAGE

ENSURE BOTH REAR CASTORS ARE CLEAN AND WITHOUT DAMAGE

ENSURE BOTH QUICK RELEASE WHEELS ARE CLEAN AND WITHOUT ANY SURFACE SCRATCHES

PACKAGING - QUALITY CHECK BEFORE SHIPMENT

ENSURE THE CARDBOARD BOX IS THE CORRECT SIZE AND CLEARLY CONFIRMS WHICH MODEL IS INSIDE THE BOX (EC1 / EC2)

ENSURE THE CARDBOARD BOX IS FREE FROM ANY SURFACE DAMAGE

ENSURE THE EVACUATION CHAIR INSIDE HAS PROTECTIVE PACKAGING AND IS A SNUG FIT INSIDE THE BOX

ENSURE THE WHITE DUST COVER – WALL BRKT – TRAINING DVD ARE PLACED INSIDE THE BOX BETWEEN THE SEAT FRAME AND REAR SQUAB ASSEMBLY