



# MiniLift Sit-to-Stand

A Unique Solution for Active Sit-to-Stand

# Extremely Comfortable, Active Sit-to-Stand With Natural Movement Pattern

The MiniLift Sit-to-Stand Lift assists the user when rising from a sitting to a standing position as gently and comfortably as possible. The handles and the side supports offer multiple hand placement options for the user. The stable anti-slip foot plate is very low, and the lower leg supports as well as the lifting arm can be adjusted to fit every user. When the MiniLift is used in combination with the appropriate lifting accessories, the user is supported under the feet, at the front of the lower legs and behind the back, which provides for a safe and secure, yet active sit-to-stand procedure. The MiniLift safely moves the user forward and upward in a natural movement pattern and, at the same time, the leg muscles and balance are exercised. It is manufactured in steel for stability and strength and has a weight capacity of 350 lb (160 kg) for the MiniLift160 and 440 lb (200 kg) for the MiniLift200.

## MiniLift160



Active sit-to-stand with natural movement pattern.



New control box with lift counter and overload indicator. Electrical emergency stop button readily available on the control box.



Ergonomic and easy-to-use hand control with 2 (EM) and 4 (EE) buttons respectively.

## MiniLift200

Easy to maneuver

Easily adjustable lift arm

Quick guide on the lift

High capacity battery with long service life

Built-in charging

Lockable rear castors

Grip-friendly handles and side supports provide multiple hand placement options for the user

Compact format

Soft, adjustable lower leg support

Footplate with low step-in height and anti-slip surface

Easy rolling, maintenance-free castors

Electrical base widening



# MiniLift Sit-to-Stand Specifications

## Part Numbers

### Material

Powder coated steel

### Safe Working Load (SWL)

MiniLift160EM/EE: 350 lb (160 kg)

MiniLift200: 440 lb (200 kg)

### Dimensions (L x W x H)

#### MiniLift160EM/EE

Length: 36.2 in (919 mm)

Length, collapsed: 43.2 in (1098 mm)

User length, min-max: 55.1–78.7 in (1400–2000 mm)

Base width, inner dimensions: 20.9–30.4 in (531–773 mm)

Base width, outer dimensions: 26.3–37.8 in (669–959 mm)

Height: 38.6 in (980 mm)

Footplate height, for user: 2.3 in (59 mm)

#### MiniLift200

Length: 38.6 in (980 mm)

User length, min-max: 55.1–78.7 in (1400–2000 mm)

Base width, inner dimensions: 17.7–34.4 in (450–875 mm)

Base width, outer dimensions: 23.2–40.6 in (590–1030 mm)

Height: 39.4 in (1000 mm)

Footplate height, for user: 3.7 in (95 mm)

### Weight

MiniLift160EM/EE: 89 lb (40 kg)

MiniLift200: 108 lb (49 kg)

### Base Widening

MiniLift160EM: Manual

MiniLift160EE: Electrical

MiniLift200: Electrical

### Castor Size

#### MiniLift160EM/EE

60300010 / 60300012: 3.9 in (100 mm)

Rear, 60300011 / 60300013: 3.9 in (100 mm)

Front, 60300011 / 60300013: 2.8 in (70 mm)

#### MiniLift200

3.9 in (100 mm)

### Lifting Speed (Without Load)

1.4 in / sec (35 mm / sec)

### Safety

Emergency lowering: Manual and electrical



45600004-008

#### 60300010: MiniLift160EM

Mobile sit-to-stand lift with manual base widening

#### 60300011: MiniLift160EM

Mobile sit-to-stand lift with manual base widening, low legs

#### 60300012: MiniLift160EE

Mobile sit-to-stand lift with electrical base widening

#### 60300013: MiniLift160EE

Mobile sit-to-stand lift with electrical base widening, low legs

#### 401100334: MiniLift200

Mobile sit-to-stand lift with electrical base widening

#### 70200033: CalfStrap

For MiniLift160EM/EE and MiniLift200

#### 70200012: SlingBarWrap MiniLift

Soft, wipeable protection for the bar

#### 70200089: Hand control HB33-6

6 buttons, with service, battery status and overload indicators, for MiniLift160EE, 60300012 and 60300013

#### 70200029: External charger Linak

External charger for all models and variants of MiniLift

#### ThoraxSling, Polyester

45500003: X-Small

45500004: Small

45500006: Medium

45500007: Large

45500008: X-Large

45500009: XX-Large

#### ThoraxSling, Wipeable

45550004: Small

45550006: Medium

45550007: Large

45550008: X-Large

45550009: XX-Large

#### ThoraxSling, Disposable,

Non-woven, 5 pcs

45590004: Small

45590006: Medium

45590007: Large

45590008: X-Large

#### ThoraxSling, with Seat Support, Polyester

45600004: Small

45600006: Medium

45600007: Large

45600008: X-Large

#### ThoraxSling, with Seat Support, Disposable, Non-Woven, 5 pcs

45690004: Small

45690006: Medium

45690007: Large

45690008: X-Large

**ThoraxSling** is a specially adapted and comfortable sit-to-stand vest that is used in combination with Handicare's MiniLift sit-to-stand lifts. It is easy for the caregiver to apply, while giving the user an added sense of security during rising and standing training. ThoraxSling also facilitates toileting. ThoraxSling is available in two models; one which provides support behind the back (ThoraxSling) and one that also gives support under the seat (ThoraxSling, with seat support), and in several materials and sizes.



45500003-009



45500006-009



45600007-008



45690004-008



Are you ready to get started  
with the MiniLift?

Call 1-888-637-8155 or visit [www.handicareusa.com](http://www.handicareusa.com)