

We provide access!

**LEHNER**  
**LIFTTECHNIK**



# Lift Collection

Platform Stairlifts

Stairlifts

Vertical Platformlifts

Stairclimbers

# Evac Skate

## The evacuation chair

The Evac Skate is an evacuation device used to move people with a disability or injury down stairways quickly and safely during an emergency. It's patented speed governor and braking system allows a small attendant to easily evacuate a larger passenger. The tracks grip the stairs, regardless of the stair construction material.

The Evac Skate is quick to fold and easy to store in a special storage cabinet to be ready in a case of emergency



### Treppentrans (motorized option)

As an option to the standard evacuation chair a powerful motor allows a small attendant to safely transport a larger person up and down stairways



The standard evacuation chair works without the need for batteries or a motor. It allows the safe and easy transport of people down the steps in case of emergencies.



## Evac-Skate - Technical Specification

### Facts

Loading capacity: 135 kg


Weight of unit: 21 kg

Dimensions 1100 x 426 x 270

Max. inclination: 40°

Lifting speed: 0,8 m/sec

Minimum landing space



# Providing access around the world!



Over 15.000 lift installations have been realized over the past 25 years which provide accessibility solutions to all kinds of public and private needs around the world!

## Contact us!

Don't hesitate to call us:  
+437278 35 14

Visit our website:  
[www.lehner-lifttechnik.at](http://www.lehner-lifttechnik.at)

Write us an email:  
[office@lehnerlifttechnik.at](mailto:office@lehnerlifttechnik.at)

Lehner Lifttechnik GmbH (© 2015)  
Salling 8, 4724 Neukirchen am Walde Austria

Call us: +65 62500555

Visit us: [www.Rehabmart.sg](http://www.Rehabmart.sg)

