

Kaye Large Bolster Chair with tray & casters

Large Bolster Chair with tray & casters

Kaye bolster chair		
Model	Small (B3)	Large (B4)
Age	2 - 6 yrs.	5 - 11 yrs.
Bolster size	15" Long x 8" Diameter	23" Long x 10" Diameter
Height from floor to seat	Adjusts 7" to 13"	Adjusts 9" to 17½"
Back height	7" above Seat	10" above Seat
Back extension height	Adds 8" to back height	Adds 11" to back height
Lateral pelvic supports	Adjusts 4½" to 6½" between cushions	Adjusts 5½" to 10" between cushions
Tray	25" Width x 23½" Depth	25" Width x 29½" Depth
User height	Up to 42"	Up to 48"
Weight capacity	55 lbs.	100 lbs.
Product weight	46 lbs.	54 lbs.

[Brochure](#)