

Kaye Bolster Chairs

Bolster Chairs are designed for children who need a wide base of support in order to achieve good pelvic position, lower extremity position, and trunk alignment for stable sitting. Two sizes are available, large and small.

The Kaye Bolster Chairs can be purchased with an adjustable Chest Support System that limits movements of the upper trunk to the range that the child can control, allowing the child to lean forward or to the sides to reach and play, yet regain the upright position with less effort.



B4C



B3C

B4CCS



B3CCS

Model Features	B3	B4
Age	2-6 years	5-11 years
Size & Weight	Up to 42" / 55 lbs (25 kg)	Up to 48" / 100 lbs (45 kg)
Bolster size	15"l x 8" dia.	23"l x 10" dia.
Height from floor to seat	Adjusts 7" to 13"	Adjusts 9" to 17.5"
Back height	7" above seat	10" above seat
Back extension height	Adds 8" to back height	Adds 11" to back height
Lateral pelvic supports	Adjusts 4.5" to 6.5" between cushions	Adjusts 5.5" to 10" between cushions
Tray	25"w x 23.5"d	25"w x 29.5"d
Chair with 4" swivel breaking casters	Model B3C	Model B4C
Chair with casters and chest support	Model B3CCS	Model B4CCS