











































































































| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | | ✓ | Description | |
|---|-------------|---------------------------------------|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|---------------|------------------------|------------------------|---|---|---------|
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c- 22 | SB6c- 26 | SB6w | SB7e | P1 | P2 | | | |
| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveller | P1 Pediatric Converter | P2 Pediatric Converter | | | |
|  | 1 | Caster Wheel | | | | | | | | | | | | | 5" (125mm) Secure Caster Wheel | \$87 |
|  | 500-265 | Caster Wheel | | | | | | | | | | | | | 4" (100mm) Robby Brand Caster Wheel Model: Kelly / ORANGE + GREY | \$94 |
|  | 500-383 | Caster Wheel | | | | | | | | | | | | | 4" (100mm) Robby Brand Caster Wheel Model: Kelly / ORANGE + GREY, identical to SB6c but given a separate part number | \$94 |
|  | 500-320 | Caster Wheel | | | | | | | | | | | | | 4" (100mm) Generic Brand Caster Wheel with brake | \$23 |
|  | 100-500 | Caster Wheel | | | | | | | | | | | | | 4" (100mm) Generic Brand Caster Wheel with brake, fitted to leg tube | \$30 |
|  | 500-367 | Rubber Foot | | | | | | | | | | | | | Rubber buffer foot | \$18 |
|  | 2 | Height Adjustment Tube (Caster Wheel) | | | | | | | | | | | | | Caster Wheel connects to this tube | \$12 |
|  | 3 | Tube Cap (25mm) | | | | | | | | | | | | | | \$2 |
|  | 100-139 | Tube Cap | | | | | | | | | | | | | | \$2 |
|  | 4 | Rolling Track | | | | | | | | | | | | | Includes: 1 x Secondary Brake (Left) (Part 5) 1 x Secondary Brake (Right) (Part 6) 1 x Primary Brake Latch (Part 7) 2 x Latch Hook (Part 8) | \$234 |
|  | 5 | Secondary Brake (Left) | | | | | | | | | | | | | Include: 2 x Aluminium Rivets | \$40 |
|  | 6 | Secondary Brake (Right) | | | | | | | | | | | | | Include: 2 x Aluminium Rivets | \$40 |
|  | 7 | Primary Brake Latch | | | | | | | | | | | | | Include: 2 x Aluminium Rivets | \$5 |
|  | 8 | Latch Hook | | | | | | | | | | | | | Include: 2 x Aluminium Rivets | \$5 |
|  | 9 | Securing Clip (Medium) | | | | | | | | | | | | | Medium length securing clips are used to attach the caster wheels on the chair rolling base unit only | \$4 |
|  | 10 | Rolling Base Leg (Left) | | | | | | | | | | | | | Includes: 4 x Leg Tube Cap (Part 12) | \$60 |
|  | 11 | Rolling Base Leg (Right) | | | | | | | | | | | | | Includes: 4 x Leg Tube Cap (Part 12) | \$60 |
|  | 12 | Leg Tube Cap | | | | | | | | | | | | | | \$2 |
|  | 13 | Rolling Base Track Wheel | | | | | | | | | | | | | | \$4 |
|  | 14 | Tube Washer (32mm) | | | | | | | | | | | | | | \$2 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c- 22 | SB6c- 26 | SB6w | SB7e | P1 | P2 | | | #VALUE! |





















| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|------------------------------|---|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|--|---------|
|  | 15 | Track Wheel Axle (Front) | | | | | | | | | | | | Includes: 2 x Stainless Steel washers 1 x Nyloc Nut | \$22 |
|  | 16 | Track Wheel Axle (Rear) | | | | | | | | | | | | Includes: 2 x Stainless Steel washers 1 x Nyloc Nut | \$22 |
|  | 17 | Primary Brake Lever | | | | | | | | | | | | Includes: 2 x Stainless Steel washers 1 x Spring + Handle + Stopper Pin | \$13 |
|  | 18 | Primary Brake Handle Sleeve | | | | | | | | | | | | | \$2 |
|  | 19 | Foot Rest (Seat Unit) | | | | | | | | | | | | Includes: 12 x Foot Rest Grips (Part 20) | \$127 |
|  | 20 | Foot Rest Grips | | | | | | | | | | | | Packed in sets of 12 pieces | \$27 |
|  | 21 | Rolling Base Height Adjustment Tube (Front) | | | | | | | | | | | | Refer to Accessories if a higher seat height is required to clear a high bath side | \$42 |
|  | 22 | Rolling Base Height Adjustment Tube (Rear) | | | | | | | | | | | | Refer to Accessories if a higher seat height is required to clear a high bath side | \$42 |
|  | 23 | Tube Washer (25mm) | | | | | | | | | | | | | \$2 |
|  | 24 | Fixing Knob | | | | | | | | | | | | | \$10 |
|  | 500-078 | Hex Bolt (replaced earlier Fixing Knob) | | | | | | | | | | | | Holds tubular seat frame to main leg frame | \$2 |
|  | 25 | Securing Clip (Large) | | | | | | | | | | | | Large securing clips are used on the chair rolling base unit only, for height adjustment. | \$4 |
|  | 26 | Seat Base | | | | | | | | | | | | Includes: 2 x Arm Rest Clip (Part 27) 4 x Cushion Connector (Part 28) 2 x Tube Cap (32mm) (Part 28) 2 x Commode Bucket Connectors (Part 31) | \$258 |
|  | 27 | Arm Latch | | | | | | | | | | | | | \$13 |
|  | 500-237 | Arm Latch | | | | | | | | | | | | This part is for SB6 series chairs only and is not interchangeable with part #27 | \$13 |
|  | 500-117 | Arm Rest Lock | | | | | | | | | | | | | \$13 |
|  | 28 | Cushion Connector | | | | | | | | | | | | | \$3 |
|  | 29 | Tube Cap (32mm) | | | | | | | | | | | | | \$2 |
|  | 30 | PU Cushion Open Front | | | | | | | | | | | | PU Foam Standard Cushion. Push-on Fit to seat plate. | \$151 |
|  | 500-020 | PU Cushion Open Front | | | | | | | | | | | | PU Cushion with drip edge around commode opening. Comes with Commode Rails factory fitted | \$161 |
|  | 500-020-B (Accessory SCSF-6) | Blow Mould Cushion Open Front | | | | | | | | | | | | Blow Mould Soft Foam cushion, screw fixed to metal supports on frame. Comes with Closure Cushion and Commode Rails factory fitted. | \$141 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |























| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | \$0 | |
|---|---------------------------------|---|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|-------------|---------|
|  | 100-505 | Cushion Open Front, complete with closure cushion | | | | | | | | | | | | \$94 | |
|  | 500-319 | Blow Mould Cushion Closed Front 18" | | | | | | | | | | | | \$94 | |
|  | 500-376 | Blow Mould Cushion Open Front 22" | | | | | | | | | | | | \$144 | |
|  | 500-381 | Blow Mould Cushion Open Front 26" | | | | | | | | | | | | \$157 | |
|  | 500-323 | Blow Mould Commode Infill Cushion | | | | | | | | | | | | \$40 | |
|  | 500-166 | P1 Seat Cushion - PU foam | | | | | | | | | | | | \$161 | |
|  | 500-166-R | P1 Seat Cushion - Standard Issue - PU Foam | | | | | | | | | | | | \$174 | |
|  | 500-167 | P2 Seat Cushion - PU Foam | | | | | | | | | | | | \$161 | |
|  | 31 | Commode Pail Holder | | | | | | | | | | | | \$13 | |
|  | 500-080 | Packer Under Commode Pail Holder | | | | | | | | | | | | \$20 | |
|  | 500-080D | OLD Packer Under Commode Pail Holder | | | | | | | | | | | | \$20 | |
|  | 500-362 | Commode Pail Holder Rails (Upgrade Models) | | | | | | | | | | | | \$40 | |
|  | 500-329 | Commode Pail Wire Holder Cradle | | | | | | | | | | | | \$40 | |
|  | 32 | Commode Pail + Lid | | | | | | | | | | | | \$20 | |
|  | 500-072 (Pail) 500-073 (Lid) | Commode Pail + Lid = 500-076 | | | | | | | | | | | | \$54 | |
|  | 100-138 | Commode Pail + Lid | | | | | | | | | | | | \$23 | |
|  | 500-074 | Pee Guard | | | | | | | | | | | | \$27 | |
|  | 32C | Commode Splash Guard | | | | | | | | | | | | \$13 | |
|  | 33 | Arm Rest (Left) | | | | | | | | | | | | \$90 | |
|  | 500-043 | Arm Rest (Left) | | | | | | | | | | | | \$90 | |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | #VALUE! | |





















| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|-------------|--------------------------------------|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|---|---------|
|  | 100-501 | Arm Rest (Left) | | | | | | | | | | | | Armrest base and arm tube supplied on SB6c chairs from serial DC-SB6c-02325 and SB6w chairs from serial DC-SB6w-01519 (but excluding DC-SB6w-01546 to 01595 inclusive) | \$47 |
|  | 500-367 | Arm Rest (Left) | | | | | | | | | | | | Complete Armrest ASSEMBLY | \$120 |
|  | 34 | Arm Rest (Right) | | | | | | | | | | | | Complete Armrest | \$90 |
|  | 500-045 | Arm Rest (Right) | | | | | | | | | | | | Complete Armrest supplied on SB6c chairs up to serial number DC-SB6c-02325, and SB6w chairs up to serial DC-SB6w-01518 (but including DC-SB6w-01546 to 01595 inclusive) | \$90 |
|  | 100-502 | Arm Rest (Right) | | | | | | | | | | | | Armrest base and arm tube supplied on SB6c chairs from serial DC-SB6c-02326 and SB6w chairs from serial DC-SB6w-01519 (but excluding DC-SB6w-01546 to 01595 inclusive) | \$47 |
|  | 500-369 | Arm Rest (Right) | | | | | | | | | | | | Complete Armrest ASSEMBLY | \$120 |
|  | 500-264 | Arm Rest Pad PU Foam (Left) | | | | | | | | | | | | PU Foam Armrest Only | \$47 |
|  | 500-333 | Arm Rest Pad Plastic (Left or Right) | | | | | | | | | | | | Plastic Armrest Only | \$27 |
|  | 500-271 | Arm Rest Pad PU (Right) | | | | | | | | | | | | PU Foam Armrest Only | \$47 |
|  | 38 | Arm Rest Cap | | | | | | | | | | | | | \$2 |
|  | 39 | Arm Rest Knob | | | | | | | | | | | | | \$12 |
|  | 40 | Back Rest Frame | | | | | | | | | | | | | \$50 |
|  | 100-515 | Back Rest Frame | | | | | | | | | | | | | \$44 |
|  | 500-042 | Back Rest Frame 18" | | | | | | | | | | | | Included: 2 x EVA Foam hand grip (Part #500-023) | \$80 |
|  | 500-372 | Back Rest Frame 22" | | | | | | | | | | | | Included: 2 x EVA Foam hand grip (Part #500-023) | \$94 |
|  | 500-377 | Back Rest Frame 26" | | | | | | | | | | | | Included: 2 x EVA Foam hand grip (Part #500-023) | \$107 |
|  | 500-351 | Back Rest Frame | | | | | | | | | | | | | \$54 |
|  | 41 | Back Rest Fabric | | | | | | | | | | | | | \$60 |
|  | 500-026 | Back Rest Fabric | | | | | | | | | | | | | \$67 |
|  | 500-334 | Back Rest Fabric | | | | | | | | | | | | Pocket on Rear of Backrest Fabric for Infill Cushion Storage | \$60 |
|  | 100-510 | Back Rest Fabric | | | | | | | | | | | | Pocket on Rear of Backrest Fabric for Infill Cushion Storage | \$60 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|-------------|---|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|---|---------|
|  | 500-375 | Back Rest Fabric 22" | | | | | | | | | | | | Pocket on Rear of Backrest Fabric for Infill Cushion Stowage | \$74 |
|  | 500-380 | Back Rest Fabric 26" | | | | | | | | | | | | Pocket on Rear of Backrest Fabric for Infill Cushion Stowage | \$80 |
|  | 42 | Hand Grips | | | | | | | | | | | | | \$13 |
|  | 500-023 | EVA foam hand grip | | | | | | | | | | | | Comes straight and has to be fed over the Back Rest Frame (Part 40B). These are longer foam grips than other models to cushion the shoulder blade area of the user | \$20 |
|  | 43 | Bridge | | | | | | | | | | | | Includes: 2 x Bridge Tongue (Long) (Part 44) 2 x Bridge Tongue (Short) (Part 45) 4 x Latch (Part 46) For custom length bridges see accessories. | \$216 |
|  | 44 | Bridge Tongue (Long) | | | | | | | | | | | | | \$40 |
|  | 45 | Bridge Tongue (Short) | | | | | | | | | | | | | \$33 |
|  | 46 | Latch | | | | | | | | | | | | Includes: 5 Rivets, 1 latch | \$40 |
|  | 47 | Shower Base | | | | | | | | | | | | Includes: 2 x Secondary Brakes (Parts 5 + 6) 1 x Primary Brake Latch (Part 7) 2 x Latch Hook (Part 8) 4 x Leg Tube Cap (Part 12) Fits into part 47. Longer adjustment tubes (4" / 100mm) may be required for shower bases that are level or below the bathroom floor, refer to Accessories / code SBR | \$345 |
|  | 48 | Standard 63mm Height Adjustment Tubes (Shower Track) | | | | | | | | | | | | | \$20 |
|  | 49 | Suction Cups | | | | | | | | | | | | Fits into part 48 | \$23 |
|  | 50 | Securing Clip (Small) | | | | | | | | | | | | Small securing clips are used to attach part 48 to the shower track (part 47) only. | \$4 |
|  | 51 | Rolling Base Legs | | | | | | | | | | | | Includes: 4 x Leg Tube Cap (Part 12) | \$284 |
|  | 52 | Rolling Base Foot Rest | | | | | | | | | | | | Includes: 12 x Foot Rest Grips (Part 20) | \$127 |
|  | 53 | Foot Rest Attachment (Nut + Bolt) | | | | | | | | | | | | Includes: 2 x Stainless Steel washers 1 x Nyloc Nut | \$5 |
|  | 54 | Height Adjustment Tubes Standard Height (Rolling Base) Pair | | | | | | | | | | | | In some situations where the tub base is higher or lower than the bathroom floor accessory unit RBE24 or TBE24 may be required to correctly adjust the height of the units. See Accessories. | \$42 |
|  | 55 | Rolling Base Track | | | | | | | | | | | | Includes: 2 x Secondary Brakes (Parts 5 + 6) 1 x Primary Brake Lever (Left) (Part 56) 1 x Primary Brake Lever (Right) (Part 57) 2 x Commode Bucket Connectors (Part | \$284 |
|  | 56 | Primary Brake Lever (Left) | | | | | | | | | | | | Complete with spring and fixing screw | \$10 |
|  | 57 | Primary Brake Lever (Right) | | | | | | | | | | | | Complete with spring and fixing screw Includes: 2 x Arm Rest Clip (Part 27) 4 x Cushion Connector (Part 28) 2 x Tube Cap (32mm) (Part 29) 4 x Track Wheel + 4 x Axle (Part 13 + 15) 2 x Primary Brake Stopper (Part 59) | \$10 |
|  | 58 | Seat Base | | | | | | | | | | | | | \$291 |
|  | 59 | Primary Brake Stopper | | | | | | | | | | | | | \$20 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |


| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|-------------|--|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|---|---------|
|  | 60 | Bridge | | | | | | | | | | | | Includes: 4 x Bridge Tongue (Short) (Part 45) 2 x Latch (Part 61) | \$191 |
|  | 61 | Latch | | | | | | | | | | | | Includes: 4 x Stainless Steel Screws 4 x Nyloc Nuts | \$33 |
|  | 62 | Bath Base Track | | | | | | | | | | | | Includes: 1 x Secondary Brake (Left) (Part 5) 1 x Secondary Brake (Right) (Part 6) 1 x Primary Brake Lever (Left) (Part 56) 1 x Primary Brake Lever (Right) (Part 57) | \$284 |
|  | 63 | Height Adjustment Tubes (Bath Base) Pair | | | | | | | | | | | | In some situations where the tub base is higher or lower than the bathroom floor accessory unit RBE24 or TBE24 may be required to correctly adjust the height of the units. See Accessories. | \$42 |
|  | 64 | Bath Base Leg (Front) | | | | | | | | | | | | This part is the front leg of the bath tub base unit. It has the footrest attachment lugs welded to it PLUS can be identified from the height adjustment holes facing the front - rear (not sideways) | \$47 |
|  | 64A | Bath Base Leg (Rear) | | | | | | | | | | | | This part is the Rear leg of the bath tub base unit. It can be identified from the height adjustment holes facing the side (not front - rear) | \$47 |
|  | 65 | Bath Base Leg (Front SB2T) | | | | | | | | | | | | This part is the front leg of the bath tub base unit. It can be identified from the height adjustment holes facing the front - rear (not sideways) | \$47 |
|  | 65A | Bath Base Leg (Rear SB2T) | | | | | | | | | | | | Includes: 2 x Leg Tube Cap (Part 13) | \$100 |
|  | 66 | Foot Rest (Bath Base) | | | | | | | | | | | | Includes: 12 x Foot Rest Grips (Part 20) | \$107 |
|  | 67 | Bath Base side Stabilising Tubes | | | | | | | | | | | | | \$54 |
|  | 68 | Rubber Feet | | | | | | | | | | | | | \$20 |
|  | 69 | Rolling Base Legs | | | | | | | | | | | | Includes: 4 x Leg Tube Cap (Part 12) | \$157 |
|  | 70 | Neck Rest attached to round tube | | | | | | | | | | | | This part is the PU foam neck rest only, fitted to a round tube. Order the support structure separately as below. | \$79 |
|  | 70A | Neck rest attached to "D" profile tube | | | | | | | | | | | | This part is the PU foam neck rest only, fitted to a "D" profile tube. Order the support structure separately as below. | \$79 |
|  | 71 | Neck Rest Height Adjustment vertical round Tube | | | | | | | | | | | | Older style round tube superceded by part #71A plus you will then need part # 72A | \$37 |
|  | 71A | Neck Rest Height Adjustment Tube with "D" profile vertical non-rotating tube | | | | | | | | | | | | NOTE: This upgrade Neck Rest Height Adjustment Tube has a D shaped flat edge and is designed to partner with part #72a which has been moulded to suit this new tube. | \$37 |
|  | 72 | Neck Rest Attachment for part #71, round tube | | | | | | | | | | | | Includes: 1 x Neck Rest Handle (Part 74). NOTE: This part does need Neck Rest Sleeve part #73 | \$94 |
|  | 72A | Neck Rest Attachment with D Shape tube aperture to suit part #71A | | | | | | | | | | | | Includes: 1 x Neck Rest Handle (Part 74) NOTE: This part does not need Neck Rest Sleeve part #73 | \$94 |
|  | 73 | Neck Rest Sleeve, for round tubes | | | | | | | | | | | | Superceded by Part #73A | \$6 |
|  | 73A | Neck Rest Sleeve, for "D" profile tubes | | | | | | | | | | | | This part is required in every Neck Rest, with part #71 and part #71A at the top where part #70A fits. | \$6 |
|  | 74 | Neck Rest Handle | | | | | | | | | | | | Includes: 1 x Stainless Steel Locking Screw 2 x Stainless Steel Washer | \$30 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|-------------|--|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|--|---------|
|  | 75 | Tilt Release Handle | | | | | | | | | | | | Includes: 1 x Stainless Steel Screw 1 x Nyloc Nut 1 x Stainless Steel Washer | \$50 |
|  | 76 | Strut Cable | | | | | | | | | | | | | \$94 |
|  | 77 | Seat Base | | | | | | | | | | | | Includes: 2 x Arm Rest Clip (Part 27) 4 x Cushion Connector (Part 28) 2 x Tube Cap (32mm) (Part 29) 4 x Track Wheel (Part 13) 4 x Track Wheel Axle (Part 15) 2 x Primary Brake Stopper (Part 59) 1 x Strut (Left) (Part 78) 1 x Strut (Right) (Part 79) 2 x Rubber Stopper (Part 80) 2 x Tilt Stopper (Part 81) 2 x Primary Brake Stopper (Part 59) | \$823 |
|  | 78 | Strut (Left) | | | | | | | | | | | | NOTE: These Struts are NOT interchangeable with struts for the SB3T Chair. Refer to Part # 93, 94 for SB3T Struts | \$117 |
|  | 79 | Strut (Right) | | | | | | | | | | | | NOTE: These Struts are NOT interchangeable with struts for the SB3T Chair. Refer to Part # 93, 94 for SB3T Struts | \$117 |
|  | 80 | Rubber Stopper | | | | | | | | | | | | | \$10 |
|  | 81 | Tilt Stopper | | | | | | | | | | | | | \$3 |
|  | 82 | Foot Rest Connector (Left) | | | | | | | | | | | | Footrest connector is used in conjunction with Part # 84 and is fitted to SB2T chairs below serial #211605 and SB3T chairs below #311666 | \$28 |
|  | 83 | Foot Rest Connector (Right) | | | | | | | | | | | | Footrest connector is used in conjunction with Part # 85 and is fitted to SB2T chairs below serial #211605 and SB3T chairs below #311666 | \$28 |
|  | 84 | Foot Rest (Left) Tool Free Adjustment | | | | | | | | | | | | This original design footrest is fitted to SB2T chairs below serial #211605 and SB3T chairs below #311666 | \$124 |
|  | 85 | Foot Rest (Right) Tool Free Adjustment | | | | | | | | | | | | This original design footrest is fitted to SB2T chairs below serial #211605 and SB3T chairs below #311666 | \$124 |
|  | 500-254 | Footrest Assembly Right | | | | | | | | | | | | This upgraded design footrest is fitted to SB2T chairs from serial #211605 and SB3T chairs from #311666. Note that locking locator pin is 80mm down from top of tube (for identification purposes) | \$107 |
|  | 500-254-6 | Footrest Assembly Right | | | | | | | | | | | | Note that locking locator pin is 140mm down from top of tube (for identification purposes) | \$107 |
|  | 500-253 | Footrest Assembly Left | | | | | | | | | | | | This upgraded design footrest is fitted to SB2T chairs from serial #211605 and SB3T chairs from #311666. Note that locking locator pin is 80mm down from top of tube (for identification purposes) | \$107 |
|  | 500-253-6 | Footrest Assembly Left | | | | | | | | | | | | Note that locking locator pin is 140mm down from top of tube (for identification purposes) | \$107 |
|  | 500-254-P | Petite-Buddy Footrest Assembly Right | | | | | | | | | | | | This footrest assembly is shorter than the standard adult footrest giving a smaller popletial length between cushion and foot plate | \$107 |
|  | 500-253-P | Petite-Buddy Footrest Assembly Left | | | | | | | | | | | | This footrest assembly is shorter than the standard adult footrest giving a smaller popletial length between cushion and foot plate | \$107 |
|  | 86 | Foot Rest Paddle (Left) | | | | | | | | | | | | | \$40 |
|  | 87 | Foot Rest Paddle (Right) | | | | | | | | | | | | | \$40 |
|  | 500-266 | Foot Rest Paddle | | | | | | | | | | | | | \$27 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | Showerbuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | \$0 | |
|---|-------------|--|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|---|-------------|---------|
|  | 500-605 | Plastic End Cap on Footrest Paddle Tube | | | | | | | | | | | | \$2 | |
|  | 88 | Foot Rest Connector Knob | | | | | | | | | | | | \$12 | |
|  | 89 | Foot Rest Height Adjustment Knob | | | | | | | | | | | | \$10 | |
|  | 500-600 | Foot Rest Height Adjustment Clip | | | | | | | | | | | | \$5 | |
|  | 90CLP | Calf Strap | | | | | | | | | | | | \$47 | |
|  | 90 | Calf Strap | | | | | | | | | | | | \$47 | |
|  | 500-413/22 | Calf Strap | | | | | | | | | | | | \$30 | |
|  | 500-413/26 | Calf Strap | | | | | | | | | | | | \$30 | |
|  | 93 | Seat Base | | | | | | | | | | | | \$258 | |
|  | 94 | Strut (Left) | | | | | | | | | | | | \$117 | |
|  | 95 | Strut (Right) | | | | | | | | | | | | \$117 | |
|  | 96 | Rolling Base Height Adjustment Tube | | | | | | | | | | | | \$120 | |
|  | SD | Shower Discs | | | | | | | | | | | | \$13 | |
|  | 97 | Bridge Support | | | | | | | | | | | | \$124 | |
|  | 500-269-R | 7.6mm Push Pin | | | | | | | | | | | | \$2 | |
|  | 500-269-S | 7.6mm Push Pin | | | | | | | | | | | | \$2 | |
|  | 500-268 | Leg Height Adjustment Fixing Knob | | | | | | | | | | | | \$10 | |
|  | 500-044 | Plastic Cap | | | | | | | | | | | | \$6 | |
|  | 500-108 | Safety Harness Clamp Bar | | | | | | | | | | | | \$9 | |
|  | 500-001 | Leg Frame 18" for "C" and "W" Model Chairs | | | | | | | | | | | | \$187 | |
|  | 500-374 | Leg Frame 22" | | | | | | | | | | | | \$218 | |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | |  | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|--------------------------------|--|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|---|---------|
|  | 500-379 | Leg Frame 26" | | | | | | | | | | | | | \$224 |
|  | 500-268 | Footrest adjust knob | | | | | | | | | | | | | \$10 |
|  | 500-255 | Splash guard clip | | | | | | | | | | | | Screw fixed to underside of cushions | \$21 |
|  | 500-021 | Seat Frame 18" | | | | | | | | | | | | | \$157 |
|  | 500-373 | Seat Frame 22" | | | | | | | | | | | | | \$174 |
|  | 500-378 | Seat Frame 26" | | | | | | | | | | | | | \$181 |
|  | 500-079 | Seat Cushion Fixing Knob | | | | | | | | | | | | | \$10 |
|  | 500-107 | Seat Cushion Fixing Screw, Self Taping | | | | | | | | | | | | | \$3 |
|  | 500-236c | Front Castor Leg for "C" Series chairs | | | | | | | | | | | | Note that the front castor leg is not interchangeable with the rear castor leg. To distinguish between C and W legs, the push-pin on C is 150mm (6") up from the bottom | \$40 |
|  | 500-236w | Front Castor Leg for "W" Series chairs | | | | | | | | | | | | Note that the front castor leg is not interchangeable with the rear castor leg. To distinguish between W and C legs, the push-pin on W is 170mm (6.75") up from the bottom | \$40 |
|  | 500-010 | Rear Castor Leg for "C" Series chairs | | | | | | | | | | | | Note that the rear castor leg is not interchangeable with the front castor leg. For SB6w - Refer 500-035 | \$40 |
|  | 500-035 | SB6w Anti Tip Rear Leg | | | | | | | | | | | | Anti-tip rear leg for SB6w only | \$53 |
|  | 500-036 | SB6w Anti Tip Leg end cap | | | | | | | | | | | | Anti-tip end cap | \$9 |
|  | 500-259 | Leg Tube | | | | | | | | | | | | This original design footrest leg tube is fitted to SB6c chairs below serial #DC-SB6c-00629 and SB6w below DC-SB6w- 00647 | \$27 |
|  | 500-259-R Single Rear Slot | Leg Tube | | | | | | | | | | | | This upgraded design footrest leg tube is fitted to SB6c chairs from serial #DC-SB6c-00629 and SB6w from DC-SB6w- 00647 and to SB2T chairs from serial #211605 and SB3T chairs from #311666 BUT IS NOT interchangeable with chairs fitted with Part 500-259-C (double slot) | \$27 |
|  | 500-259-C Slot Front & Rear | Leg Tube | | | | | | | | | | | | New improved part with a slot front and rear for the locking pin to engage into. This part is not interchangeable with chairs initially fitted with part 500-259-R or part 500-259. Fitted to SB2T from DC212063 onwards, SB3T from DC312290 onwards, SB6c from 988 onwards and SB6w from 636-654 plus 680 onwards. | \$27 |
|  | 500-005 | Leg Plastic Insert | | | | | | | | | | | | | \$25 |
|  | 500-270 | Safety Restraints x 1 set. Can be used as either Chest or Lap (additional sets available as accessory) | | | | | | | | | | | | Fabric belts (webbing). Fixed to push handle for waist belt or underside of cushion rear tube for lap belt using 2 x part #500-108 (order separately) | \$117 |
|  | 118 | Safety Restraints (set) Chest and Lap | | | | | | | | | | | | New Webbing Fabric belts. Screw Fixed lap belt with chest belt being passed through loops in rear of fabric back rest | \$117 |
|  | 500-059 | SB6w Hand Brake | | | | | | | | | | | | Price PER PAIR. Fitted to all chairs to serial DC-SB6w-01518 (and also including DC-SB6w-01546 to 01595 inclusive) | \$84 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lite | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveler | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|-------------|--|------------------------------|----------|---------------|-------------------|-------------------|------------------|-------------------|-------------------|--------------|------------------------|------------------------|--|---------|
|  | 100-503 | SB6w Hand Brake | | | | | | | | | | | | Price PER PAIR. Fitted to SB6w chairs from serial DC-SB6w-01519 (but excluding DC-SB6w-01546 to 01595 inclusive) | \$84 |
|  | 500-054 | SB6w 24" wheel axle RH | | | | | | | | | | | | Axle for 24" wheel comes complete as kit with fixing bolts | \$54 |
|  | 500-053 | SB6w 24" wheel axle LH | | | | | | | | | | | | Axle for 24" wheel comes complete as kit with fixing bolts | \$54 |
|  | 500-053-B | Stainless Steel Bolt Secures 500-053/54 Axle Housing | | | | | | | | | | | | | \$3 |
|  | 500-275 | SB6w 24" Wheel | | | | | | | | | | | | | \$64 |
|  | 500-052 | SB6w 24" wheel axle only | | | | | | | | | | | | | \$23 |
|  | 500-149 | P1 RH Arm Assembly | | | | | | | | | | | | Complete arm that comes with arm pad and plastic caps | \$122 |
|  | 500-142 | P1 LH Arm Assembly | | | | | | | | | | | | Complete arm that comes with arm pad and plastic caps | \$122 |
|  | 500-152 | P1 Back Cushion | | | | | | | | | | | | Smallest Pediatric = P1 size | \$174 |
|  | 500-153 | P1 - P2 Seat Cushion support bracket | | | | | | | | | | | | | \$40 |
|  | 500-155 | P1 - P2 Pediatric seat support bracket recess panel | | | | | | | | | | | | | \$14 |
|  | 500-263 | P1 - P2 Seat lock Knob | | | | | | | | | | | | For adjustment of the seat back angle and position. | \$13 |
|  | 500-157 | P1 - P2 Pediatric seat Spine Collar | | | | | | | | | | | | | \$10 |
|  | 500-158 | P1 - P2 Pediatric seat spine end cap | | | | | | | | | | | | | \$10 |
|  | 500-161 | P2 back Cushion | | | | | | | | | | | | | \$194 |
|  | 500-163 | P1 - P2 M6 x 30 Bolt | | | | | | | | | | | | Attaches Seat Cushion Support bracket 500-153 to the PU Back Cushion 500-161 | \$3 |
|  | 500-169 | P2 LH Arm Assembly | | | | | | | | | | | | | \$122 |
|  | 500-173 | P2 RH Arm Assembly | | | | | | | | | | | | | \$122 |
|  | 500-178 | P1 Head Support | | | | | | | | | | | | | \$35 |
|  | 500-179 | P2 Head Support | | | | | | | | | | | | | \$40 |
|  | 500-480 | P1 - Pediatric Push Handle / Backrest Support Frame | | | | | | | | | | | | New Design from October 2017 | \$127 |
| 19/10/20 | Part Number | Part Description | This Part Fits Models Marked | | | | | | | | | | ✓ | Description | #VALUE! |
| Model Order Code | | | SB1 | SB2 | SB2T | SB3T | SB6c | SB6c-22 | SB6c-26 | SB6w | SB7e | P1 | P2 | | #VALUE! |

| Model Names | | | ShowerBuddy | TubBuddy | TubBuddy Tilt | Roll-InBuddy Tilt | Roll-InBuddy Lie | Roll-In Buddy XL | Roll-In Buddy XXL | Roll-InBuddy Solo | Eco-Traveller | P1 Pediatric Converter | P2 Pediatric Converter | | \$0 |
|---|---------|---|-------------|----------|---------------|-------------------|------------------|------------------|-------------------|-------------------|---------------|------------------------|------------------------|------------------------------|-------|
|  | 500-481 | P2 - Pediatric Push Handle / Backrest Support Frame | | | | | | | | | | | | New Design from October 2017 | \$127 |