



Supporting you

1300 362 923
www.healthtec.co

WALKING AND STANDING RANGE

Electric Parallel Bars

We pride ourselves on designing and manufacturing quality products that make our clients working lives easier. Our electric parallel bars are a perfect example of this.

The electric function enables the practitioner to adjust the bars to fit the patients required height. If for example the patient is in a wheel chair, the bars can be used to lift the patient into a standing position removing any OH&S issues and potential of injuries.

As the unit bolts straight to the floor there are no obstacles for the client to negotiate as they walk, allowing them to totally focus on their practitioners instructions.

Available in standard lengths of 4, 6 & 8 metres we can also do custom lengths if required. 4 metre set also available with width adjustment.

FEATURES

Durable powder coated rust resistant frame

No tripping hazards along walkway

Australian designed & manufactured

SPECIFICATIONS

LIFT CAPACITY	200kg (per side — moving)
HEIGHT RANGE	600mm to 1000mm with infinite adjustment through total range of movement
OPERATION	Battery operated



The Healthtec Warranty

Healthtec offers a lifetime warranty on all our manufactured tables and a five year warranty on motors and foot switches from date of purchase. Original suppliers warranty applies.

ARTG 132671

© August 2018





Supporting you

PRICE SHEET


www.healthtec.com.au

Electric Parallel Bars

Product Description	CODE	QTY
3M Electric Parallel Bars (600mm to 1000mm)	8032	
4M Electric Parallel Bars (600mm to 1000mm)	8038	
4M Electric Parallel Bars (600mm to 1000mm) w/Electric Width Adjustment	8041	
6M Electric Parallel Bars (600mm to 1000mm)	8040	
8M Electric Parallel Bars (600mm to 1000mm)	8039	
Universal Parallel Bar (single side) — 3M Electric Hi-Lo (600mm to 1000mm)	8031	
Universal Parallel Bar (single side) — 4M Electric Hi-Lo (600mm to 1000mm)	8035	
Universal Parallel Bar (single side) — 6M Electric Hi-Lo (600mm to 1000mm)	8034	
Universal Parallel Bar (single side) — 8M Electric Hi-Lo (600mm to 1000mm)	8033	
Universal Parallel Bar (single side) — 4M Electric Hi-Lo Electric Width Adjustment (600mm to 1000mm)	8044	
4M Manual Height Parallel Bars (600mm to 1100mm)	8036	

Freight

Freight charges are not included in our prices.
Please contact your Healthtec representative for freight costs and options.



FREECALL
1300 362 923
info@healthtec.com.au



Healthtec Pty Ltd ABN 94 224 518 379 ARTG 132671
© August 2018

Disclaimer: Due to our policy of constant product improvement,
specifications may change without notice.
Check with us for the latest specifications and pricing.