

Bed Pressure-release | Floor Stepping Fall Prevention Sensor Pad

The SBM-8000 is a Bed Pressure-release | Floor Stepping Fall Prevention Sensor Pad mat that has pressure sensors which are activated when pressure is released or applied, depending on usage preference. When used in conjunction with STS-100FM, a wireless transmitter, an alarm alert can be sent out.

It is ideal for caring for patients in hospitals or at home as it allows you to be aware when the person gets out of bed and thus offer immediate assistance.

In the same way, it is an excellent device for security uses since it can go completely unnoticed by possible intruders, probably pending other more conventional alarm systems and not stepping on the carpet that gave them away.

Features and Specification

- Frequency type: FSK (accuracy and reliability)
- Battery: 2 x CR2032
- CPU TI: which has its own communication protocol
- Button ID: 16 million ID's available
- Call indicator: red LED on the top
- 2 colors available: Black / White
- Size: 100x60x6.5mm
- Weight: 50G
- Appropriate for: Hospitals, Home and various uses in security aspect
- Can be paired with STS-100FM, a wireless transmitter