



A Better Wheelchair, A Better Fit

**KARMA MEDICAL PRODUCTS CO., LTD.**

康揚股份有限公司

T +886 5 206 6688

F +886 5 206 7788

www.karmamedical.com / E globalsales@karma.com.tw

BRANCHES: U. K. / TAIWAN / GERMANY / THAILAND / CHINA / INDIA / U.S.A.

## KARMA KM-BT10 Bariatric wheelchair Order Form

P.O.# \_\_\_\_\_ Date: \_\_\_\_\_

Karma Agent: \_\_\_\_\_

### Product Information

- The KM-BT10 was designed for user with a 550lbs (250kgs) limit.
- All measurements are without cushion.

### Frame Finish

- Diamond Black

STD

### Seat Width & Depth (Selection 1 or 2)

#### Selection 1 :

STD S\$ 4350/set Seat

Width can grow from 22" to 26"

Seat Depth can grow from 18" to 20"

set at:

Seat Width \ Seat Depth	22"	24"	26"
18"	●	●	●
20"	●	●	●

- ☐ 22" Seat Width x 18" Seat Depth \_\_\_\_\_ set
- ☐ 22" Seat Width x 20" Seat Depth \_\_\_\_\_ set
- ☐ 24" Seat Width x 18" Seat Depth \_\_\_\_\_ set
- ☐ 24" Seat Width x 20" Seat Depth \_\_\_\_\_ set
- ☐ 26" Seat Width x 18" Seat Depth \_\_\_\_\_ set
- ☐ 26" Seat Width x 20" Seat Depth \_\_\_\_\_ set

#### Selection 2 :

+S\$ 550

Seat Width can grow from 26" to 30"

Seat Depth can grow from 22" to 24"

set at:

Seat Width \ Seat Depth	26"	28"	30"
22"	●	●	●
24"	●	●	●

- ☐ 26" Seat Width x 22" Seat Depth \_\_\_\_\_ set
- ☐ 26" Seat Width x 24" Seat Depth \_\_\_\_\_ set
- ☐ 28" Seat Width x 22" Seat Depth \_\_\_\_\_ set
- ☐ 28" Seat Width x 24" Seat Depth \_\_\_\_\_ set
- ☐ 30" Seat Width x 22" Seat Depth \_\_\_\_\_ set
- ☐ 30" Seat Width x 24" Seat Depth \_\_\_\_\_ set

## Seat to Floor Height - Front

Selection 1: 17" \_\_\_\_, 18" \_\_\_\_, 19" \_\_\_\_ STD  
Based on 7" \*1 3/4 castor and 24\*2.125 rear wheels


Selection 2: 17.5" \_\_\_\_, 18.5" \_\_\_\_, 19.5" \_\_\_\_ +S\$ 180  
Based on 200\*50 castor and 24\*2.125 rear wheels

## Backrest Height

☐ Selection 1: 16" ~ 18" Adjustable Backrest Height, set at: NC  
16" \_\_\_\_, 17" \_\_\_\_, 18" \_\_\_\_

☐ Selection 2: 18" ~ 20" Adjustable Backrest Height, set at: NC  
18" \_\_\_\_, 19" \_\_\_\_, 20" \_\_\_\_

## Center-of-Gravity Position



Seat Depth	Backrest Adjustable Plate-Outside	Rear Wheel Bracket Plate			
		Active Position		Stable Position	
	Your Selection	Please tick the rear wheel position correspondent to the seat depth			
18"	① . ③ . ⑤	【A B】 <input type="checkbox"/>	【B C】 <input type="checkbox"/>	【D E】 <input type="checkbox"/>	【E F】 <input type="checkbox"/>
20"	③ . ⑤ . ⑦	【D E】 <input type="checkbox"/>	【E F】 <input type="checkbox"/>	【F G】 <input type="checkbox"/>	【G H】 <input type="checkbox"/>
22"	⑤ . ⑦ . ⑨	【F G】 <input type="checkbox"/>	【G H】 <input type="checkbox"/>	【H I】 <input type="checkbox"/>	
24"	⑦ . ⑨ . ⑪	【H I】 <input type="checkbox"/>		【H I】 <input type="checkbox"/> + Rear wheel extension plate L <input type="checkbox"/> or M <input type="checkbox"/>	

## Back Angle

−5° ~ 15° Back Angle Adjustment, set at:

NC

−5° \_\_\_\_\_ −0° \_\_\_\_\_ 5° \_\_\_\_\_ 10° \_\_\_\_\_ 15° \_\_\_\_\_

## Front Castor

☐ Front Casters are hard solid tires (7" \*1 3/4)

STD

☐ Front Casters are hard solid tires (200\*50)

NC

## Brake

☒ Parking Brake

STD

☒ Attendant (Wire) Brake

STD

## Rear Wheels

☐ Rear Wheels Pneumatic tires 24" \*2.125

STD

## Rear Axle

☐ Fixed

STD

☐ Quick Release

+S\$ 120

## Handrim

☐ Aluminum Handrims

STD

## Space Between The Handrim and Rear Wheel

☐ 15mm

STD

☐ 5mm

+S\$ 90

## Armrest Type

☐ Flip-Back & Detachable armrest

STD

☐ Flip-Back & Detachable armrest with Height Adjustable

+S\$ 105

## Footrest

☐ Swing-away footrest

STD

☐ Elevating and Swing away footrest

+S\$ 240

## Footplate

3D Adjustable Aluminum Footplate

STD

## Footrest Accessories

☐ Calf Strap

STD

☐ Adjustable Length Heel Loop

+S\$ 36

## Seat Positioning Belt

■ Hip Belt with Full Length 1570mm x 2

STD

## Seat Accessories

☐ 3cm Foam Cushion

STD

## Back Accessories

■ Foldable Push Bar

STD

## Miscellaneous Accessories

■ Anti-tipper

STD

☐ IV Pole

+S\$ 120

☐ Oxygen Tank Holder

+S\$ 130

(need DIY assemble the accessories)