

LATERAL TURNING BED

Motorized

KS-803



Manual

KS-1032



Description : 4-functions controlled electronically – lateral turning (30 degrees left & right), head-raising (0-75 degrees), leg-raising (0-40 degrees), bed-height adjustment (55-85 cm, inclusive of standard 10cm mattress). Steel bed frame. Mesh bed base. Aluminium alloy bedrail. 5" castors.

Dimension : 212X970X55-85cm(inclusive of mattress)

Kangshen Medical Equipment Co. Ltd is a Chinese-German joint venture located in Foshan City, China. It has obtained ISO13485-2003 Quality System certification.

Rehab Mart