

Oxford® Classic Mobile Lifts

Spares

Effective 12th February 2018

Oxford® Mini 140

PART NUMBER	ITEM DESCRIPTION	LIST PRICE
39001256	CLASSIC FRONT CASTORS (PAIR - INC FIXINGS) (75mm castor x2, M8 Nyloc nut x2, M8 csk screw x2 & 19mm screw cover x2)	\$113
39001257	CLASSIC REAR CASTORS (PAIR - INC FIXINGS) (100mm castor x2, M10 Nyloc nut x2, M10 csk screw x2 & 25mm screw cover x2)	\$112
31009213	MINI 140 LEG	\$219
39001192	MINI140/MIDI180 CROSS MEMBER INC FIXINGS (M12 csk screw x2, M6 hex screw x1, Anti-reverse bar x1, Brass bearing washer x4, leg compression sleeve x2, 25mm screw cover x2, Acetal wear washer x1, M16 Binx nut x1, Tie rod sleeve x2, M12 nylon washer x4 & M10 bttm screw x2)	\$458
39001191	MINI 140/MIDI 180 TIE RODS WITH FIXINGS (PAIR) (Tie rod LH x1, tie rod RH x1, M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4)	\$88
39001188	MINI 140/MIDI 180 FOOT PEDAL ASSEMBLY (Pedal assembly x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4)	\$394
39001393	LEG MECHANISM ASSEMBLY (Leg mechanism x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4) <i>Electric leg models only</i>	\$183
39001189	CLASSIC ADHESIVE FOOT PEDAL PADS (NOT SHOWN)	\$82
39001190	CLASSIC FOOT PEDAL FASTENING KIT (M16 binx nut x1, Acetal wear washer x1, M12 Nylon washer x8, Mast locking knob x1, M10 button hd screw x2, M8 button hd screw x2, Foot pedal cover x1, Tie rod sleeve x4) (NOT SHOWN)	\$67
39001181	CLASSIC FOOT PEDAL PIVOT COVER (NOT SHOWN)	\$141
39001260	MINI 140 MAST (INC FIXINGS) [Mini mast fabrication x1, Fulcrum pin x1, 10mm external circlip x2, No8 pozi pan/ tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2 & Mast decals]	\$405
39001402	MINI 140 MAST (INC FIXINGS) [Mini mast fabrication x1, Fulcrum pin x1, 10mm external circlip x2, No8 pozi pan/ tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2, actuator bracket & Mast decals] <i>Electric leg models only</i>	\$482
39001258	MINI 140 BOOM (INC FIXINGS) (Mini Boom fabrication x1, Boom retainer washer x2, Boom fulcrum pin end cap x2, Boom fulcrum bush x1, Boom fulcrum pin x1, M8 csk screw x2, 19mm screw cover x2, Boom end cap x1, plastic cap x2, 10mm external circlip x4, Fulcrum pin x1, Cap retaining washer x2 & Boom decals)	\$277
39001200	MINI 140/MIDI 180 TOP FIXING KIT (Boom retainer washer x2, Boom fulcrum bush x1, Boom fulcrum pin x1, M8 csk screw x2, Boom fulcrum pin end cap x2, Fulcrum pin x3, cap retaining washer x2, 10mm circlip x6, Plastic cap x2, Push handle end stop x2, 19mm screw cover x2, Boom end cap x1)	\$84
39001199	MINI 140/MIDI 180 BASE FIXING KIT (M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4, M16 binx nut x1, Acetal wear washer x1, Mast locking knob x1, Foot pedal cover x1, M10 csk screw x2, 19mm screw cover x2, 25mm screw cover x2, M10 Nyloc nut x2, M8 Nyloc nut x2, Brass leg bearing washer x4, Leg compression sleeve x2, M12 csk screw x2, M10 Nyloc nut x2)	\$113
39001273	MINI 140 HYDRAULIC PUMP (INC FIXINGS) (Mini Hydraulic Pump x1, Fulcrum pin x2 & 10mm external circlip x4)	\$1,721
39001270	MINI 140 ACTUATOR LA31 (INC FIXINGS) (LA31 300mm stroke x1, Fulcrum pin x2, 10mm external circlips x4)	\$474

39001396	LEG ACTUATOR (INC FIXING) (Set pin, nylon spacer, nylon washer, spring washer) Electric leg models only	\$408
39001398	LOWER ACTUATOR FIXING KIT Electric leg models only	\$229
39001274	SMART MONITOR RETRO-FIT KIT (Smart monitor x1, Smart monitor handset (1 channel only) x1, Battery pack x 1, M5 hex cap screw x1)	\$1,122
39001167	SMART MONITOR CONTROL BOX - 1 CHANNEL	\$711
39001395	SMART MONITOR CONTROL BOX - 2 CHANNEL Electric leg models only	\$745
39001180	SMART MONITOR USER HANDSET - 1 CHANNEL	\$184
39001394	SMART MONITOR USER HANDSET - 2 CHANNEL Electric leg models only	\$179
39001178	SMART MONITOR ACTUATOR RESET HANDSET (NOT SHOWN)	\$343
39001179	SMART MONITOR CALIBRATION HANDSET (NOT SHOWN)	\$237
39001297	SERVICE PERIOD RE-SET HANDSET (SMART™ MONITOR ONLY) (NOT SHOWN)	\$340
11012180	BATTERY PACK 24V RECHARGEABLE *	\$471
39001201	CLASSIC MOUNTING BRACKET (INC FIXINGS) (Mounting bracket x1, No8 * Pozi pans/tap screw x3)	\$59
39001183	CLASSIC CHARGER & STAND ASSEMBLY WITH UK LEAD	\$431
11012183	CHARGER LEAD - UK (NOT SHOWN)	\$141
11012184	CHARGER LEAD - EU (NOT SHOWN)	\$123
11013804	CHARGER LEAD - AU (NOT SHOWN)	\$157
31009237	CLASSIC SMALL SPREADER BAR *	\$151
31009196	CLASSIC STANDARD SPREADER BAR *	\$157
31009197	CLASSIC LARGE SPREADER BAR *	\$170
11014480	CLASSIC SMALL SPREADER BAR COVER (NOT SHOWN)	\$72
11014250	CLASSIC STANDARD SPREADER BAR COVER (NOT SHOWN)	\$76
11014263	CLASSIC LARGE SPREADER BAR COVER (NOT SHOWN)	\$91
39001171	MINI 140/MIDI 180 SPREADER BAR FIXING KIT (Midi/Mini Boom end cap x1, Fulcrum pin x1, 10mm circlips x2, Cap retaining washer x2, Plastic cap x2, Hook shaft washer x5 & O-ring x5) (NOT SHOWN)	\$54
39001304	CLASSIC SPREADER BAR SLING RETAINING WASHER KIT (Retaining washer x2, m6 screw x2, spring x2, spacing post x2)	\$46
39001279	CLASSIC SPREADER BAR FULCRUM PIN CAP (NOT SHOWN)	\$22
11014467	MINI 140 ELECTRIC DECAL SET (NOT SHOWN)	\$33
11014482	MINI 140 HYDRAULIC DECAL SET (NOT SHOWN)	\$33
39001198	CLASSIC CIRCLIPS (NOT SHOWN)	\$47
39001182	CLASSIC PUSH HANDLE END CAPS	\$29
39001276	MINI 140/MIDI 180 BOOM END CAPS	\$143
39001277	CLASSIC 19MM GREY SCREW COVERS (NOT SHOWN)	\$30
39001278	CLASSIC 25MM GREY SCREW COVERS (NOT SHOWN)	\$27
39001397	RECTANGULAR TUBE INSERT FOR LEG ACTUATOR BRACKET - WHITE (NOT SHOWN) Electric leg models only	\$123
39001399	19MM UNI CAP - WHITE (NOT SHOWN) Electric leg models only	\$21

Oxford® Midi 180

PART NUMBER	ITEM DESCRIPTION	#VALUE!
39001256	CLASSIC FRONT CASTORS (PAIR - INC FIXINGS) (75mm castor x2, M8 Nyloc nut x2, M8 csks screw x2 & 19mm screw cover x2)	\$113
39001257	CLASSIC REAR CASTORS (PAIR - INC FIXINGS) (100mm castor x2, M10 Nyloc nut x2, M10 csk screw x2 & 25mm screw cover x2)	\$112
31009203	MIDI 180 LEG	\$230
39001192	MINI140/MIDI180 CROSS MEMBER INC FIXINGS (M12 csk screw x2, M6 hex screw x1, Anti-reverse bar x1, Brass bearing washer x4, leg compression sleeve x2, 25mm screw cover x2, Acetal wear washer x1, M16 Binx nut x1, Tie rod sleeve x2, M12 nylon washer x4 & M10 bttm screw x2)	\$458
39001191	MINI 140/MIDI 180 TIE RODS WITH FIXINGS (PAIR) (Tie rod LH x1, tie rod RH x1, M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4)	\$88

39001188	MINI 140/MIDI 180 FOOT PEDAL ASSEMBLY (Foot pedal assembly x1, M10 Button head screw x2, Tie rod sleeve x2, Nylon washer x4)	\$394
39001393	LEG MECHANISM ASSEMBLY (Leg mechanism x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4)	\$183
	Electric leg models only	
39001189	CLASSIC ADHESIVE FOOT PEDAL PADS (NOT SHOWN)	\$82
39001190	CLASSIC FOOT PEDAL FASTENING KIT (M16 binx nut x1, Acetal wear washer x1, M12 Nylon washer x8, Mast locking knob x1, M10 button hd screw x2, M8 button hd screw x2, Foot pedal cover x1, Tie rod sleeve x4) (NOT SHOWN)	\$67
39001181	CLASSIC FOOT PEDAL PIVOT COVER (NOT SHOWN)	\$141
39001176	MIDI 180 MAST (INC FIXINGS) (Midi mast fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, Push handle end cap x2, Boom retainer washer x2, Boom fulcrum bush x1, Mounting bracket x1, Boom fulcrum pin end cap x2, Boom fulcrum pin x1, 19mm screw cover x2 & Mast decals)	\$405
39001401	MIDI 180 MAST (INC FIXINGS) (Midi mast fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, Push handle end cap x2, Boom retainer washer x2, Boom fulcrum bush x1, Mounting bracket x1, Boom fulcrum pin end cap x2, Boom fulcrum pin x1, 19mm screw cover x2, actuator bracket & Mast decals) Electric leg models only	\$471
39001175	MIDI 180 BOOM (INC FIXINGS) (Midi boom fabrication x1, boom retainer washer x2, boom fulcrum pin end cap x2, boom fulcrum bush x1, boom fulcrum pin x1, M8 csk screw x2, 19mm screw cover x2, boom end cap x1, Plastic cap x2, 10mm external circlip x4, Fulcrum pin x1, Cap retaining washer x2 & boom decals)	\$276
39001200	MINI 140/MIDI 180 TOP FIXING KIT (Boom retainer washer x2, Boom fulcrum bush x1, Boom fulcrum pin x1, M8 csk screw x2, Boom fulcrum pin end cap x2, Fulcrum pin x3, cap retaining washer x2, 10mm circlip x6, Plastic cap x2, Push handle end stop x2, 19mm screw cover x2, Boom end cap x1)	\$84
39001199	MINI 140/MIDI 180 BASE FIXING KIT (M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4, M16 binx nut x1, Acetal wear washer x1, Mast locking knob x1, Foot pedal cover x1, M10 csk screw x2, 19mm screw cover x2, 25mm screw cover x2, M10 Nyloc nut x2, M8 Nyloc nut x2, Brass leg bearing washer x4, Leg compression sleeve x2, M12 csk screw x2, M10 Nyloc nut x2)	\$113
39001272	MIDI 180 HYDRAULIC PUMP (INC FIXINGS) (Midi Hydraulic Pump x1, Fulcrum pin x2 & 10mm external circlip x4)	\$1,771
39001177	MIDI 180 ACTUATOR LA44 (INC FIXINGS) (LA44 250mm stroke x1, Fulcrum pin x2, 10mm external circlip x4, LA44 Lead x1, emergency descent decal x1)	\$941
11013922	ACTUATOR POWER LEAD (LA44) (NOT SHOWN)	\$28
39001396	LEG ACTUATOR (INC FIXING) (Set pin, nylon spacer, nylon washer, spring washer) Electric leg models only	\$408
39001398	LOWER ACTUATOR FIXING KIT Electric leg models only	\$229
39001274	SMART MONITOR RETROFIT KIT (Smart monitor x1, Smart monitor handset (1 channel only) x1, Battery pack x1, M5 hex cap screw x1)	\$1,122
39001167	SMART MONITOR CONTROL BOX - 1 CHANNEL	\$711
39001395	SMART MONITOR CONTROL BOX - 2 CHANNEL Electric leg models only	\$745
39001180	SMART MONITOR USER HANDSET - 1 CHANNEL	\$184
39001394	SMART MONITOR USER HANDSET - 2 CHANNEL Electric leg models only	\$179
39001178	SMART MONITOR ACTUATOR RESET HANDSET (NOT SHOWN)	\$343
39001179	SMART MONITOR CALIBRATION HANDSET (NOT SHOWN)	\$237
39001297	SERVICE PERIOD RE-SET HANDSET (SMART™ MONITOR ONLY) (NOT SHOWN)	\$340
11012180	BATTERY PACK 24V RECHARGEABLE *	\$471
39001201	CLASSIC MOUNTING BRACKET (INC FIXINGS) (Mounting bracket x1, No8 * Pozi pans/tap screw x3)	\$59
39001183	CLASSIC CHARGER & STAND ASSEMBLY WITH UK LEAD	\$431
11012183	CHARGER LEAD - UK (NOT SHOWN)	\$141
11012184	CHARGER LEAD - EU (NOT SHOWN)	\$123
11013804	CHARGER LEAD - AU (NOT SHOWN)	\$157
31009237	CLASSIC SMALL SPREADER BAR *	\$151
31009196	CLASSIC STANDARD SPREADER BAR *	\$157
31009197	CLASSIC LARGE SPREADER BAR *	\$170

11014480	CLASSIC SMALL SPREADER BAR COVER (NOT SHOWN)	\$72
11014250	CLASSIC STANDARD SPREADER BAR COVER (NOT SHOWN)	\$76
11014263	CLASSIC LARGE SPREADER BAR COVER (NOT SHOWN)	\$91
39001171	MINI 140/MIDI 180 SPREADER BAR FIXING KIT (Midi/Mini Boom end cap x1, Fulcrum pin x1, 10mm circlips x2, Cap retaining washer x2, Plastic cap x2, Hook shaft washer x5 & O-ring x5) (NOT SHOWN)	\$54
39001304	CLASSIC SPREADER BAR SLING RETAINING WASHER KIT (Retaining washer x2, m6 screw x2, spring x2, spacing post x2)	\$46
39001279	CLASSIC SPREADER BAR FULCRUM PIN CAP (NOT SHOWN)	\$22
11014468	MIDI 180 ELECTRIC DECAL SET (NOT SHOWN)	\$33
11014483	MIDI 180 HYDRAULIC DECAL SET (NOT SHOWN)	\$38
39001198	CLASSIC CIRCLIPS (NOT SHOWN)	\$47
39001182	CLASSIC PUSH HANDLE END CAPS	\$29
39001276	MINI 140/MIDI 180 BOOM END CAP	\$143
39001277	CLASSIC 19MM GREY SCREW COVERS (NOT SHOWN)	\$30
39001278	CLASSIC 25MM GREY SCREW COVERS (NOT SHOWN)	\$27
39001397	RECTANGULAR TUBE INSERT FOR LEG ACTUATOR BRACKET - WHITE (NOT SHOWN) <i>Electric leg models only</i>	\$123
39001399	19MM UNI CAP - WHITE (NOT SHOWN) <i>Electric leg models only</i>	\$21

Oxford® Major 200

PART NUMBER	ITEM DESCRIPTION	#VALUE!
39001256	CLASSIC FRONT CASTORS (PAIR - INC FIXINGS) (75mm castor x2, M8 Nyloc nut x2, M8 csks screw x2 & 19mm screw cover x2)	\$113
39001257	CLASSIC REAR CASTORS (PAIR - INC FIXINGS) (100mm castor x2, M10 Nyloc nut x2, M10 csk screw x2 & 25mm screw cover x2)	\$112
31009208	MAJOR 200/MAXI 180 LEG	\$240
39001265	MAJOR200/MAXI180 CROSSMEMBER INC FIXINGS (M12 csk screw x2, M6 hex screw x1, Anti-reverse bar x1, Brass bearing washer x4, leg compression sleeve x2, 25mm screw cover x2, Acetal wear washer x1, M16 Binx nut x1, Tie rod sleeve x2, M12 nylon washer x4 & M10 bttm screw x2)	\$482
39001264	MAJOR 200/MAXI 180 TIE RODS WITH FIXINGS (PAIR) (Tie rod LH x1, tie rod RH x1, M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4)	\$91
39001275	MAJOR 200/MAXI 180 FOOT PEDAL ASSEMBLY (Foot pedal assembly x1, M10 Button head screw x2, Tie rod sleeve x2, Nylon washer x4)	\$402
39001405	LEG MECHANISM ASSEMBLY (Leg mechanism x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4) <i>Electric leg models only</i>	\$186
39001189	CLASSIC ADHESIVE FOOT PEDAL PADS (NOT SHOWN)	\$82
39001190	CLASSIC FOOT PEDAL FASTENING KIT (M16 binx nut x1, Acetal wear washer x1, M12 Nylon washer x8, Mast locking knob x1, M10 button hd screw x2, M8 button hd screw x2, Foot pedal cover x1, Tie rod sleeve x4) (NOT SHOWN)	\$67
39001181	CLASSIC FOOT PEDAL PIVOT COVER (NOT SHOWN)	\$141
39001262	MAJOR 200 MAST (INC FIXINGS) (Major Mast Fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2 & Mast decals)	\$465
39001403	MAJOR 200 MAST (INC FIXINGS) (Major Mast Fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2, actuator bracket & Mast decals) <i>Electric leg models only</i>	\$542
39001259	MAJOR 200/MAXI 180 BOOM (INC FIXINGS) (Major/Maxi Boom fabrication x1, Boom retainer washer x2, Boom fulcrum pin end cap x2, Boom fulcrum bush x1, boom fulcrum pin x1, M8 csk screw x2, 19mm screw cover x2, Boom end cap x1, plastic cap x2, 10mm external circlip x4, Fulcrum pin x1, Cap retaining washer x2 & Boom decals)	\$320

39001268	MAJOR 200/MAXI 180 TOP FIXING KIT (Boom retainer washer x2, Boom fulcrum bush x1, Boom fulcrum pin x1, M8 csk screw x2, Boom fulcrum pin end cap x2, Fulcrum pin x3, cap retaining washer x2, 10mm circlip x6, Plastic cap x2, Push handle end stop x2, 19mm screw cover x2, Boom end cap x1)	\$114
39001269	MAJOR 200/MAXI 180 BASE FIXING KIT (M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4, M16 binx nut x1, Acetal wear washer x1, Mast locking knob x1, Foot pedal cover x1, M10 csk screw x2, 19mm screw cover x2, 25mm screw cover x2, M10 Nyloc nut x2, M8 Nyloc nut x2, Brass leg bearing washer x4, Leg compression sleeve x2, M12 csk screw x2, M10 Nyloc nut x2)	\$113
39001271	MAJOR 200 HYDRAULIC PUMP (INC FIXINGS) (Major Hydraulic Pump x1, Fulcrum pin x2 & 10mm external circlip x4)	\$1,766
39001263	MAJOR200/MAXI180 ACTUATOR LA44 (INC FIXINGS) (LA44 300mm stroke x1, Fulcrum pin x2, 10mm external circlip x4, LA44 Lead x1, emergency descent decals x1)	\$970
11013922	ACTUATOR POWER LEAD (LA44) (NOT SHOWN)	\$28
39001396	LEG ACTUATOR (INC FIXING) (Set pin, nylon spacer, nylon washer, spring washer) Electric leg models only	\$408
39001398	LOWER ACTUATOR FIXING KIT Electric leg models only	\$229
39001274	SMART MONITOR RETROFIT KIT (Smart monitor x1, Smart monitor handset (1 channel only) x1, Battery pack x1, M5 hex cap screw x1)	\$1,122
39001167	SMART MONITOR CONTROL BOX - 1 CHANNEL	\$711
39001395	SMART MONITOR CONTROL BOX - 2 CHANNEL Electric leg models only	\$745
39001180	SMART MONITOR USER HANDSET - 1 CHANNEL	\$184
39001394	SMART MONITOR USER HANDSET - 2 CHANNEL Electric leg models only	\$179
39001178	SMART MONITOR ACTUATOR RESET HANDSET (NOT SHOWN)	\$343
39001179	SMART MONITOR CALIBRATION HANDSET (NOT SHOWN)	\$237
39001297	SERVICE PERIOD RE-SET HANDSET (SMART™ MONITOR ONLY) (NOT SHOWN)	\$340
	BATTERY without pack	\$78
11012180	BATTERY PACK 24V RECHARGEABLE *	\$471
39001201	CLASSIC MOUNTING BRACKET (INC FIXINGS) (Mounting bracket x1, No8 * Pozi pans/tap screw x3)	\$59
39001183	CLASSIC CHARGER & STAND ASSEMBLY WITH UK LEAD	\$431
11012183	CHARGER LEAD - UK (NOT SHOWN)	\$141
11012184	CHARGER LEAD - EU (NOT SHOWN)	\$123
11013804	CHARGER LEAD - AU (NOT SHOWN)	\$157
31009237	CLASSIC SMALL SPREADER BAR *	\$151
31009196	CLASSIC STANDARD SPREADER BAR *	\$157
31009197	CLASSIC LARGE SPREADER BAR *	\$170
11014480	CLASSIC SMALL SPREADER BAR COVER (NOT SHOWN)	\$72
11014250	CLASSIC STANDARD SPREADER BAR COVER (NOT SHOWN)	\$76
11014263	CLASSIC LARGE SPREADER BAR COVER (NOT SHOWN)	\$91
39001172	MAJOR 200/MAXI 180 S / BAR FIXING KIT (Major/Maxi Boom end cap x1, Fulcrum pin x1, 10mm circlips x2, Cap retaining washer x2, Plastic cap x2, Hook shaft washer x5 & O-ring x5) (NOT SHOWN)	\$60
39001304	CLASSIC SPREADER BAR SLING RETAINING WASHER KIT (Retaining washer x2, m6 screw x2, spring x2, spacing post x2)	\$46
39001279	CLASSIC SPREADER BAR FULCRUM PIN CAP (NOT SHOWN)	\$22
11014470	MAJOR 200 ELECTRIC DECAL SET (NOT SHOWN)	\$33
11014484	MAJOR 200 HYDRAULIC DECAL SET (NOT SHOWN)	\$33
39001198	CLASSIC CIRCLIPS (NOT SHOWN)	\$47
39001182	CLASSIC PUSH HANDLE END CAPS	\$29
39001280	MAJOR 200/MAXI 180 BOOM END CAP	\$143
39001277	CLASSIC 19MM GREY SCREW COVERS (NOT SHOWN)	\$30
39001278	CLASSIC 25MM GREY SCREW COVERS (NOT SHOWN)	\$27
39001397	RECTANGULAR TUBE INSERT FOR LEG ACTUATOR BRACKET - WHITE (NOT SHOWN) Electric leg models only	\$123
39001399	19MM UNI CAP - WHITE (NOT SHOWN) Electric leg models only	\$21
		\$0
		\$0

PART NUMBER	ITEM DESCRIPTION	#VALUE!
39001256	CLASSIC FRONT CASTORS (PAIR - INC FIXINGS) (75mm castor x2, M8 Nyloc nut x2, M8 csks screw x2 & 19mm screw cover x2)	\$113
39001257	CLASSIC REAR CASTORS (PAIR - INC FIXINGS) (100mm castor x2, M10 Nyloc nut x2, M10 csk screw x2 & 25mm screw cover x2)	\$112
31009208	MAJOR 200/MAXI 180 LEG	\$240
39001265	MAJOR200/MAXI180 CROSSMEMBER INC FIXINGS (M12 csk screw x2, M6 hex screw x1, Anti-reverse bar x1, Brass bearing washer x4, leg compression sleeve x2, 25mm screw cover x2, Acetal wear washer x1, M16 Binx nut x1, Tie rod sleeve x2, M12 nylon washer x4 & M10 bttn screw x2)	\$482
39001264	MAJOR 200/MAXI 180 TIE RODS WITH FIXINGS (PAIR) (Tie rod LH x1, tie rod RH x1, M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4)	\$91
39001275	MAJOR 200/MAXI 180 FOOT PEDAL ASSEMBLY (Foot pedal assembly x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4)	\$402
39001405	LEG MECHANISM ASSEMBLY (Leg mechanism x1, M10 Button head screw x2, Tie rod sleeve x2, nylon washer x4) <i>Electric leg models only</i>	\$186
39001189	CLASSIC ADHESIVE FOOT PEDAL PADS (NOT SHOWN)	\$82
39001190	CLASSIC FOOT PEDAL FASTENING KIT (M16 binx nut x1, Acetal wear washer x1, M12 Nylon washer x8, Mast locking knob x1, M10 button hd screw x2, M8 button hd screw x2, Foot pedal cover x1, Tie rod sleeve x4) (NOT SHOWN)	\$67
39001181	CLASSIC FOOT PEDAL PIVOT COVER (NOT SHOWN)	\$141
39001261	MAXI 180 MAST (INC FIXINGS) (Maxi Mast Fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2 & Mast decals)	\$465
39001404	MAXI 180 MAST (INC FIXINGS) (Maxi Mast Fabrication x1, Fulcrum pin x1, 10mm external circlip x4, No8 pozi pan/tap screw x3, M8 csk screw x2, push handle end cap x2, boom retainer washer x2, boom fulcrum bush x1, Mounting bracket x1, boom fulcrum pin end cap x2, boom fulcrum pin x1, 19mm screw cover x2, actuator bracket & Mast decals) <i>Electric leg models only</i>	\$531
39001259	MAJOR 200/MAXI 180 BOOM (INC FIXINGS) (Major/Maxi Boom fabrication x1, Boom retainer washer x2, Boom fulcrum pin end cap x2, Boom fulcrum bush x1, boom fulcrum pin x1, M8 csk screw x2, 19mm screw cover x2, Boom end cap x1, plastic cap x2, 10mm external circlip x4, Fulcrum pin x1, Cap retaining washer x2 & Boom decals)	\$320
39001268	MAJOR 200/MAXI 180 TOP FIXING KIT (Boom retainer washer x2, Boom fulcrum bush x1, Boom fulcrum pin x1, M8 csk screw x2, Boom fulcrum pin end cap x2, Fulcrum pin x3, cap retaining washer x2, 10mm circlip x6, Plastic cap x2, Push handle end stop x2, 19mm screw cover x2, Boom end cap x1)	\$114
39001269	MAJOR 200/MAXI 180 BASE FIXING KIT (M10 button hd screw x2, M8 button HD screw x2, M12 Nylon washer x8, Tie rod sleeve x4, M16 binx nut x1, Acetal wear washer x1, Mast locking knob x1, Foot pedal cover x1, M10 csk screw x2, 19mm screw cover x2, 25mm screw cover x2, M10 Nyloc nut x2, M8 Nyloc nut x2, Brass leg bearing washer x4, Leg compression sleeve x2, M12 csk screw x2, M10 Nyloc nut x2)	\$113
39001263	MAJOR200/MAXI180 ACTUATOR LA44 (INC FIXINGS) (LA44 300mm stroke x1, Fulcrum pin x2, 10mm external circlip x4, LA44 Lead x1, emergency descent decals x1)	\$970
11013922	ACTUATOR POWER LEAD (LA44) (NOT SHOWN)	\$28
39001396	LEG ACTUATOR (INC FIXING) (Set pin, nylon spacer, nylon washer, spring washer) <i>Electric leg models only</i>	\$408
39001398	LOWER ACTUATOR FIXING KIT <i>Electric leg models only</i>	\$229
39001274	SMART MONITOR RETROFIT KIT (Smart monitor x1, Smart monitor handset (1 channel only) x1, Battery pack x1, M5 hex cap screw x1)	\$1,122
39001167	SMART MONITOR CONTROL BOX - 1 CHANNEL	\$711
39001395	SMART MONITOR CONTROL BOX - 2 CHANNEL <i>Electric leg models only</i>	\$745
39001180	SMART MONITOR USER HANDSET - 1 CHANNEL	\$184
39001394	SMART MONITOR USER HANDSET - 2 CHANNEL <i>Electric leg models only</i>	\$179
39001178	SMART MONITOR ACTUATOR RESET HANDSET (NOT SHOWN)	\$343

39001179	SMART MONITOR CALIBRATION HANDSET (NOT SHOWN)	\$237
39001297	SERVICE PERIOD RE-SET HANDSET (SMART™ MONITOR ONLY) (NOT SHOWN)	\$340
11012180	BATTERY PACK 24V RECHARGEABLE *	\$471
39001201	CLASSIC MOUNTING BRACKET (INC FIXINGS) (Mounting bracket x1, No8 * Pozi pans/tap screw x3)	\$59
39001183	CLASSIC CHARGER & STAND ASSEMBLY WITH UK LEAD	\$431
11012183	CHARGER LEAD - UK (NOT SHOWN)	\$141
11012184	CHARGER LEAD - EU (NOT SHOWN)	\$123
11013804	CHARGER LEAD - AU (NOT SHOWN)	\$157
31009237	CLASSIC SMALL SPREADER BAR *	\$151
31009196	CLASSIC STANDARD SPREADER BAR *	\$157
31009197	CLASSIC LARGE SPREADER BAR *	\$170
11014480	CLASSIC SMALL SPREADER BAR COVER (NOT SHOWN)	\$72
11014250	CLASSIC STANDARD SPREADER BAR COVER (NOT SHOWN)	\$76
11014263	CLASSIC LARGE SPREADER BAR COVER (NOT SHOWN)	\$91
39001172	MAJOR 200/MAXI 180 SPREADER BAR FIXING KIT (Major/Maxi Boom end cap x1, Fulcrum pin x1, 10mm circlips x2, Cap retaining washer x2, Plastic cap x2, Hook shaft washer x5 & O-ring x5) (NOT SHOWN)	\$60
39001304	CLASSIC SPREADER BAR SLING RETAINING WASHER KIT (Retaining washer x2, m6 screw x2, spring x2, spacing post x2)	\$46
39001279	CLASSIC SPREADER BAR FULCRUM PIN CAP (NOT SHOWN)	\$22
11014469	MAXI 180 DECAL SET (NOT SHOWN)	\$33
39001198	CLASSIC CIRCLIPS (NOT SHOWN)	\$47
39001182	CLASSIC PUSH HANDLE END CAPS	\$29
39001280	MAJOR 200/MAXI 180 BOOM END CAP	\$143
39001277	CLASSIC 19MM GREY SCREW COVERS (NOT SHOWN)	\$30
39001278	CLASSIC 25MM GREY SCREW COVERS (NOT SHOWN)	\$27
39001397	RECTANGULAR TUBE INSERT FOR LEG ACTUATOR BRACKET - WHITE (NOT SHOWN) <i>Electric leg models only</i>	\$123
39001399	19MM UNI CAP - WHITE (NOT SHOWN) <i>Electric leg models only</i>	\$21

Oxford® Standaid 140

PART NUMBER	ITEM DESCRIPTION	#VALUE!
39001195	FRONT CASTOR KIT: 75mm nylon castor, M8 CSK screw (Pair)	\$67
39001194	REAR CASTOR KIT: 100mm nylon braked castor M10 button head bolt (Pair)	\$105
417641.11001	LEG FABRICATION: Suitable for left & right config (not inc. castors)	\$198
808003.10000	CROSS MEMBER ASSEMBLY, INCLUDING NUTSERT	\$526
EM-2679	LINK ROD	\$38
232388.10014	PIVOT PLATE ASSY: inc. belville washers, circlip, cap, & locking knob	\$114
457641.11001	MAST COMPLETE	\$596
808003.10094	BOOM COMPLETE	\$532
172001.10002	LINAK CONTROLLER	\$654
11012180	LINAK BATTERY *	\$471
39001201	LINAK BATTERY AND CONTROLLER BRACKET AND FIXINGS	\$59
39001183	LINAK BATTERY CHARGER AND DESK TOP UNIT (UK) (NOT SHOWN)	\$431
11012183	CHARGER LEAD - UK (NOT SHOWN)	\$141
11012184	CHARGER LEAD - EU (NOT SHOWN)	\$123
11013804	CHARGER LEAD - AU (NOT SHOWN)	\$157
293388.10000	LINAK HAND CONTROL	\$193

182001.10003	LINAK ACTUATOR (250MM STROKE)	\$935
232388.10015	LEG ADJUSTER HANDLE	\$171
533388.10001	FOOT TRAY FIXINGS AND NON-SLIP PADS	\$674
SA/2569	NON-SLIP PAD - LH	\$36
SA/2570	NON-SLIP PAD - RH	\$36
426388.10000	KNEE PAD AND FIXINGS	\$138
232388.10009	MAST & BOOM FIXING KIT: Nylon washers (x2), Nylon pivot tube (x1), Steel pivot pin (x1), Washers (x2), M8 CSK bolts (x2)	\$39
232388.10011 *Includes item EM-2679	LEG FIXING KIT: Washers (x2), Washer springs (x2), Circlips (x2), Rod link assembly (x2), Leg pivot fabrication (x2)	\$125
232388.10004	ACTUATOR FULCRUM FIXING KIT: Fixing pins (x2), Circlips (x4)	\$22
412.10006	DECAL SET (NOT SHOWN)	\$33
655530.10000	HANDLE GRIP: leg adjuster handle	\$23
384.11000	KNOB LEVER HANDLE FIXING SCREW (NOT SHOWN)	\$5
SA/2692-K	SLING RETAINING CLIP (PAIR)	\$22
EM/3683	STRAIGHT LINE STEERING DEVICE (NOT SHOWN)	\$79
323-12060	HAND GRIP	\$14

