



## VERO PRC

POWER RECLINING CHAIR



### KEY FEATURES

- Sealed two-button power switch embedded on the side of the chair (dual switches and/or wired pendant options available)
- Easy to clean and service
- Durable construction - thoroughly tested to withstand the demands of healthcare environments
- Power-return from Trendelenburg (Trendelenburg models only)
- **4-year limited warranty** on the power recline control module and actuator motor

## Switch to the new Verō PRC

The Verō Power Reclining Chair (PRC) is the newest addition to our popular line of infinitely configurable Verō recliners. Verō PRC is comfortable and convenient for the patient while offering exceptional benefits to the caregiver and added value for facility managers. Caregivers can now assist patients up to 500 lbs at the push of a button with no fear of strain. Trendelenburg models are thoughtfully designed with the capability to **return patients upright from the trendelenburg position**, minimizing potential for caregiver injury.

A sealed power switch is embeded on the patient right side of the recliner for easy access. Add an optional second switch, at patient left side and/or a wired hand pendant for additional patient and caregiver convenience. Verō PRC powers the entire line of Verō recliners offering unmatched configurations, durability, patient comfort and caregiver peace of mind.

Contact Winco today to arrange for a demonstration of the Verō PRC or to request a quote!

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



## SPECIFICATIONS: VERÒ PRC

FOOT OPTION	PEDESTAL FEET		3" DUAL CASTERS		5" CASTERS	
RECLINER SIZE	STANDARD	X-LARGE	STANDARD	X-LARGE	STANDARD	X-LARGE
SEAT HEIGHT TO FLOOR	20.25 in (51.4 cm)		20.75 in (52.7 cm)		23.75 in (60.3 cm)	
SEAT WIDTH	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)	22.5 in (57.2 cm)	26.5 in (67.3 cm)
SEAT DEPTH	23.5 in (59.7 cm)					
OVERALL HEIGHT	44 in (111.8 cm)		44.5 in (113 cm)		47.5 in (120.7 cm)	
OVERALL WIDTH (FIXED ARMS / NO TABLES)*	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)	29.5 in (11.6 cm)	34 in (13.4 cm)
SEAT BACK HEIGHT	28.75 in ( 73 cm)					
OVERALL LENGTH	Upright: 41 in (104.1 cm) Reclined: 69.25 in (175.9 cm)					
ARM HEIGHT	25.5 in (64.8 cm)		26 in (66 cm)		29 in ( 73.7 cm)	
WEIGHT (NO OPTIONS)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)	140 lb (63.5 kg)	158 lb (71.7 kg)
WEIGHT CAPACITY	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)	350 lb (158.8 kg)	500 lb (226.8 kg)

## OPTIONS: VERÒ PRC

- Wired Hand Pendant
- Dual 180 Degree Swing Arms
- Trendelenburg Positioning
- Standard or XL Seat Widths
- 3" Dual Casters, 5" Casters or Pedestal Feet
- Feet Heat Only or Heat and Massage
- Packages Sliding, Weight Bearing Footplate



- Side Tables and/or Pivot Table
- Battery Back-Up for Recline Functions
- IV Pole
- Push Handle
- Accessory Hook
- Custom Armrest and Headrest Covers
- Urethane Arm Caps

+ Many more clinically useful options and accessories - Visit our website for more information

**WARRANTY** Contact Winco or visit our website for warranty details.

**COLORS** Contact Winco or visit our website for available color selections.

**NOTE** Winco reserves the right to change specifications without notice.

### CERTIFICATIONS

- 1) Winco recliners meet ANSI/BIFMA X5.4 for reliability testing
- 2) Winco foam meets California Technical Bulletin 117, Flammability Test Procedure for Resilient Filling Materials used in Upholstered Furniture
- 3) NR or TB options are available. NR meets AB2998 regulations and TB meets NFPA 101 and ASTM E1537-16 Standards for Upholstered Furniture
- 4) PRC models and Heat and Massage, Heat Only Options: Meet electrical safety standards, UL 1647 and CAN/ CSA-C22.2 No.68-92.
- 5) All individual electrical components for Heat and Massage / Heat Only options and power recline certified to IEC 60601-1.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**