

# Tumble Forms 2<sup>®</sup> Feeder Seat<sup>®</sup> System Order Form



Please mark clearly with an X the items you wish to purchase. Fax form to Patterson Medical at **800-547-4333** or call in your order to our Pediatric Helpline at **800-499-0285**. Prices effective from January 1, 2014 through December 31, 2014.

Account No. \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Buyer: \_\_\_\_\_  
Attention: \_\_\_\_\_

Ship to: \_\_\_\_\_  
Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
Phone: (     ) \_\_\_\_\_

## Feeder Seat Positioner Features

- Accommodates infants through adolescents
- Shoulder harness slots allow 4" of vertical adjustment of the shoulder straps to accommodate children of different heights
- Contoured interior has a 90° seat-to-back relationship to provide posturally correct seating
- A 45° hip strap and quick-release H-belt help maintain child's position
- Seamless coverings are washable, odor, urine and stain-resistant and nontoxic

## Mobile Floor Sitter (includes both the Feeder Seat and the Mobile Base)

<u>4444M</u>	Medium, for children up to 48"	<b>\$843.30</b>
<u>4542M</u>	Large for children up to 60"	<b>\$1,017.55</b>
<u>4542P</u>	Extra-Large for adolescents up to 72"	<b>\$1,506.25</b>

## Mobile Base (may be ordered separately)

<u>4444MC</u>	Medium Mobile Base	<b>\$539.95</b>
<u>4542N</u>	Large Mobile Base	<b>\$534.50</b>
<u>4542S</u>	Extra-Large Mobile Base	<b>\$831.70</b>

## Tray for Feeder Seat

<u>2795T</u>	Tray for Small/Medium/Large Feeder Seats	<b>\$1,099.00</b>
<u>2795TXL</u>	Tray for Extra Large Feeder Seats	<b>\$1,099.00</b>

## Choose One Color

Blue	Red	Teal	Purple	Product	Price
4542BB	4542BR	4542BT	4542BP	Deluxe Floor Sitter Small	<b>\$492.30</b>
4542AB	4542AR	4542AT	4542AP	Deluxe Floor Sitter Medium	<b>\$582.00</b>
4542DB	4542DR	4542DT	4542DP	Deluxe Floor Sitter Large	<b>\$751.10</b>
4542EB	4542ER	4542ET	4542EP	Deluxe Floor Sitter Extra Large	<b>\$1,452.50</b>

## Or Mix & Match

Blue	Red	Teal	Purple	Product	Price
2795NB	2795NR	2795NT	2795NP	Small Feeder Seat Positioner	<b>\$376.35</b>
2795MB	2795MR	2795MT	2795MP	Medium Feeder Seat Positioner	<b>\$499.10</b>
2795XB	2795XR	2795XT	2795XP	Large Feeder Seat Positioner	<b>\$649.10</b>
4542WB	4542WR	4542WT	4542WP	Floor Sitter Wedge, fits Small, Medium & Large Feeder Seats	<b>\$339.50</b>
2795YB	2795YR	2795YT	2795YP	Extra Large Feeder Seat Positioner	<b>\$1,055.35</b>
4542YB	4542YR	4542YT	4542YP	Floor Sitter Wedge, fits Extra Large Feeder Seats	<b>\$744.50</b>

## Feeder Seat System Sizing Chart

Dimension	Small Feeder Seat	Medium Feeder Seat	Large Feeder Seat	Extra-Large Feeder Seat
Overall Height	21" (53cm)	27" (68.5cm)	34" (86cm)	41" (104cm)
Outside Width	10½" (26.5cm)	14" (35.5cm)	17" (43cm)	21" (53cm)
Inside Depth	8" (20cm)	9½" (24cm)	13¼" (33cm)	16½" (42cm)
Inner Width	7¾" (18.7cm)	9" (22.8cm)	12½" (31.7cm)	15" (38cm)
Inner Height	18" (46cm)	23" (58.5cm)	29" (73cm)	35" (89cm)
Child's Height	up to 36" (91cm)	up to 48" (122cm)	up to 60" (150cm)	up to 72" (183cm)

**Tumble Forms 2<sup>®</sup>**  
a Patterson Brand

**With Anti-Microbial Protection**

Odor Resistance, Increased Durability

