

# Pool and Spa Access Hoist

This versatile product was invented by Ken Wyborn at the same time he invented the Wymo Wheelchair Hoist. We have very positive feedback relating to the installation and use of the Pool and Spa Access Hoist.



This product can help people who have severe disabilities, aged swimmers and those who have minor mobility problems. The carer positions the seat so the user of the hoist can transfer into the seat, then the carer raises the seat and pushes the seat around the post over the pool or spa, the carer then lowers the seat into the pool or spa then at the correct depth and then user of the hoist the swims off the seat.

## Features

- Specially designed carriage with 450mm Clearance
- Seat carriage rotates 360<sup>0</sup>
- 24 volt system using three different type of power systems
- Choice of 3 different types of post bases to suit your requirements
- Choice of heavy duty post & winch system which brings rating to 200Kg
- Sturdy winch unit construction
- Light grey powdercoated stainless steel which blends into pool or spa environment
- Small freight size (0.22m<sup>3</sup>)
- 1 year Warrantee

# The *wymo* Pool & Spa Access Chair Hoist

- \* Affordable for private homes
- \* Comfortable easy access to pool or spa
- \* Contains own power source
- \* Built-in long life batteries (24volt)
- \* Will lift any person with ease (up to 135kg)
- \* 24 volt battery charger
- \* Carer operated
- \* Powder coated stainless steel construction
- \* Can be adapted for spa use with shorter post



Proudly manufactured  
and supplied by :-

**wymo** LIFTING AIDS  
24 Swallow Road  
South Grafton 2460

Phone (066) 43 1606  
Fax (066) 42 6717

See us on the internet  
[www.wymo.com.au](http://www.wymo.com.au)  
Email [wymo@nor.com.au](mailto:wymo@nor.com.au)

Retail  
Balestier: 12 Boon Teck Rd, #01-01 Shaw Corner T: 62500555  
Outram Park MRT : 100 Eu Tong Sen St #02-34 Pearls Ctr T: 63247077  
Paya Lebar MRT: 10 Eunos Rd 8 #B1-02 Sing Post Ctr T:6455 1551  
Upp Thomson: 23 Tagore Lane #04-19 Tagore 23 Warehouse  
T:250 0555  
Website : [www.RehabMart.sg](http://www.RehabMart.sg) Email : [rehab@pacific.net.sg](mailto:rehab@pacific.net.sg)

AVAILABLE FROM

*Rehab Mart*