




ReOrder# 912428



ReOrder# 912441

Re-Order #	Description	Size	Unit
912428	Spiral Gel Activity Pad	20" x 22"	Each
912441	Spiral Gel Activity Pad used with Light Box	15" x 15"	Each



Purpose

User guides plastic disks around gel swirl for tactile and visual stimulation, enhanced eye-hand coordination, increased finger strength, improved hand dexterity, and greater attention span. The gel surface and foam base protect against injury caused by uncontrolled motor activity.

Easy Care


These vinyl pads are heat sealed and can be cleaned with a mild, non-chlorinated detergent and can be disinfected with a mild, non-chloride-containing disinfectant. Although the gel not toxic, medical assistance must be invoked if the gel is taken through the mouth. If the product is a crack or damage, it must be discarded.

Caution

Possible choking hazard. Do not use for persons cognitively under three (3) years old.

Other Gel Sensory Activity Pads from **SkiL-Care**TM CORPORATION

Gel-Wave Pad



Reorder# 912431

User moves plastic disks through gel waves

Tic-Tac-Toe Gel Pad



Reorder# 912435

A fun way to fine-tune coordination and increase attention span. Allows two people to participate in activity.

Manufactured by: SkiL-Care Corporation • Made in the U.S.A