




FOLDABLE LIGHTWEIGHT POWER CHAIR

SQPC-07B LIGHTWEIGHT LITHIUM BATTERY ELECTRIC CHAIR		Picture
Function Anodized 6063 aluminum frame Flip-up armrest and folded footrest, 7"x2" shock absorbing front caster, with brushless motor and controller, with lithium battery, detachable back cushion, detachable seat cushion, washable seat cushion cover. The chair could be folded within 3 seconds without any tools.		
Specification		
Seat Width:	410mm	
Seat Height:	410mm	
Backrest Height:	470mm	
Net Weight:	23.4kgs with three battery	
Weight Capacity:	120kgs	
Caster:	7"x2" shock absorbing rubber caster	
Rear Wheel:	8" x2" width caster	
Motor:	180WX2 Brushless motor	
Battery:	Lithium 24V 7AH x 1, 2 or 3	
Max speed:	6kmph	
Gradeability:	12degree	
Charging period:	5 hours	
Endurance:	18 / 36 / 54 KM (1, 2 or 3 batteries)	
Overall size:	H940mm X W565mm X L920mm	
Folded size:	655 x 560 x 300 mm	
Carton Size:	710 x 610 x 350mm	
20'FCL	159PCS	
40'HQ FCL	408PCS	
Extra Battey		



Retail Outlets (Rehab Mart Homecare): **Balestier** 12 Boon Teck Rd, #01-02 Shaw S 329586 Tel: 6.2500555 **Chinatown** 100 Eu Tong Sen St #02-34 Pearls Ctr T: 63247077 **Paya Lebar MRT** 10 Eunos Rd 8 #B1-02 Sing Post Ctr T:6455 155 **North** 23 Tagore Lane #04-19 Tagore 23 Warehouse S 787601 T:6 2500555 **Website** : www.RehabMart.sg