


Parts Manual

Oxford® Calibre



NO	PART NUMBER	DESCRIPTION	QTY
1	0Y0066	100mm rear castor (braked)	2
2	0Y0917	Straight line steering device	1
3	0Y0073	100mm front castor (non-braked) (including wrench)	2
4	0Y0871	Front castor bracket assembly (including castors & fixings)	1
5	0Y0889	Front castor bracket end caps	2
6	0Y0861	Mast extrusion (including decals)	1
7	0Y0863	Boom extrusion (including decals)	1
8	0Y0865	Push handle (including fixings)	1
9	0Y0866	Casted boom to mast pivot joint (including fixings)	1
10	0Y0867	Cradle assembly (no scale) (including decals)	1
11	OXF-SCALE-LP	Oxford scale - UK	1
	OXF-SCALE-LP-EU	Oxford scale - EU	1
12	0Y0868	Base assembly (including base casting, base plate & fixings)	1
13	0Y0890	Base plate	1
14	0Y0888	Base fixings kit	1
15	0Y0869	Left leg extrusion (including fixings)	1
16	0Y0870	Right leg extrusion (including fixings)	1
17	39001370	Control box (Smart™ Monitor)	1
	0Y0873	Control box (Standard)	1
18	0Y0067	Battery pack	1
19	39001309	User handset (Smart™ Monitor only)	1
	0Y0075	User handset (Standard control box)	1
-	39001368	Calibration handset (Smart™ Monitor only) (not shown)	1
-	39001369	Actuator reset handset (Smart™ Monitor only) (not shown)	1
20	39001316	Smart™ Monitor retrofit kit (inc. control box, battery and user handset)	1
21	0Y0053	Charger and stand assembly (exc. charger lead)	1
22	0Y0071	Charger lead (UK)	1
23	0Y0072	Charger lead (EU)	1
24	0Y0078	Foot push pad assembly	2
25	0Y0875	Actuator (including power lead - not shown)	1
-	0Y0893	Actuator power lead (not shown)	1
26	0Y0884	Actuator fixings kit (upper & lower)	1
27	0Y0069	Battery and controller bracket (inc. fixings)	1
28	0Y0877	Leg opening actuator (inc. fixings)	1
29	0Y0886	Leg pivot assembly - left hand	1
30	0Y0887	Leg pivot assembly - right hand	1
31	0Y0880	Pivot joint caps	2
32	0Y0882	Boom to mast pivot - plastic covers	2
33	0Y0883	Mast to mast pivot - plastic covers	2
34	0Y0885	Boom to scale quick-release pin kit	1
35	0Y0891	Boom end plastic caps - black	2
-	0Y0881	Fixings kit - complete (not shown)	1
-	0Y0892	Cradle wear washer kit (not shown)	1
-	0Y0894	Leg actuator limit switch (not shown)	2
-	0Y0878	Decal set (not shown)	1
-	0Y0128	Castor tightening wrench (not shown)	1

