



” The Roomer S makes it easy to transfer a patient from one room to another without making any modifications to doorways or headers.

Roomer S

Room-to-room transfers without modification to doorways

The Roomer S overhead lift is a portable overhead lift that offers the possibility to transfer a patient from one room to another in one single motion without making any modifications to doorway headers. Its portability ensures that asset investment is greatly reduced as the same lift can easily be used in multiple rooms. The lift has a weight capacity of 220kg and can be used with either a 2, 3 or 4-point suspension.

A comfortable experience

The Roomer S offers optimum comfort for the patient when used with our sling portfolio. Our room-to-room function reduces the number of lifts and minimizes strains on both patients and caregivers.

The Roomer S is suitable for rooms with sliding doors or fire doors as no modification to the doorway is required, and eliminates the need for turntables and rail switches. It offers soft-start, soft-stop and constant lifting speed regardless of patient weight and direction, making lifting a comfortable experience. Charging can be done conveniently through the hand control for optimal efficiency. The Roomer S also



features a unique telescoping hanger bar that offers different widths in one product to ensure any patient experiences a comfortable sitting position.

Lift intelligence

The Roomer S has the possibility to store usage and performance data. With our data and service application and a PC or tablet device, it lets you find out more about the usage and the condition of your lift. Detailed information can be viewed directly in the application or exported to other media for additional analysis.

Technical specifications

Description	Part no.	Max Load	Dimensions	Weight	Battery	Battery Charger
Roomer S 220 kg incl. lift, charger, charging dock, Europe	55300-EU	220kg (484lbs)	[WxLxH] 400x160x390mm [15.7x6.2x15.4"]	10.9kg (24 lbs)	2,3 Ah 2x12V	Input 100-240 V, AC 50-60 Hz Output 27,8V, 0,8A DC
Roomer S 220 kg incl. lift, charger, charging dock, Australia	55300-AU					
Roomer S 220 kg incl. lift, charger, charging dock, Great Britain	55300-GB					
Roomer S 220 kg incl. lift, charger, charging dock, North america	55300-NA					
Roomer S 220 kg incl. lift, charger, charging dock, Japan	55300-JP					
Emergency Lowering	Hanger Bar Option	IP-class	Lifting Distance	Lifting Speed		
Electronic and mechanical	2, 3, 4-point 35, 45 and 55 cm [13.8", 17.7" and 21.7"]	IPX4 (lift) IPX5 (hand control)	2080mm (81,9")	1.8m/ min - 3m/min 70" /min - 118" / min	ROHS-2 approved Tested and approved according to ISO 10535:2006 ETL-listed	