

Altair

Reduce heavy investments

Our new Altair overhead lift is a lightweight portable lift with a discreet and modern design. Its portability ensures greatly reduced asset investments, as the same lift is easily used in multiple rooms. The lift is available in three weight capacities from 150kg to 220kg and up to 300kg. The Altair can be used with either a 2, 3 or 4-point suspension.

Unbeatable comfort

Together with the Human Care sling range, the Altair offers unbeatable comfort. The unique telescopic hanger bar makes it possible to adjust the width according to patient needs. The 3-point suspension offers an even weight distribution and elongated sitting position. The soft-start, soft-stop and constant lifting speed and direction makes lifting a comfortable experience for patients regardless of their weight. The Altair offers two-speed lifting, high battery capacity and convenient battery charging through the hand control.



Lift intelligence

The Altair has the ability to store usage and performance data. With our data and service application and a PC or tablet device, it lets you find out more about the usage and the condition of your lift. Detailed information can be viewed directly in the application or exported to other media for additional analysis.

Technical specifications

Description	Part no.	Safe Working Load	Dimensions	Weight	Battery	
Altair 150 kg incl. lift, charger, charging dock, Europe	55200-EU	150, 220, 300kg (330, 484, 661lbs)	[WxHx L] 400x140x300 mm (15.7x5.5x11.8")	9,9kg (22lbs)	2,3 Ah 2x12V	
Altair 150 kg incl. lift, charger, charging dock, Australia	55200-AU					
Altair 150 kg incl. lift, charger, charging dock, Great Britain	55200-GB					
Altair 150 kg incl. lift, charger, charging dock, North America	55200-NA					
Altair 150 kg incl. lift, charger, charging dock, Japan	55200-JP					
Altair 220 kg incl. lift, charger, charging dock, Europe	55210-EU					
Altair 220 kg incl. lift, charger, charging dock, Australia	55210-AU					
Altair 220 kg incl. lift, charger, charging dock, Great Britain	55210-GB					
Altair 220 kg incl. lift, charger, charging dock, North America	55210-NA					
Altair 220 kg incl. lift, charger, charging dock, Japan	55210-JP					
Altair 300 kg incl. lift, charger, charging dock, Europe	55220-EU					
Altair 300 kg incl. lift, charger, charging dock, Australia	55220-AU					
Altair 300 kg incl. lift, charger, charging dock, Great Britain	55220-GB					
Altair 300 kg incl. lift, charger, charging dock, North America	55220-NA					
Altair 300 kg incl. lift, charger, charging dock, Japan	55220-JP					
Battery charger	Emergency lowering	Hanger bar option	IP-class	Lifting Distance	Lifting speed	
Input 100-240 V, AC 50-60 Hz Output 27,8V, 0,8A DC	Electronic and mechanical	2, 3, 4-point 35, 45 and 55 cm (13.8", 17.7" and 21.7")	IPX4 (lift) IPX5 (hand control)	2080mm (81,9")	1.8m/ min - 3m/min 70"/min - 118"/min	ROHS-2 approved Tested and approved accord- ing to ISO 10535:2006 ETL-listed

Portable Lift Systems



Ensures greatly reduced asset investments, as the same lift can be used in multiple rooms.