



healthtec

Supporting you

1300 362 923

www.healthtec.com.au

## REHABILITATION

# Sliding Top Tilt Table

**Ideal for Rehabilitation Wards, Diagnostic Clinics doing venous testing or Physiotherapy the sliding top allows patients to practice exercises actively.**

Suitable for a broad spectrum of patients where decreased leg strength is a problem, the advantage of having a sliding top allows the patients to increase their strength and endurance safely.

Patients are able to practice independently or with carers and allows active practice for patients undertaking any rehabilitation programme.



### Features

- Electric Height & Tilt Adjustment with hand switch
- Plantar/Dorsi Flexion Footplates
- Patient Lifter Access
- Removable & Adjustable shoulder horns for premium posture support
- Friction Free Slide Top can be locked in various positions
- 3 x Patient Safety Straps supplied as standard
- Footplates can foldaway converting the Tilt Table into an Examination Table when required
- Footplates can allow for single leg squats be performed
- Independently operated Tilt and height adjustment – can be operated at the same time or individually from any position
- Inclination Indicator to measure tilt angle
- Heavy duty industrial lockable castors
- 5 year warranty on frame & electrics
- Australian designed and manufactured

### Specifications

- 200kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 550mm Minimum Height
- 1000mm Maximum Height
- 710mm Table Width
- Infinite tilt adjustment to +85°
- +/- 10° Plantar/Dorsi Flexion Footplates
- Hoist access to allow for patient lifter



### The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© October 2013




Supporting you

www.healthtec.com.au

# PRICE SHEET

## Sliding Top Tilt Table

	CODE	QTY
<b>Product Description</b>		
LynX Sliding Top Tilt Table (710 wide)	53832	
<b>Width Options</b>		
770mm Width Option	-75	
Custom Width Option	6014	
<b>Foam &amp; Vinyl Options</b>		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
Leather Upgrade	6008	
*please refer to our colour chart for vinyl colour options		
<b>Nose Hole Options</b>		
Nose Hole	6011	
Nose Hole Plug	6010	
<b>Table Protection</b>		
Paper Towel Roll Holder	6017	
<b>Extended Warranty</b>		
10 years Frame & Electrics	7800	
<b>Freight</b>		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	



**FRECALL**  
**1300 362 923**  
 info@healthtec.com.au

Healthtec Pty Ltd ABN 94 224 518 379 TGA AR-00106892-6  
 © October 2013

Disclaimer: Due to our policy of constant product improvement, specifications may change without notice. Check with us for the latest specifications and pricing.