



Apex Health Care Mfg., Inc.



No. 6 Industrial 2nd Road, Tou Chiao Industrial Park,
 Min Hsiung Shiang, Chia Yi, Taiwan
 Tel : 886-5-2133889 Fax : 886-5-2133899
 E-Mail : sales@apexcare.com.tw
<http://www.apexcare.com.tw>



Model no: APC-80873 Hydraulic Shower bed

Features:



- . Powered by hydraulic jack.
 - . Reinforced Aluminum Alloy Frame.
 - . Scissor cross frame chassis.
 - . 7° Tilt adjustment on head/foot ends.
 - . Off-set chassis for easy approaching.
 - . Hi-Lo adjustments by foot-pedal.
 - . Medical grade Non-P PVC Mattress.
 - . Safe Working Load: 181kgs (400lbs)
 - . Steering system included.
 - . foldable side rails.
 - . Castors: 360° turning.
- Front: 3" lockable castors
 Rear: 4" lockable castors.

Specifications: unit: mm, () as reference value.



Turning radius: around 2067mm.

