

## Air Mattress (Bubble) 3 inch



APP-B01+APP-P08



XIAMEN SENYANG CO.,LTD.


*Certificated: FDA, ISO13485, CE*

*Rehab Mart*

## Mattress (APP-B01)

- Medical class PVC 0.32mm
- Size : 200cm(L) x 90cm(W) x 7.0cm(H), end-flap 40cm
- Laser-punctured air vents to improve ventilation,
- Mattress weight 1.95 KG
- Weight capacity 135 kg
- Repair kit and hose included
- Suitable for stage one pressure sore

## Air Pump (APP-P08)



## Features

Aluminium pump case inside , 6-8L/min

- static mode for patient transfer and upright position
- low and normal pressure indicator
- pressure adjustable
- easily accessible filter for easy cleaning
- hook for hanging on bed
- easily accessible fuse for easy changing
- Warranty : one year warranty for pump

## Specifications

Power supply: 110V-220V/50-60Hz

Air output: 6-7 liter/min.

Pressure range: 35mmHg - 130 mmHg

Cycle time: 6 min

Fuse: 1A

Power consumption: 7W

Size of pump: 27cm(L) x 13cm(W) x 9cm(H)

Weight: 1.6kg