

 **Netti** VISION FAMILY
COMFORT

Netti[®]
BY ALU REHAB

MAKING EVERYDAY
LIFE BETTER
THE NORWEGIAN WAY



Enable joy of life


My-Netti.com

IT'S THE MOMENTS OF JOY THAT MAKE THE DIFFERENCE

NETTI VISION, BRINGING MOMENTS OF JOY

More moments of happiness for the individual means a general sense of well-being. Healthier users have a better appetite, need less medication, and have a stronger physical condition. More moments of happiness for family, loved ones and carers.

That is what Netti Vision is all about.




THE NETTI SEATING SYSTEM IS SIMPLY MORE COMFORTABLE

User comfort comes from careful design to ensure an optimum interaction between the wheelchair, the cushions and the different accessories. The thoroughly thought out design creates synergies so that the separate solutions support each other and secures a comfortable sitting position in accordance with the 7 principles of seating comfort.

For example pelvic stability is assured by a number of advanced design features that work together:


- the design of the seating cushion
- the shape of the seating plate beneath the cushion
- the pivot points for the leg supports and the back rest
- the back cushion
- the individually selected accessories, such as comfort pads, abduction block or thigh supports

It's this attention to detail that makes a Netti chair feel different. It's what makes the Netti seating experience so good.



MORE ACTIVE AND HEALTHY

Netti Vision enables the wheelchair user to better digestion, more activities and to have closer contact with people who loves them. The Netti Seating System is unique and secures the optimal sitting position throughout the day. Easy to adjust for variation of pressure.




GET MORE OUT IN THE CHAIR

Netti Vision gives the user more golden moments because it is possible – and safe – to travel by car on a sunny day – together with family members or care givers. The daily life is easier for the family members or care givers because it is easy to organize transportation from home to a care centre etc.



COMPACT AND EASY TRANSFER

Netti Vision has innovative features and enable the wheelchair users to an easy transfer to and from a bed, a toilet or an ordinary chair. The foldable frame of Netti Vision gives you more space in the living room and the kitchen. It is easy to fold and store the wheelchair.



THE NEW VISION SERIES COMPACT COMFORT

Tilt range between
-2 and +30 degrees.
This prevents sliding and
offers the opportunity
for variation.

TAILOR MADE INNOVATION


Combining the competence of the leading comfort wheelchair designer from Norway with seating experts from Japan, and the result is Netti Vision; Making everyday life better the Norwegian way.

Netti Vision is a dream to work with for the therapists, the care givers and the care managers. It is lighter, has innovative and simple ways to adjust the features and it has a modern design.

MORE COMPACT



Folds easily for storage and transport - a unique feature for a comfort wheelchair.



COLOUR YOUR LIFE

The range of patterns and colors can vary and be changed due to availability and order size.


CLINICAL TESTING AND USERS WITH LONG TERM SEATING

Netti Vision can be a tailor made solution for users with complex deceases and the need for long term seating.

Clinical studies shows that users of Netti Vision can have less pain than before, and the users achieve longer seating time with less frequency of sliding.

Report from one Netti Vision user:

"I see the Netti Vision as a therapeutic aid that will enable me to improve my autonomy."



- + Ideally located pivot points enable ergonomic positioning and good posture in all positions.
- + The height of the back cushion can be adjusted for correct lumbar and side support.

- + Rigid frame properties provides stability for the whole body.
- + Seat depth adjustable from 35 cm to 42,5 cm for maximum pressure distribution and stability.
- + Arm support can be lowered to seat cushion level for easy transfer.
- + Cushions provides maximum pressure distribution, comfort and relief.
- + Integrated side supports increase stability.

Foot operated parking brake and anti tippers.

Seat width	Seat depth	Back height	Seat height	Max load	Weight	Tilt	Back recline
40 cm	35-42,5 cm	45 cm	38,5 cm	100 kg	22 kg	-2°-30°	90°-110°

COMFORTABLE AND SUPPORTIVE SEATING

Netti Vision is available with the entire range of Netti Cushions. The Netti seat and backrest cushions are ideal for users who have to sit in a chair for long periods because they assure good pressure distribution and micro climate, which is essential for a healthy skin.

Below is an extract from our seat and back cushion range. The entire product portfolio and technical specifications is at My-Netti.com



Netti Uno Back

The comfortable entry model. One foam piece with basic ergonomic and integrated side support.



Netti Smart

High quality multilayer back cushion with moderate side supports and integrated lumbar support.



Netti Stabil

A backrest cushion similar to the Netti Smart but with larger side supports. For user with tendency to tilt to the side.



Netti Kyphotic

Special designed multilayer back cushion for people with kyphosis.



Netti Uno Seat

Multilayer seat cushion with high pressure distribution properties.



Netti Sit

Multi layer seat cushion with integrated side supports, protection cover and high pressure distribution properties.

EVEN OUR
STANDARD
QUALITY
IS AMONG
THE BEST



Netti Nordic

RELIEF AND STABILITY

Pressure distribution and side support properties with visual design elements to make positioning easy. Multilayer back specially fitted to be used together with Netti Vision.

Netti Contour

POSITIONS HIPS FOR BETTER POSTURE

Anatomic shaped seat cushion with high pressure distribution properties and low seat height.



DESIGNED IN NORWAY FOR ULTIMATE COMFORT

Our range of ergonomic comfort chairs is the most comprehensive available. Netti wheelchairs epitomise advanced Norwegian technology with a rigid frame, high quality cushions and advanced fabrics. This is what gives them their outstandingly high performance. Naturally, all our components are manufactured to exacting Norwegian quality standards, they meet all applicable European norms and ISO 9000 and 13485.



This is a selection of our other solutions.
Learn more at My-Netti.com.



Netti[®] 4U FAMILY COMFORT

VARIETY OF ADJUSTMENT POSSIBILITIES

The range of the Netti 4U series starts with entry models up to medium advanced comfort wheelchairs.



Netti[®] MOBILE FAMILY COMFORT

THE COMFORTABLE POWER WHEELCHAIR

Finally, a power wheelchair offering excellent seating and posture possibilities.



Netti[®] III FAMILY COMFORT

HIGHEST ADAPTABILITY FOR INDIVIDUAL COMFORT

Netti III Comfort is the most adaptable and adjustable choice for achieving maximum comfort.



Netti[®] DYNAMIC SYSTEM COMFORT

MAXIMUM FREEDOM OF MOVEMENT

Accommodates the extension patterns in order to provide the best seating solution for users having unvolunteer muscle contractions.

UK

TOGETHER WE CREATE
ELEGANT SOLUTIONS TO
ENABLE JOY OF LIFE

FEEL FREE TO CONTACT US:



Alu Rehab AS
Bedriftsvegen 23
4353 Klepp Stasjon
Norway
T: +47 51 78 62 20
post@My-Netti.com