



The MULTI-LIFT lifts a driver or passenger in and out of an ordinary car, van truck or SUV. The MULTI-LIFT can also be used as a scooter/ wheelchair lift in certain vans.

The lift is powered by the electrical system of your car. A handy magnetic lift control box, attached to a spiral cord, allows flexible entry and exit and easy storage of the box on any metal surface. Sideways movements are done manually.

Available in various types and sizes, our "Backsaver" sling design requires no lifting in order to put the sling on.

If the lift is permanently removed, the car remains virtually unaffected, maintaining maximum resale value of your vehicle.

Installation can be done by over 2000 certified Access Unlimited dealers across North America and around the world.



• **Multi-Lift is adaptable to most cars, vans, and trucks.**

• **The lift arm weighs 23 lbs. and is easily lifted off its mounting bracket, when not in use. You can use the same lift for several cars or even home stations!**

• **The lift can remain on the bracket while driving or riding.**

• **Lifting capacity is maximum 275 lbs.**



570 Hance Rd Binghamton NY 13903

1-800-849-2143

www.accessunlimited.com

moreinfo@accessunlimited.com