



HD

SPECIFICATIONS MODEL 65000

TRAVEL SURFACES

- Grass, gravel, sidewalks, linoleum, carpeting, indoor and outdoor surfaces

SPEED AND POWER

- 0-5 mph (8.0 kph)
- 30-mile (48.3 km) battery capacity*
- 24V electrical system for two 33/35 amp, airline-safe, sealed batteries**
- 24V onboard battery charger with retractable AC cord reel

HEAVY DUTY UNIPAC® & VARI-SPEED DRIVE SYSTEM

- Sealed, lifetime-lubricated, drive system
- No belts/no chains

MANEUVERABILITY

- 38-inch (96.5 cm) turning radius
- 3.5-inch (8.9 cm) underbody ground clearance

DIMENSIONS

- 47.8 inches (121.4 cm) long
- 23.5 inches (59.7 cm) wide (rear wheel width)
- 39.5 inches (100.3 cm) high
- 186 pounds (84.4 kg) (with batteries)

COLOR

- Cinnamon Red ■

SAFETY

- Electromagnetic braking with manual release
- Extended seat lock with bend
- Flip-up, adjustable, removable arms
- Black rear safety bumper with stability wheels
- ASIC MPC electronics controller
- EZ Tilt handle with HI/LO speed controls
(TouchTron handle available)

PROVIDES

- Steel, painted, one-piece extended platform
- Dynamic, regenerative & electromagnetic brakes
- Pneumatic 9 inch front & 10.5 inch rear tires
- 20-inch (50.8 cm) vinyl folding seat with heavy duty seat post
- 360-degree seat rotation
- 500-pound (227 kg) weight capacity
- 25-inch (63.5 cm) adjustable black handle
- Front universal basket option
- Amigo Limited Lifetime Warranty

* Range based on a 180-pound (81 kg) driver on a flat surface. Range depends on battery type, driver weight, terrain traveled, ambient temperature, and user's charging routine.

** Batteries not included.

AMIGO MOBILITY INTERNATIONAL, INC.

6693 Dixie Highway • Bridgeport, MI 48722-9725
800 248-9130 • www.myamigo.com