

Stylish Design.

Light in Weight & Budget

EB103-A is the latest model designed for the European market. Made for everyday use, indoor & outdoor. Lightest dismantled part weighs only 23 Kg **\$3270**



■ Specifications

Item no	LY-EB103-A
Seat width	A:18";J:16"
Seat Length	17"
Seatback Height	17"
Seat Base Height	18"
Width	A:27";J:25"
Height	36"
Length	41.5"
Motor Output	150W×24V×2pcs
Battery	21AH×12V×2pcs
Batteries Weight	15kgs (7.5kgs×2)
Charger	24V×3mp
N.W.(Without Battery)	30kgs
N.W.(With Battery)	45kgs
Brakes	Electric-Magnetic Brakes
Max loading	120kgs
Slop Gradeability	10 °
Range	25km
Front Wheels	200×50 PU castor
Rear Wheels	2.5" ×6" pneumatic tyre
Maximum Speed	8.5 km/h
Brand of Controller	Dynamic

Rehab Mart

One Stop Home Care

www.RehabMart.com.sg

Tel: 62500555 Fax: 65-62507576 Email: rehab@pacific.net.sg

Retail:

Balestier: 12 Boon Teck Rd, #01-02 Shaw Corner T: 62500555

Outram Park MRT : 100 Eu Tong Sen St #02-34 Pearls Ctr
T: 63247077

Paya Lebar MRT: 10 Eunos Rd 8 #B1-02 Singapore Post Ctr
T:6455 1551