

# Far Infrared Heating Pad

The Far Infrared Heating Pad is used for superior Far Infrared(FIR)Technology- a proven form of heat therapy that works more efficiently and effectively than traditional heating techniques,90% R.H. Best for pain/stress relief-muscle relaxation, Improvement of blood circulation and promote healing.

## Features:



- Comes complete with different braces, Temperature Controller, Heating Pad, AC/DC Adaptor.
- There are Shoulder, Elbow, Hand & Wrist, Back, Thigh, Calf, Knee, Ankle, Foot pads, which is ergonomically designed to fit the body part.
- Designed from 100% carbon fiber fabric which in itself is an electrically conductive material that naturally emits infrared energy to support therapeutic application.
- Ergonomic, breathable design provides support, comfort and 100% coverage from a second injury. Remove heating pad before washing, hand wash the brace by neutral cleaning solution and hang dry, do not bleach, dry clean, iron, wring or tumble dry.
- For the back/hip/rib brace, the USB port can be supplied, heating therapy by rechargeable battery, more convenient.

### 3-Stage Temperature Cycle

12V, 1.5A, 30 minute auto shut-off for safety

- Low Temperature (40-45°C)
- Medium Temperature (45-50°C)
- High Temperature (50-55°C)

## Body Part



## Indication

Heat is generally used for chronic injuries or injuries that have no inflammation or swelling and it helps relax tight muscles or muscles spasms.

## Far Infrared

Far infrared is part of the spectrum of light emanating from the sun. Far infrared is not visible to the eye, but you feel it as heat penetrating into your skin. It believed to retain body heat and emit safe far-infrared rays that penetrate deep into the joint and surrounding tissues which increase blood flow, reduces inflammation and decreases pain. Therefore, FIR was called as *Light Life*.