

User Information	Dealer Information
Order No.	Order No.
Order Date:	Order Date:
Client Name:	Client Name:
Phone:	Phone:

Flexx

Versatile design for Heavy Duty users



Standard Features

- 170 kg Maximum Weight Capacity
- Self Propulsion with 22" or 24" Q-Release Rear Wheels
- Transit Propulsion with 14" Q-Release Rear Wheel
- 6", 7" or 8" Caster Wheels
- Seat Widths: 505/ 555/ 605mm
- Seat Depths: 405/ 455/ 505mm
- Angle Adjustable Backrest: -4°/ 0°/ 4°/ 8°/ 12°/ 16°
- Backrest Heights: 430mm/ 510mm
- 80° Swing in/Swing out Legrest with Flip-up Alu. Footplates
- Push to Lock Brakes with Brake Extension

◆ Standard feature ○ Free of Charge

Seat				
Seat Width and Seat Depth				
Width Depth	505mm(20")	555mm(22")	605mm(24")	
405mm(16")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
455mm(18")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
505mm(20")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Seat-to-Floor Height				
Front Caster Rear Wheel	6"	7"	8"	Caster Positions
14"	<input type="checkbox"/> 18"/ L1 <input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2	<input type="checkbox"/> 20"/ L3	S= Short Fork L= Long Fork H= for 21" seat height
22"	<input type="checkbox"/> 18"/ L1 <input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2 <input type="checkbox"/> 20"/ L3	<input type="checkbox"/> 20"/ L2	
24"	<input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2 <input type="checkbox"/> 20"/ L3	<input type="checkbox"/> 20"/ L2 <input type="checkbox"/> 21"/ H2	
Seat Upholstery				
<input type="checkbox"/> Padded Seat Upholstery				◆
Backrest				

Backrest	
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 430mm Back Height Fixed Backrest	◆
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 510mm Back Height Fixed Backrest	○
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 430mm Back Height Foldable Backrest	\$
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 510mm Back Height Foldable Backrest	\$
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Height Adjust 380mm to 460mm Fixed Backrest	\$
Backrest Upholstery	
<input type="checkbox"/> Without Backrest Upholstery	○
<input type="checkbox"/> Padded Backrest Upholstery	◆
<input type="checkbox"/> Backrest Upholstery for Height Adjustable Backrest	\$
<input type="checkbox"/> Strap Tension Adjustable Backrest (Including Detachable Padded Cover)	\$
Legrests/Footplates	
Legrest RH	
<input type="checkbox"/> 80° RH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) RH Swing in/Swing out	\$
<input type="checkbox"/> Amputee Legrest RH	\$
Legrest LH	
<input type="checkbox"/> 80° LH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) LH Swing in/Swing out	\$
<input type="checkbox"/> Amputee Legrest LH	\$
Footplate RH (Aluminum Footplate is recommended while 80° Swing in/Swing out Legrest is ordered)	
<input type="checkbox"/> Flip-up RH Aluminum Footplate	◆
<input type="checkbox"/> Flip-up RH Tube-in-center Footplate	○
<input type="checkbox"/> Flip-up RH Angle Adjustable Footplate	\$
Footplate LH (Aluminum Footplate is recommended while 80° Swing in/Swing out Legrest is ordered)	
<input type="checkbox"/> Flip-up LH Aluminum Footplate	◆
<input type="checkbox"/> Flip-up LH Tube-in-center Footplate	○
<input type="checkbox"/> Flip-up LH Angle Adjustable Footplate	\$
Armrests	
Armrest RH	
<input type="checkbox"/> Swing-away Ergo Armpad RH Adj. in Height	◆
<input type="checkbox"/> Swing-away Boat Armpad RH Adj. in Height	\$
<input type="checkbox"/> Swing-away Hemi Armpad RH Adj. in Height	\$
Armrest LH	
<input type="checkbox"/> Swing-away Ergo Armpad LH Adj. in Height	◆
<input type="checkbox"/> Swing-away Boat Armpad LH Adj. in Height	\$
<input type="checkbox"/> Swing-away Hemi Armpad LH Adj. in Height	\$
Front Castors/Forks	
Castors	
<input type="checkbox"/> 150mm x 30mm Solid Tire (6")	○
<input type="checkbox"/> 175mm x 30mm Solid Tire (7")	◆
<input type="checkbox"/> 200mm x 30mm Solid Tire (8")	○
<input type="checkbox"/> 200mm x 50mm Solid Tire (8")	\$
Forks	
<input type="checkbox"/> Stylish Fork	◆
<input type="checkbox"/> 21" Seat Height Long Fork	○
<input type="checkbox"/> 2" Wide Long Fork (While 200mm x 50mm Tire is ordered)	\$
Rear Wheels	
Rear Wheel	

<input type="checkbox"/> 14" Solid Tire Plastic Rim	○
<input type="checkbox"/> 14" Pneumatic Tire Plastic Rim	\$
<input type="checkbox"/> 22" Pneumatic Tire Spoked Rim with Plastic Handrim	\$
<input type="checkbox"/> 22" Pneumatic Tire Spoked Rim with Aluminum Handrim	\$
<input type="checkbox"/> 24" Solid Tire Spoked Rim with Aluminum Handrim	◆
<input type="checkbox"/> 24" Pneumatic Tire Spoked Rim with Ergo Handrim	\$
<input type="checkbox"/> 24" Pneumatic Tire Spoked Rim with Aluminum Handrim	\$
Rear Wheels Axle	
<input type="checkbox"/> Fixed Axle (Only available for 14" rear wheel)	○
<input type="checkbox"/> Quick Release Axle (Karma version)	◆
Brakes	
<input type="checkbox"/> Push to Lock Brakes with Brake Extension	◆
<input type="checkbox"/> Attendant Brakes	\$
<input type="checkbox"/> Attendant Hub Brakes	\$
Safety Options	
<input type="checkbox"/> 4-Point Tie-down Brackets	\$
<input type="checkbox"/> Swing-away Height Adjustable Anti-tippers Pa	\$
Accessories	
<input type="checkbox"/> Lap Belt (Safety Belt)	\$
<input type="checkbox"/> Cane Holder	\$
<input type="checkbox"/> IV Pole Holder	\$
<input type="checkbox"/> Oxygen Tank Holder	\$
<input type="checkbox"/> Air Pump	\$
Price	
Total price	