

RollerSlide

SystemRoMedic.

*For supine
transfers*

Bridges gaps

*Moves with the user
during the transfer*



RollerSlide is used for supine transfers of inert users between a bed and an x-ray or operating table, shower trolley, ambulance stretcher or similar. RollerSlide consists of a padded sliding board with a cover made from ULF (Ultra Low Friction), our own sophisticated sliding material that slides round the board during transfer. RollerSlide is available in several versions and sizes.

RollerSlide

How can a bridge move and still be safe?

Differences in height can pose a big problem. And sometimes it is difficult to get close enough to ensure smooth transfers for users and staff. RollerSlide, which is used for supine transfers, bridges gaps and moves with the users during the transfer.

RollerSlide consists of a padded sliding board with a cover made from ULF (Ultra Low Friction), our own sophisticated sliding material that slides round the board during transfer. The sliding material reduces the friction under the user, which facilitates the transfer and avoids unnecessary lifting for the staff. RollerSlide also ensures very gentle transfers for users who are injured or very sensitive to pain.

RollerSlide is both easy to use and simple to maintain. The replaceable sliding cover can easily be cleaned and disposable covers provide for even further improved hygiene and less laundry.



RollerSlide is available in several versions and sizes. RollerSlide, flexi is a lighter and bendable version of RollerSlide. And RollerSlide, light is a light-weight version of RollerSlide, flexi. All three versions of RollerSlide are available in two lengths, one which covers the entire body and one which covers the trunk. The longer versions are foldable in the middle.



RollerSlide



RollerSlide, flexi

- A rigid inner board guarantees a stable transfer even when there is a gap and a difference in height.
- RollerSlide bridges gaps of up to 20 cm/8" and evens out any height differences.
- The sliding cover made from ULF, our unique sliding material, rotates around the board and gives a very simple and smooth transfer.
- RollerSlide is foldable which makes it easy to handle and transport.
- RollerSlide has handles at the short ends.
- RollerSlide, flexi can be folded lengthwise along the centerline of the user or the bed and it effortlessly adapts to the surface it rests on.
- RollerSlide, flexi does not require a flat surface and therefore it is very practical for transfers from ambulance stretchers or other uneven surfaces.
- In all other respects, RollerSlide, flexi is similar to the basic RollerSlide model.
- RollerSlideFlexi is patented by DE patent no.103 34 270 and EP patent applied for.

QUICK REFERENCE

RollerSlide is available in several versions and sizes. RollerSlide, flexi is a lighter and bendable version of RollerSlide. RollerSlide, light is a light-weight version of RollerSlide, flexi. RollerSlide consists of a plywood or plastic board, foam padding, nylon and a sliding cover made of ULF (Ultra Low Friction), our own sophisticated sliding material that slides round the board during transfer

Art. No.	Product	Material	Description	Weight
9145	RollerSlide	Plywood	175 x 50 cm/70 x 20", foldable	4,1 kg / 9,0 lbs
9155	RollerSlide	Plywood	85 x 50 cm/34 x 20"	2,2 kg / 4,8 lbs
9148	RollerSlide, flexi	Plastic	175 x 50 cm/70 x 20", foldable	3,9 kg / 8,6 lbs
9158	RollerSlide, flexi	Plastic	85 x 50 cm/34 x 20"	2 kg / 4,4 lbs
9120	RollerSlide, light	Plastic	175 x 50 cm/70 x 20", foldable	2 kg / 4,4 lbs
9110	RollerSlide, light	Plastic	85 x 50 cm/34 x 20"	1,1 kg / 2,4 lbs

Accessories

9136	Sliding cover	For 9120, 9145 and 9148
9156	Sliding cover	For 9110, 9155 and 9158
9138	Disposable cover	For 9120, 9145 and 9148, 2 x 50 pcs
9168	Disposable cover	For 9110, 9155 and 9158, 2 x 115 pcs
1085	WallRack	For convenient storing of sliding mats and disposable covers for EasySlide and RollerSlide

