



Date _____ Customer code C- _____

Invoice address _____ Contact person _____

_____ Reference _____

_____ Other delivery address _____

Invoice address = delivery address

Quote

Basic price: 3.495,- €

excl. VAT

Seat width
D2025 25 cm
D2030 30 cm
D2035 35 cm

Seat depth
A 32,5 cm
B 35,0 cm
C 37,5 cm

Capacity Max. 90 kg
Seat height See page 5

D22152 Side upholstery (reduces SW with 2 x 2 cm) 129,- €

= Standard equipment

FRAME

Item nr.	Description	Price in €
Standard	<input checked="" type="checkbox"/> Push handles 30 cm, height adjustable and swingable	0,-
Standard	<input checked="" type="checkbox"/> Frame brake, brake lever 9 cm	0,-
Standard	<input checked="" type="checkbox"/> Drum brakes with quick release	0,-
D21808	<input type="checkbox"/> Push handles 500 mm, height adjustable and swingable	189,-
D83349	<input type="checkbox"/> Pushbow, height and angle adjustable	437,-
D20824	<input type="checkbox"/> Frame extension incl. brakes on lower frame tube	269,-
D41074	<input type="checkbox"/> Car fastener	158,-

SOLUTIONS FOR RESPIRATION MACHINES AND OTHER MEDICAL DEVICES

D21186	<input type="checkbox"/> Universal device for different medical equipment (incl. frame extender D20824) – only SW 35 cm	1697,-
D21187	<input type="checkbox"/> Cover plate for universal device (additional level)	449,-
D42900	<input type="checkbox"/> Extra set of 4 mounting brackets, incl. screws	234,-


BACK UNIT

Standard	<input checked="" type="checkbox"/>	Back cushion Netti Smart - side support 12 cm - 3D cover	0,-
Standard	<input checked="" type="checkbox"/>	Back cushion cover black	0,-
D30530	<input type="checkbox"/>	Back cushion cover black - red	0,-
Standard	<input checked="" type="checkbox"/>	Back height 40 cm (incl. adjustable velcro back)	0,-
D21006	<input type="checkbox"/>	Back height 50 cm for SW 35 cm (incl. back extender + longer cushion – only with black cover)	294,-
D2080M	<input type="checkbox"/>	Back height 50 cm for SW 25-30 cm (incl. back extender + longer cushion – only with black cover)	On request
D24051		Side support correction mini for Netti SMART <input type="checkbox"/> 1 pc. <input type="checkbox"/> 2 pcs.	319,-
D21023	<input type="checkbox"/>	Back cushion Netti STABIL - with 17 cm side support - 3D cover	497,-
D26830	<input type="checkbox"/>	Side support stable left – for Netti STABIL	219,-
D26843	<input type="checkbox"/>	Side support stable right – for Netti STABIL	219,-
D27541	<input type="checkbox"/>	Upholstery for side support stable, pair	134,-
D20881	<input type="checkbox"/>	H-belt, with upholstery and car belt lock	227,-
D41645	<input type="checkbox"/>	Belt cross harness (we recommend D86890 H-harness adaptor kit)	198,-
D86890	<input type="checkbox"/>	H-harness adaptor kit (from March 2016)	250,-
Libertas 15° back flexibility incl. long anti tips and frame extender (not intended for transport in car) <input type="checkbox"/> D82414 80N <input type="checkbox"/> D82415 120N <input type="checkbox"/> D82416 160N			
Looking for a wheelchair with flexible back, seat, head and leg support – please contact us for further information on Netti Dynamic System.			

SEAT UNIT

Standard	<input checked="" type="checkbox"/>	Netti Sit anti-decubitus seat cushion with incontinence layer under the 3D cover	0,-
Standard	<input checked="" type="checkbox"/>	Solid seat plate	0,-
D21048	<input type="checkbox"/>	Incontinence cover for seat cushion	212,-
D21017	<input type="checkbox"/>	Netti SIT °C - anti-decubitus seat cushion with temperature regulating features and 3D cover	497,-
D63400	<input type="checkbox"/>	Vital Base Venti - anti-decubitus cushion with 3D cover	439,-
D20845	<input type="checkbox"/>	Hip belt, with upholstery and car belt lock	107,-
D86282	<input type="checkbox"/>	Belt bracket kit	98,-
D24297	<input type="checkbox"/>	Abduction block mini (W: 6 cm x L: 11 cm x H: 5 cm)	193,-
D20830		Removable knee and thigh support <input type="checkbox"/> 1 pc. <input type="checkbox"/> 2 pcs.	188,-

LEG RESTS AND ACCESSORIES

Standard	<input checked="" type="checkbox"/>	Universal leg rests with knee angle adjustment and footboard	0,-
D21011	<input type="checkbox"/>	Manual adjustable leg rests with automatic length adjustment, angle adjustable foot plates, calf supports and knee upholstery	383,-
D21081	<input type="checkbox"/>	Extra abduction for leg rests +5 cm	380,-
D2082M	<input type="checkbox"/>	Foot plate holder + 9 cm	56,-
D21041	<input type="checkbox"/>	Foot box upholstered	469,-
D23014	<input type="checkbox"/>	Foot fixation, pair	132,-
		Ankle huggers	
	<input type="checkbox"/>	D44863 S (L: 17-20 cm)	152,-
	<input type="checkbox"/>	D44864 M (L: 19-23 cm)	152,-
	<input type="checkbox"/>	D44865 L (L: 22-29 cm)	159,-
	<input type="checkbox"/>	D44904 XL (L: 28-33 cm)	174,-

ARM SUPPORTS AND ACCESSORIES FOR ARM REST

Standard	<input checked="" type="checkbox"/>	Arm support height and depth adjustable	0,-
Standard	<input checked="" type="checkbox"/>	Arm pads 330 x 55 mm	0,-
D21096	<input type="checkbox"/>	Arm pads 385 x 58 mm	112,-
D21077	<input type="checkbox"/>	Arm pads 420 x 80 mm	112,-
D21078	<input type="checkbox"/>	Arm pads 445 x 70 mm	112,-
D21079	<input type="checkbox"/>	Arm pads 525 x 80 mm	112,-
D21927	<input type="checkbox"/>	Arm pads curved at the back 525 x 80 mm	148,-
D41649	<input type="checkbox"/>	Arm pad upholstery extra soft (for arm pads 385 x 58 mm) – 3D cover (not to be combined with a tray)	99,-
D41650	<input type="checkbox"/>	Arm pad upholstery extra soft (for arm pads 420 x 80 mm) – 3D cover (not to be combined with a tray)	123,-

MAIN WHEEL

Standard	<input checked="" type="checkbox"/>	Drum brake wheels 22 x 1" PU with quick release	0,-
Standard	<input checked="" type="checkbox"/>	4° camber	0,-
D22006	<input type="checkbox"/>	Main wheel size: 1 3/8, flexel	59,-
D22008	<input type="checkbox"/>	Main wheel size: 1 3/8, pneumatic	59,-
D22002	<input type="checkbox"/>	Drum brake wheels 20 x 1" PU	137,-
D22018	<input type="checkbox"/>	Drum brake wheels 24 x 1" PU	137,-
D22009	<input type="checkbox"/>	Hand and spoke protector, transparent, pair	138,-
D22013	<input type="checkbox"/>	Push rim cover, pair	112,-

FRONT WHEELS

Standard	<input checked="" type="checkbox"/>	Front wheels: 6" – 150 x 30 PU	0,-
D22003	<input type="checkbox"/>	Front wheels: 5" – 125 x 27 PU	79,-
D22017	<input type="checkbox"/>	Front wheels: 7" – 175 x 33 PU	79,-
D21021	<input type="checkbox"/>	Front wheels: 7" – 175 x 45 flexel	79,-
D21012	<input type="checkbox"/>	Front wheels: 7" – 175 x 45 pneumatic	79,-
D21026	<input type="checkbox"/>	Front wheels: 8" – 200 x 50 flexel	79,-
D20664	<input type="checkbox"/>	Front wheel distance busher to simplify foot propelling	158,-

HEAD SUPPORTS & ACCESSORIES

D86133	<input type="checkbox"/>	Head support A mini with side support, 350 mm vertical bar	388,-
D23973	<input type="checkbox"/>	Head support A mini, with side support, 200 mm vertical bar	382,-
D20574	<input type="checkbox"/>	Head support B, small, oval	210,-
D21027	<input type="checkbox"/>	Head support adjustment tube (horizontal) + 10 cm, length 25 cm	79,-
D21028	<input type="checkbox"/>	Head support adjustment tube (horizontal) + 23 cm, length 38 cm	89,-
D21054	<input type="checkbox"/>	Head support adjustment tube 500 mm	79,-
D23012	<input type="checkbox"/>	Easy Care cover for head support type A + B – removable	79,-
D23015	<input type="checkbox"/>	Grey head support cushion with Mars cover	214,-

TRAYS

D21050	<input type="checkbox"/>	Tray reinforced for SW 25 - 35 cm	356,-
D21052	<input type="checkbox"/>	Upholstery for tray	113,-

VARIOUS ACCESSORIES



D21087	<input type="checkbox"/>	Clinic-Kit (tools and star wheels)	106,-
D21239	<input type="checkbox"/>	Toolkit	89,-
D42269	<input type="checkbox"/>	Wheelchair bag	140,-
D41987	<input type="checkbox"/>	Intravenous rack	298,-

TECHNICAL DATA

Tilt (seat)	-7° to +23°
Recline (back)	87° to 120°
Total width of the wheelchair	Seat width + 21 cm
Total length	80 cm with leg rests
Total height	62 cm
Weight of wheelchair	App. 25 kg (without cushion)
Frame colour	Silver
Cover variants	Stretchable & flexible cover 3D

Seat height – WITHOUT seat cushion (measured from floor to seat plate)

Main wheel size	20"		22"		24"	
Upper position/front wheel size			41,5 cm	5 & 6"	44 cm	5,6 & 7"
Lower position/front wheel size	42 cm	5,6 & 7"	44,5 cm	6,7 & 8"	47 cm	7 & 8"



Measurements

A Hip width/Seat width	_____ cm
B Thigh length / seat depth	_____ cm
C Arm rest height	_____ cm
D Leg length / seat height	_____ cm
E Backrest height	_____ cm
F Trunk width	_____ cm

Leg rests	Min. calf length	Max. calf length
Seat width 25-35 cm Manual adjustable leg rests	29 cm	39 cm
Seat width 25-30 cm Universal leg rests	22 cm	32 cm
Seat width 35 cm Universal leg rests	27 cm	37 cm