



SOUTHERN CROSS



Supporting you

1300 362 923

www.healthtec.com.au

Examination Table

The SC electric lift Examination Table makes switching to an electric height adjustable table more affordable than ever before.

Offering a generous lifting capacity of 250kg with a lift-up backrest enabling a comfortable seated or semi-reclined position.



Features

- Electric Height Adjustment with footswitch
- X-Frame Lift action minimises pinch and crush points
- Table can be moved similar to moving a wheelbarrow
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 250kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 540mm Minimum Height
- 860mm Maximum Height
- 710mm Table Width
- 1930mm Table Length
- Manual backrest adjustment +70°
- Footswitch Control



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© August 2013




Supporting you

www.healthtec.com.au

PRICE SHEET

Examination Table

Product Description	CODE	Qty
Product Description		
SC Examination Table (710 wide) – 1/3-2/3 Two Section	56121	
Foam & Vinyl Options		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
*please refer to our colour chart for vinyl colour options		
Nose Hole Option		
Nose Hole in Backrest	6011-2-7	
Nose Hole Plug	6010	
Armrest Option		
Procedure Arm	5102	
Gynae Option		
Gynae Leg Supports	5114	
Table Protection		
Paper Towel Roll Holder	6017	
Disposable Fitted Bed Sheet – 10 per pack	ANWFIT	
Extended Warranty		
10 years Frame & Electrics	7800	
Freight		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	



FRECALL
1300 362 923
 info@healthtec.com.au

Healthtec Pty Ltd ABN 94 224 518 379 TGA AR-00106892-6
 © August 2013

Disclaimer: Due to our policy of constant product improvement, specifications may change without notice. Check with us for the latest specifications and pricing.