

02 Rehabilitation class

HW60022F

180W Handy & Feet Indoor Rehabilitation Exerciser Training Bike

Function Description :

voltage: AC 220~ 240V/ 100~120V

HZ: 50Hz

PRODUCT WEIGHT: 15KGS



With Spasm Mode

1. Power : Adjustable Power 80W/100W/120W/180W
2. 51 different Adjustable stepless Speed from 15 RPM ~ 65RPM
3. Spasm Mode
4. Active Training
5. Automatic Training
6. Automatic Active training Speeds : 4 different speeds 15RPM/30RPM/45RPM /60RPM
7. With 2 direction: Forward / Reverse
8. With Big LCD display with Blue Back light on Machine
9. Hands/Arms / Legs/ Foot training
10. Timers function: 5min 10 min 15min 20min 25min
11. Product weight: 19kgs
12. The best electronics upper and lower limbs physiotherapy and rehabilitation exerciser product with adjustable base Rack for hemiplegia and stro
13. Leg foot arm hand could do training together , or could choose to do training separately
14. Machine Height could be adjustable
15. with Manual training as well
17. With Electric protection function.
18. The upper hand part is with Adjustale function from 1~15cm FROM front to back
19. The metal role part is with 1~15 cm adjustable for the height for different height people use.
20. The product has the function of anti falling
21. The rubber foot is with anti-slip function.
22. The leg support could be wrapped safety with feet together when put on the pedal.
23. The hand pedal is soft and anti-slipped
24. Product Size: 61cm * 37 width * 101~116cm height (adjustable)



Big LCD Screen



Hand Arm Exercise



Feet leg Exercise