

Available from May 2022



Date _____ Customer code C- _____
 Invoice address _____ Contact person _____
 _____ Reference _____
 _____ Other delivery address _____

- Invoice address = delivery address Quote
 The chair will be used as seat for car transport

Basic price: 5.795,- € Excl. VAT

*Some accessories will only be available from summer 2022

Seat width	D30035	<input type="checkbox"/>	350 mm	Seat depth	A	<input type="checkbox"/>	355 mm*
	D30038	<input type="checkbox"/>	380 mm		B	<input type="checkbox"/>	380 mm*
	D30040	<input type="checkbox"/>	400 mm		C	<input type="checkbox"/>	405 mm
				D	<input type="checkbox"/>	430 mm	
				E	<input type="checkbox"/>	455 mm	

Seat width	D30043	<input type="checkbox"/>	430 mm	Seat depth	A	<input type="checkbox"/>	405 mm*
	D30045	<input type="checkbox"/>	450 mm		B	<input type="checkbox"/>	430 mm*
	D30050	<input type="checkbox"/>	500 mm		C	<input type="checkbox"/>	455 mm
				D	<input type="checkbox"/>	480 mm	
				E	<input type="checkbox"/>	505 mm	

Capacity Max. 135 kg

Seat height See page 6

*Please take note of the possible conflict between the Netti Grandis leg support and front castors when the seat plate is horizontal and the leg supports are adjusted in 90°.

■ = Standard equipment ■ = included in introduction campaign

Special offer

FRAME

Item nr.	Description	Price in €
Standard ■	Wheel footprint – frame length and anti-tippers adjust automatically when chair is tilted	0,-
Standard ■	Push bow height adjustable and foldable	0,-
Standard ■	Anti-tippers with pedal	0,-
Standard ■	Car fastener	0,-
Standard ■	User brake	0,-
Standard ■	Drum brakes	0,-
D30824* □	Frame extension bracket for main wheels (main wheels moved 60 mm backwards)	290,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
 DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300

Available from May 2022

D21022	<input type="checkbox"/>	Easy Care cover for seat and back – wipeable (not available for Netti Sit °C, Kyphotic, Ride cushions and Netti Relaxo head support)	0,-
--------	--------------------------	--	-----

BACK UNIT

Standard	<input checked="" type="checkbox"/>	Back height 500 (incl. adjustable velcro back)	0,-
Standard	<input checked="" type="checkbox"/>	Built-it back extender for back height 610 mm (longer cushion to be ordered separately)	0,-
Standard	<input checked="" type="checkbox"/>	Back cushion Netti Smart - side support 120 mm - 3D cover in black	0,-
D37376*	<input type="checkbox"/>	Side support correction with upholstery left - for Netti Smart	235,-
D37378*	<input type="checkbox"/>	Side support correction with upholstery right - for Netti Smart	235,-
D21023	<input checked="" type="checkbox"/>	Back cushion Netti Stabil - with 170 mm side support and 3D cover in black/gold (value 497,-)	0,-
D21123	<input type="checkbox"/>	Back cushion Netti Super Stabil – with 170 mm reinforced side support and 3D cover in black/gold (black version available now)	579,-
D26831	<input type="checkbox"/>	Side support stable left – for Netti Stabil, Super Stabil & Kyphotic Please note: Can only be used when back extender is in top position or removed.	240,-
D26844	<input type="checkbox"/>	Side support stable right – for Netti Stabil, Super Stabil & Kyphotic Please note: Can only be used when back extender is in top position or removed.	240,-
D31006	<input type="checkbox"/>	Back cushion delivered in back height 600 mm – according to cushion choice: Netti Smart, Stabil or Super Stabil delivered with black 3D cover	135,-
D20630	<input type="checkbox"/>	Back cushion Netti Kyphotic - black 3D/Mars cover - only back height 500 mm (back height 600 mm in special adaptation)	667,-
	<input type="checkbox"/>	Ride Designs Back solution* (See separate order form) Article number: _____ / Accessories: _____	
D21042	<input type="checkbox"/>	Lumbar support	164,-
D21084	<input type="checkbox"/>	Back wedge	164,-



Netti Dynamic back cylinder

Allows up to 30° trunk movements backwards. The standard back recline can be locked in any position as the dynamic movement is blocked simultaneously, which can be appropriate during transport or positioning. Including Bodypoint Evoflex pelvic stabilizer 2-point belt: Please choose size: XS M L

For more information and choice of strength of gas spring, feel free to contact us. If combined with other options, please send us a request as some solutions are only possible as special adaptations.

D32231	<input type="checkbox"/>	300N
D32232	<input type="checkbox"/>	450N
D32233	<input type="checkbox"/>	600N
D32234	<input type="checkbox"/>	750N

Price on request

Available from May 2022



For a complete Dynamic System with dynamic seat unit, please see the Netti Dynamic AdaptPro order form.

SEAT UNIT

Standard	<input checked="" type="checkbox"/>	Seat plate with pressure reduction properties. Easy to adjust for correct positioning	0,-
Standard	<input checked="" type="checkbox"/>	Rail for mounting of pelvic belts	
Standard	<input checked="" type="checkbox"/>	Seat cushion Netti Sit – pressure relief seat cushion with incontinence layer under the 3D cover	0,-
D21017	<input type="checkbox"/>	Seat cushion cover °C for Netti Sit cushion with temperature regulating features – only 3D cover	530,-
D21048	<input type="checkbox"/>	Incontinence cover for seat cushion	212,-
	<input type="checkbox"/>	Ride Designs Seat cushion (See separate order form) Article number: _____ / Accessories: _____	
D37490*	<input type="checkbox"/>	Abduction block small	250,-
D37491*	<input type="checkbox"/>	Abduction block medium	250,-
D37492*	<input type="checkbox"/>	Abduction block large	250,-

LEG SUPPORTS

Netti Grandis angle adjustable leg support

Standard	<input checked="" type="checkbox"/>	Netti Grandis angle adjustable positioning leg supports for complex individual configurations. Reduces pressure and shear forces. For users with e.g. contractions, leg asymmetry, elephantiasis.	0,-
			
		The foot support is characterized by unique ergonomics and user-friendliness due to lifting handles, scale displays and conical extension pieces. Angle adjustable incl. angle limiter and adjustable calf- and foot supports.	
D88602	<input type="checkbox"/>	Calf support upholstery extra soft - pair - EC or 3D cover depending on the fabric for the whole chair	190,-
D88607	<input type="checkbox"/>	Knee upholstery for the upper lower leg (3D or EC cover depending on the fabric for the whole chair)	102,-
D21745	<input type="checkbox"/>	Calf support bracket upholstery - pair - skai cover	49,-
D88014	<input type="checkbox"/>	Foot fixation, pair incl. mounting kit	212,-
D88603	<input type="checkbox"/>	Foot board, lockable and upwards swingable	347,-
D88605	<input type="checkbox"/>	Footbox, upholstered with Skai cover (approved for car transport)	654,-
D21293	<input type="checkbox"/>	Calf strap	79,-
D86900	<input type="checkbox"/>	Bodypoint calf strap M	162,-
D86901	<input type="checkbox"/>	Bodypoint calf strap L	162,-
			
		Bodypoint ankle huggers with upholstery	
	<input type="checkbox"/>	D88863 S (L: 170-200 mm)	293,-
	<input type="checkbox"/>	D88864 M (L: 190-230 mm)	293,-
	<input type="checkbox"/>	D88865 L (L: 220-290 mm)	303,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300


Available from May 2022

<input type="checkbox"/>	D88866 XL (L: 280-330 mm)	324,-
--------------------------	---------------------------	-------



Netti Grandis amputation support

D88611	<input type="checkbox"/> Amputation support left (amputation support on both sides only together with frame extension D30824) incl. Netti Grandis leg support right side	453,-
D88610	<input type="checkbox"/> Amputation support right (amputation support on both sides only together with frame extension D30824) incl. Netti Grandis leg support left side	453,-

Netti Grandis leg support with fixed angle

D88609	<input type="checkbox"/> Netti Grandis leg supports with fixed angle for complex individual configurations. Reduces pressure and shear forces. For users with e.g. contractions, leg asymmetry. The foot support is characterized by unique ergonomics and user-friendliness due to scale displays and conical extension pieces. Adjustable in a fixed angle – also below 90°. Including heel strap.	117,-
		
D88603	<input type="checkbox"/> Foot board, lockable and upwards swingable	347,-
D88605	<input type="checkbox"/> Footbox, upholstered with Skai cover (approved for car transport)	654,-

Manually angle adjustable leg supports



D21111	<input type="checkbox"/> Manually angle adjustable leg supports with -automatic length adjustment -angle adjustable foot plates -calf supports and knee upholstery -reinforced extension pieces	117,-
		
D88211	<input type="checkbox"/> Leg support upper parts with lever instead of star wheel for angle adjustment	80,-
	Bodypoint ankle huggers with upholstery	
	<input type="checkbox"/> D44863 S (L: 170-200 mm)	213,-
	<input type="checkbox"/> D44864 M (L: 190-230 mm)	213,-
	<input type="checkbox"/> D44865 L (L: 220-290 mm)	223,-
	<input type="checkbox"/> D44904 XL (L: 280-330 mm)	244,-
D80087	<input type="checkbox"/> Depth adjustable hinge for footplates (50 mm in all) – mounted	444,-
D21008	<input type="checkbox"/> Gas spring	878,-
D21044	<input type="checkbox"/> Foot board, adjustable	287,-
D21041	<input type="checkbox"/> Foot box upholstered (not approved for car transport, if this is needed, please use the Netti Grandis solution which is approved)	459,-
D24000	<input type="checkbox"/> Calf support, swingable, pair	165,-
D21745	<input type="checkbox"/> Calf support bracket upholstery - pair - skai cover	49,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300

Available from May 2022



D41617	<input type="checkbox"/>	Calf support upholstery extra soft - pair - EC or 3D cover depending on the fabric for the whole chair	190,-
D28453	<input type="checkbox"/>	Foot plate upholstery - pair - skai cover	269,-
D86900	<input type="checkbox"/>	Bodypoint calf strap M	162,-
D86901	<input type="checkbox"/>	Bodypoint calf strap L	162,-

Universal leg supports

D24006	<input type="checkbox"/>	Universal leg supports with knee angle adjustment, adjustable below 90° for contract leg situations with seat unit tilted when chair is pushed. Without calf supports.	117,-
			
		Bodypoint ankle huggers with upholstery incl. mounting kit	
	<input type="checkbox"/>	D44863 S (L: 170-200 mm)	213,-
	<input type="checkbox"/>	D44864 M (L: 190-230 mm)	213,-
	<input type="checkbox"/>	D44865 L (L: 220-290 mm)	223,-
	<input type="checkbox"/>	D44904 XL (L: 280-330 mm)	244,-
D80087	<input type="checkbox"/>	Depth adjustable hinge for footplates (50 mm in all) – mounted	444,-
D21044	<input type="checkbox"/>	Foot board, adjustable without heel strap	287,-
D21041	<input type="checkbox"/>	Foot box upholstered (not approved for car transport, if this is needed, please use the Netti Grandis solution which is approved)	459,-
D86900	<input type="checkbox"/>	Bodypoint calf strap M	162,-
D86901	<input type="checkbox"/>	Bodypoint calf strap L	162,-

Netti Dynamic angle adjustable leg rests



		Netti Dynamic angle adjustable leg rests DUAL with calf supports and footboard, 100N.	Price on request
D82296	<input type="checkbox"/>	Adjustment piece 300 mm (lower leg length 450-550 mm)	
D82298	<input type="checkbox"/>	Adjustment piece 214 mm (lower leg length 350-450 mm)	
			
		<u>Adjustable for different thigh lengths:</u>	
D82318	<input type="checkbox"/>	Adjustment piece 300 mm (lower leg length 450-550 mm)	
D82319	<input type="checkbox"/>	Adjustment piece 214 mm (lower leg length 350-450 mm)	
D45245	<input type="checkbox"/>	Reinforced automatic footboard lock - Not possible to combine with D82318/9 (extension bracket). With this lock the footboard is not angle adjustable, but this can be offered as special adaptation.	Price on request
D88211	<input type="checkbox"/>	Leg support upper parts with lever instead of star wheel for angle adjustment	80,-
D32296	<input type="checkbox"/>	Extra-long knee upholstery 550 mm – 3D cover	276,-
		Bodypoint ankle huggers with upholstery	
	<input type="checkbox"/>	D44863 S (L: 170-200 mm)	213,-
	<input type="checkbox"/>	D44864 M (L: 190-230 mm)	213,-
	<input type="checkbox"/>	D44865 L (L: 220-290 mm)	223,-
	<input type="checkbox"/>	D44904 XL (L: 280-330 mm)	244,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300

Available from May 2022

D21745	<input type="checkbox"/>	Calf support bracket upholstery - pair - skai cover	49,-
D41617	<input type="checkbox"/>	Calf support upholstery extra soft - pair - EC or 3D cover depending on the fabric for the whole chair	190,-
D86900	<input type="checkbox"/>	Bodypoint calf strap M	162,-
D86901	<input type="checkbox"/>	Bodypoint calf strap L	162,-

ARM SUPPORTS AND ACCESSORIES FOR ARM SUPPORTS

Standard	<input checked="" type="checkbox"/>	Arm support swingable to the back for easy side transfer	0,-
Standard	<input checked="" type="checkbox"/>	Arm support hinged at the wheelchair back	0,-
Standard	<input checked="" type="checkbox"/>	Arm supports height and depth adjustable with tools	0,-
Standard	<input checked="" type="checkbox"/>	Arm pads 385 x 58 mm	0,-
D21077	<input type="checkbox"/>	Arm pads 420 x 80 mm	112,-
D21078	<input type="checkbox"/>	Arm pads 445 x 70 mm	112,-
D21079	<input type="checkbox"/>	Arm pads 525 x 80 mm	112,-
D41649	<input type="checkbox"/>	Standard arm pad upholstery extra soft, EC or 3D cover depending on the fabric for the whole chair (not to be combined with a tray)	99,-
D41650	<input type="checkbox"/>	Arm pad upholstery extra soft for D21077, EC or 3D cover depending on the fabric for the whole chair (not to be combined with a tray)	123,-
D22185	<input type="checkbox"/>	Arm support hemiplegic left (not to be combined with a tray) Please note: When large amounts of tilt are used, a conflict between arm support and main wheel are possible.	320,-
D22186	<input type="checkbox"/>	Arm support hemiplegic right (not to be combined with a tray) Please note: When large amounts of tilt are used, a conflict between arm support and main wheel are possible.	320,-
D23003	<input type="checkbox"/>	Hemiplegic cushion	196,-

MAIN AND FRONT WHEEL

Standard	<input checked="" type="checkbox"/>	Main wheels 24 x 1 PU black edition with drum brakes	0,-
Standard	<input checked="" type="checkbox"/>	Front wheels: 7" 175 x 45 PU with infill Seat height <input type="checkbox"/> 417 mm <input type="checkbox"/> 447 mm	0,-
D32035*	<input type="checkbox"/>	Main wheels 16" PU with drum brakes and frame extension Seat height <input type="checkbox"/> 417 mm <input type="checkbox"/> 447 mm	Price on request
D22009	<input type="checkbox"/>	Hand and spoke protectors for 24" wheels, transparent, pair	138,-
D22013	<input type="checkbox"/>	Push rim cover for 24" wheels, pair	112,-

HEAD SUPPORTS & ACCESSORIES FOR HEAD SUPPORTS



Relaxo

A

B

C

D

E


F

D30572 **Head support Relaxo, only available with 3D cover (value 454,-)**

0,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300

Available from May 2022

D30575	<input type="checkbox"/>	Head support A, with side support	305,-
D30566	<input type="checkbox"/>	Head support A with side support and head band	454,-
D82568	<input type="checkbox"/>	Extra comfort pads for the head band of D30566, pair	32,-
D30574	<input type="checkbox"/>	Head support B, small, oval	305,-
D30573	<input type="checkbox"/>	Head support C, large, square	305,-
 D32746	<input type="checkbox"/>	Netti Dynamic head support C – If combined with other options, please send us a request, as some solutions are only possible as special adaptations. For more information, feel free to contact us.	Price on request
D30576	<input type="checkbox"/>	Head support D, comfort (pressure distributive)	385,-
D37575	<input type="checkbox"/>	Head support E with adjustable side supports	465,-
D37576	<input type="checkbox"/>	Head support F with cheek support	611,-
D37579	<input type="checkbox"/>	Headrest E with shorter side support right	598,-
D37580	<input type="checkbox"/>	Headrest E with shorter side support left	598,-
D21027	<input type="checkbox"/>	Headrest adjustment tube (horizontal) + 100 mm, length 250 mm	90,-
D23004	<input type="checkbox"/>	Comfort cushion for head rest type C – removable with Mars cover	162,-
D23012	<input type="checkbox"/>	Easy Care cover for head rest type A, B, C – removable	79,-
D23013	<input type="checkbox"/>	3D cover for head rest type A, B, C – removable	79,-
D23015	<input type="checkbox"/>	Grey head rest cushion with Mars cover	214,-



D35372		[A] Savant head support incl. headband, small – height 130 mm User forehead 460-560 mm User ½ circumference behind ears 150-180 mm Inside distance along top support 300 mm, along middle support 180 mm (not approved for car transport – if this is needed please order an additional Netti headrest)	Price on request
D35374		[A] Savant head support incl. headband, adult – height 150 mm User forehead 560-610 mm User ½ circumference behind ears 200-300 mm Inside distance along top support 360 mm, along middle support 200 mm (not approved for car transport – if this is needed please order an additional Netti headrest)	Price on request
D45373	<input type="checkbox"/>	[B] Side pads for Savant head support. For extra lateral support. For the upper arms of the head support. Cannot be combined with the headband	Price on request
D45375	<input type="checkbox"/>	Headrest adjustment tube (horizontal) + 100 mm, length 250 mm	Price on request

TRAYS

D28826	<input type="checkbox"/>	Tray swingable	242,-
D28828	<input type="checkbox"/>	Tray lockable	242,-
D21052	<input type="checkbox"/>	Upholstery for tray	113,-

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300







Available from May 2022

VARIOUS ACCESSORIES






D21239	<input type="checkbox"/>	Tool kit	89,-
--------	--------------------------	----------	------

BELT OPTIONS FROM BODYPOINT AND NETTI

Belts and accessories for seat unit

				
[A]	[B]	[C]	[D]	[E]
Standard	<input checked="" type="checkbox"/>	Rail for mounting of pelvic belts	0,-	
D20845	<input type="checkbox"/>	Hip belt, with upholstery and car belt lock [A]	107,-	
		Bodypoint 2-point pelvic stabiliser Evoflex incl. mounting kit [B]	183,-	
D45021	<input type="checkbox"/>	XS (overall pad length: 380 mm)		
D45022	<input type="checkbox"/>	M (overall pad length: 460 mm)		
D45023	<input type="checkbox"/>	L (overall pad length: 620 mm)		
		Bodypoint 2-point hip belt incl. mouting kit [E]	213,-	
D45024	<input type="checkbox"/>	XS (overall pad length: 360 mm)		
D45025	<input type="checkbox"/>	S (overall pad length: 400 mm)		
D45026	<input type="checkbox"/>	M (overall pad length: 460 mm)		
D45040	<input type="checkbox"/>	Security cover for type [B], [C] and [E] in XS/M/L	70,-	
		Bodypoint 4-point pelvic stabiliser incl. mounting kit [C]	376,-	
D36889	<input type="checkbox"/>	XS (overall pad length 360 mm)		
D36893	<input type="checkbox"/>	M (overall pad length: 460 mm)		
D36894	<input type="checkbox"/>	L (overall pad length: 620 mm)		
D60001	<input type="checkbox"/>	Abduction belt [D]	497,-	
		Bodypoint supportive universal elastic strap	170,-	
D45084	<input type="checkbox"/>	S (Height 80 mm - 1420 mm long, easily cut-to-length)		
D45085	<input type="checkbox"/>	L (Height 110 mm - 1420 mm long, easily cut-to-length)		
D45086	<input type="checkbox"/>	XL (Height 150 mm - 1420 mm long, easily cut-to-length)		

Belts and accessories for back unit

				
[A]	[B]	[C]	[D]	[E]
D36890*	<input type="checkbox"/>	H-harness adaptor kit (recommended for all upper body harnesses) [A]	310,-	
D20881	<input type="checkbox"/>	H-belt, with upholstery and car belt lock [B]	227,-	
D41645	<input type="checkbox"/>	Belt cross harness [C]	175,-	
		Bodypoint Stayflex chest support without zipper [D]	246,-	
D89805	<input type="checkbox"/>	XS (Shoulder width 240-280 mm)		
D89898	<input type="checkbox"/>	S (Shoulder width 280-330 mm)		

NORWAY: Alu Rehab AS | Bedriftsvegen 23 | N-4353 Klepp stasjon | post@My-Netti.com | Phone +47 51 78 62 20
DENMARK: Alu Rehab ApS | Kløftehøj 8 | DK-8680 Ry | info@My-Netti.com | Phone +45 87887300

Available from May 2022

D86895	<input type="checkbox"/>	M (Shoulder width 330-410 mm)	
D86896	<input type="checkbox"/>	L (Shoulder width 410-480 mm)	
Bodypoint Pivofit shoulder harness [E]			246,-
D89802	<input type="checkbox"/>	XS (Shoulder width 240-280 mm)	
D89803	<input type="checkbox"/>	S (Shoulder width 280-330 mm)	
D86898	<input type="checkbox"/>	M (Shoulder width 330-410 mm)	
D86899	<input type="checkbox"/>	L (Shoulder width 410-480 mm)	

NETTI CUSTOMIZED

Please describe your current case and how a possible solution could look. If you have any picture or video material, we would also like to see this. To offer the best solution, please describe the user in the best possible way, e.g. such as any given diagnosis, user weight and/or the exact pattern movement of the user.

User weight:


Please inform us if there are any special circumstances in the living area that has an influence for the chair: e.g. the total possible width _____ mm and length _____ mm of the chair.

TECHNICAL DATA

Tilt (seat)	-3° to +30°
Recline (back)	87° to 133°
Total width of the wheelchair	Seat width + 232 mm
Total length	Min: 1280 mm with vertical leg support, vertical push bow and horizontal seat. Max. with 30° tilt is 1600 mm.
Distance arm support to seat plate	Min: 265 mm/ max. 377 mm
Weight of wheelchair	SW 450 – app. 53 kg incl. wheels, cushions, and head support
Frame colour	Pearl gold and black
Cover variants	Stretchable & flexible cover 3D Easy Care antibacterial cover, washable
Crash test	ISO 7176-19 (135 kg)

Available from May 2022

Leg supports	Min. calf length	Max. calf length
Manual adjustable leg rests	320 mm	680 mm
Universal leg rests	360 mm	500 mm
Netti Grandis leg rests	360 mm	640 mm
Netti Grandis fixed leg rests	380 mm	550 mm
Netti Dynamic leg rests	350 (450) mm	450 (550) mm



Measurements

- A Hip width/Seat width _____ mm
- B Thigh length / seat depth _____ mm
- C Arm rest height _____ mm
- D Leg length / seat height _____ mm
- E Backrest height _____ mm
- F Trunk width _____ mm