


FOAM WHEELCHAIR CUSHIONS

Reorder #	Cover	Cushion Size (inches)
□707052	Vinyl	18" x 16" x 6"
□753210	Nylon	16" x 18" x 2"
□753230	Vinyl	16" x 18" x 2"
□753240	No Cover	16" x 18" x 2"
□753250	Nylon	16" x 18" x 3"
□753270	Vinyl	16" x 18" x 3"
□753280	No Cover	16" x 18" x 3"
□753285	Nylon	16" x 18" x 4"
□753290	Vinyl	16" x 18" x 4"
□753295	No Cover	16" x 18" x 4"
□		

Maximum cushioning with minimum compression make Skil-Care High Density Foam Cushions your best choice for wheelchair and geri-chair comfort.



COVER CARE

Nylon and vinyl covers may be cleaned with a spray-and-wipe cleaner and damp cloth. Laundering the nylon cover will cause it to lose its water repellant quality. If the nylon cover is laundered, do not wash or dry it at temperatures over 180°F (82°C). Bleaching is not recommended.