

SafeSit Rollator

with forward-facing Seat

Patent
Registered
Design

Description

This ergonomically-design rollator allows user to be seated facing the same direction as when he is walking with the equipment. It allows tired user to be seated as quickly as possible while eliminating guesswork as the user can see the seat before seating.

The seat is designed as an attachment module for ease of handling.

Benefits

1. Reduce fall-related injuries among the elderly while enhancing their mobility as the user will have a firm grip of the rollator before he sits down.
2. Promote active and independent lifestyle among the elderly as he no longer needs to look for seat when he uses the rollator during a stroll.

