



Scoop Stretcher YDC-4B



Scan for product details



Technical features:

- * It's a separable-type emergency stretcher for ambulance to transfer fractured & serious patients.
- * At both ends of the stretcher, a hinged separate and joined arrangement is set in the middle, which can separate the stretcher into two parts,
- * The stretcher could strengthen the patients at origin place which could reduce the further hurt to the patients.
- * The patient can be quickly & easily put on the stretcher and also the stretcher can be taken out from the back of patient without moving the patient.
- * The length of the stretcher can be adjusted according to patient's body.
- * The foot-end is designed as a narrow-frame structure.
- * It is made of high-strength aluminum alloy materials.
- * It characterized by its being light-weight, small-sized, easy-carried, use-safely and easy for sterilization and clean.
- * It's mainly used for hospitals, sports, ambulance carrying patients and wounded person.

Model	Unfolded Size (L×W×H)	Folded Size (L×W×H)	Self-weight	Load Bearing	Packing		
					Quantity	Size (L×W×H)	G.W
YDC-4B	214×42×7cm	167×42×7cm	8Kg	160Kg	1	170×43×8cm	9Kg