


Stretcher 4-Fold Aluminum YDC-1A10



Scan for product details



Technical features:

- * The stretcher with outrigger, it can be folded contracted by transverse and longitudinal, and with handbag.
- * It is made of high-strength aluminum alloy and oxford leather materials.
- * It characterized by its being light-weight, small-sized, easy-carried, use-safely and easy for sterilization.
- * It's mainly used for hospitals, sports, ambulance, battle fields and outdoors carrying patients and wounded person.



Model	Unfolded Size (L×W×H)	Folded Size (L×W×H)	Self-weight	Load Bearing	Packing		
					Quantity	Size (L×W×H)	G.W
YDC-1A10	210×53×12cm	61×30×12cm	6Kg	159Kg	1	63×25×18cm	7Kg