



LIFE IS TIME™
 When minutes count
EMERGENCY OXYGEN
 can be the difference between
LIFE OR DEATH



*This unit
 could save
 your life.*

*Only the very best and complete
EMERGENCY OXYGEN UNIT
 provides supplemental oxygen
 for First Aid in **TWO WAYS:***

- **INHALATOR**
 ... for a breathing victim
 and
- **CPR RESUSCITATOR adjunct**
 ... for a non-breathing victim



LIFE® OxygenPac features a LIFE® mask for "mouth-to-mask" resuscitation with a one-way valve. This eliminates direct contact with the victim's mouth and nose while resuscitating a non-breathing person. Maximum hygiene is provided by the one-way valve, which eliminates transmittal of fluids and contagious diseases. LIFE® mask fits both adult and child.
 (Patented)

LIFE® OxygenPac MEETS OR EXCEEDS ALL GOVERNMENT AND FDA REQUIREMENTS: A NON-PRESCRIPTION DEVICE INTENDED FOR EMERGENCY USE.

www.LIFECorporation.com

5 YEAR WARRANTY

LIFE® is a Registered Trademark of LIFE Corporation

Patents Approved and Pending

LIFE® OxygenPac

should be part of your first aid program.

- **SIMPLE and EASY TO USE**, with external Two-Step Instructions and Clear Cover
- **Wall Mounted for Recognition and Accessibility** always available on display (20" x 6 3/8" x 9")
- **Easily Portable with Convenient Handle** (16 lbs.)
- **Simple "On-Off" Control Lever**
- **Constant Reading Supply Gauge**
 90 Minute Supply @ 6 LPM (566 liters O₂)
- **Resuscitation LIFE® mask with One-Way Valve**
 Fits both adult and child
- **REFILLABLE**, low cost, by any compressed gas distributor with universal CGA #870 fitting
- **Tight Seal Protects Unit from Dust and Moisture**
- **Units come Ready to Use - Full of oxygen** with wall brackets and screws

MODEL FEATURES
 LIFE-6FF 6 LPM (Liters Per Minute) Fixed Flow

LIFE-612 6 & 12 LPM - "**NORM**" & "**HIGH**" (2 settings)
 provides up to recommended 100% inspired O₂

LIFE-025 0 to 25 LPM Variable Flow (for EMT's)
 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25

Made in the U.S.A. ©2005 LIFE Corporation, Milwaukee

*Incontestable trademarks of LIFE Corporation

**IDEAL LOCATIONS FOR THE LIFE® OxygenPac
... WHEREVER THERE ARE PEOPLE!**



INDUSTRIAL/FACTORY



PUBLIC FACILITIES



ENTERTAINMENT/RECREATION



OFFICES

LIFE® OxygenPac is attractive. Prominently displayed, it makes the simple but strong statement, "This Organization Cares For Its Employees And Visitors."



**American Medical Association
GUIDELINES ON OXYGEN EQUIPMENT IN CPR ...**
as published by the American Heart Association

- "... supplemental oxygen should be used as soon as it is available ..."
- "... oxygen should not be withheld for fear of suppressing respiration ..."
- "The highest possible oxygen concentration (preferable 100%) should be administered as soon as possible to all patients with cardiac or pulmonary arrest or other patients with suspected hypoxemia regardless of cause."
- "Masks equipped with one-way valve ... are recommended for mouth to mask ventilation."

**EMERGENCY OXYGEN
INSTRUCTIONS**

1. Move lever → ON
2. Place mask over nose and mouth



IF BREATHING HAS STOPPED, START MOUTH-TO-MASK RESUSCITATION (BLOW UNTIL CHEST RISES, EVERY 5 SECONDS)

612 MODEL LIFE® OxygenPac

- SAME SUPERIOR FEATURES of standard LIFE® OxygenPac
- PLUS: NORM & HIGH FLOW RATES
Regulator with "Click-Stop" two settings
 - "NORM" - 6 LITERS PER MINUTE for GENERAL emergency situations delivers approximately 45% inspiration concentration to the average adult breathing victim as an INHALATOR
 - "HIGH" - 12 LITERS PER MINUTE for EXTREME emergency situations delivers oxygen concentrations approaching the recommended 100% inspiration concentration as an INHALATOR for breathing victims and as a RESUSCITATOR for oxygen supplemented mouth-to-mask resuscitation for a non-breathing victim.
- SIMPLE OPERATION IS EASY-TO-USE
LIFE® OxygenPac 612 MODEL cannot be turned to below 6 LPM (Lower flow rates are not advised when used with mask delivery.)

5 YEAR WARRANTY

Our warranty includes provisions that all parts and materials are warranted against failure due to defective material or poor workmanship for FIVE YEARS from date of purchase.

LIFE® OxygenPac TECHNICAL SPECIFICATIONS

<p>REGULATOR</p> <p>Flow length.....90 minutes ± U.S.P oxygen Liter flow.....6 LPM constant Style.....Piston</p> <p>OXYGEN PAC UNIT</p> <p>H x W x D.....20" x 6 3/8" x 9" Weight.....16 lbs.</p> <p>Component material:</p> <p>Case.....Polyethylene Viewing cover.....Lexan® polycarbonate Tight seal protects unit from dust and moisture</p> <p>GAUGE</p> <p>Readout.....Constant Cylinder contents "Full" indicator.....Green wedge Maximum pressure rating...4000 PSI</p> <p>SHUT-OFF VALVE</p> <p>"On-Off" lever.....Adjustable ratchet Thread.....1/2" - 14 NPT taper Safety burst relief valve.....3025 PSI Refill adaptor.....C.G.A. #870</p>	<p>CYLINDER (included)</p> <p>Fill capacity.....566 liters "Full" pressure.....2015 PSI Test pressure.....3358 PSI Material.....REFILLABLE seamless steel</p> <p>Dimensions:</p> <p>Outside diameter.....5 1/4" Overall length.....14" Minimum wall thickness.....0.108"</p> <p>MASK</p> <p>Style.....Mouth-to-mask resuscitator with one-way flow valve Material.....Clear vinyl</p> <p>TUBING</p> <p>Length.....7 ft. Clamp-proof star diameter</p> <p>Brackets and Fasteners for Wall Mounting Included with Every Unit</p> <p>REFILLABLE by Any Medical Industrial Gas Distributor with Universally Standard C.G.A. #870 Fitting</p>
--	--

ORDERING INFORMATION

MODEL #	FEATURES
LIFE-6FF	6 LPM (Liters Per Minute) Fixed Flow
LIFE-612	6 & 12 LPM - "NORM" & "HIGH" (2 settings) provides up to recommended 100% inspired O2
LIFE-025	0 to 25 LPM Variable Flow (for EMT's) 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25
LIFE-101	Spare Cylinder Assembly (cylinder, valve, gauge, O2)
LIFE-102	LIFE® mask & 1-way valve Fits both adult and child (mask, valve, strap, O2 npl)
LIFE-BAG	LIFE® Bag/mask manual Resuscitator with O2 inlet

(All units complete with Cylinder, Oxygen, Brackets and Fasteners for Wall Mounting)

Authorized Distributor:

Uncompromisingly manufactured with quality materials, technology and workmanship



*Incontestable trademarks of LIFE Corporation
LIFE® is a Registered Trademark of LIFE Corporation



PATENTS APPROVED AND PENDING Made in the U.S.A.

© 2005 LIFE Corporation, Milwaukee

www.LIFECorporation.com