

Anti-Thrust Cushion



Reorder #	Description	Size
757110	Anti-Thrust Cushion, Foam-Only, Soft Base	18" x 16"
757115	Anti-Thrust Cushion, w/Gel Pad, Soft Base	18" x 16"
757120	Anti-Thrust Cushion, Foam-Only, Firm Base	18" x 16"
757125	Anti-Thrust Cushion, w/Gel Pad, Firm Base	18" x 16"

PURPOSE

Provides restraint-free positioning for wheelchair users with forward thrust and sacral sitting problems. Maintains pelvic neutrality. Designed as a comfortable alternative to pommel-style cushions, especially for males.

Anti-Thrust barrier diminishes problems associated with forward thrust and sacral seating while maintaining pelvic neutrality.

Safety straps secure cushion to wheelchair.

Available gel pad offers maximum pressure relief in the critical coccyx area.

Low-shear cover prevents skin-damaging friction. This incontinent-proof cover removes for easy laundering.

High-density foam provides comfort and pressure relief

Contoured base eliminates sling-seat hammocking and minimizes seat rise. Choose soft or firm base.

