

# DEBUG TILT IN SPACE RECLINING BEACH WHEELCHAIR



# DEBUG TILT-IN-SPACE RECLINING BEACH WHEELCHAIR

debug

Two custom features combined into one ultimate beach wheelchair. Easily switch from a Tilt-in-Space chair to a Reclining chair in seconds without tools. With so much flexibility in seating positions the Debug Tilt-in-Space beach wheelchair keeps you comfortable and properly seated. You will be ready for a fun day out on the beach, on a nature trail or just for a relaxing day in your own backyard.

LONG LASTING #316L STAINLESS STEEL FRAME / REMOVE-ABLE ARTICULATING TELESCOPIC LEG REST / FITS THROUGH A STANDARD ADA DOORWAY / LOW MAINTENANCE DEL-RIN BUSHINGS / POLYURETHANE PNEUMATIC WHEELS ROLL ACROSS SAND, SNOW AND GRAVEL / DUAL STITCHED SLING WITH TENERA 130 THREAD / PHIFERTEX PLUS SLING MATERIAL RESISTS STRETCHING, MILDEW AND BACTERIA, COLD CRACK AND ULTRAVIOLET DAMAGE AND IS MACHINE WASHABLE.

## SPECIFICATIONS

SEAT SIZES	ENGLISH (INCH)	METRIC (CM)
Widths	14, 18, 22	36, 46, 56
Average depth	17	43
Seat-to-floor height	17 – 21	43 – 53
Backrest height	18	46
Elevating leg rest <i>knee to heel</i>	15 – 22	38 – 56
Telescope footrest <i>knee to heel</i>	15 – 22	38 – 56

## WEIGHTS

Capacity	250 – 400 lbs	114 – 182 kg
Chair weight	40 – 65 lbs	18 – 30 kg

\* Shown with optional headrest and custom two piece cushion.

Manufactured in the USA,  
Debug has a 10 year warranty on stainless steel frames and a 30 day money back guarantee.

© Copyright Deming Designs Inc., 2001–2016. All rights reserved.

ORDER ONLINE [WWW.DEBUGMOBILITYPRODUCTS.COM](http://WWW.DEBUGMOBILITYPRODUCTS.COM)