

Name:	
School/Organization:	
Dealer/Rep:	

Phone:			
E-mail:			
Sales Rep:		Quote:	

PARALLEL BARS

ELECTRIC PLATFORM MOUNT

STANDARD FEATURES:

- Load capacity: 500 lbs
- 15-29"W x 84" or 120"L x 27-52"H
- Platform: 43"W x 120"L x 15-29"W between columns
- Stainless steel hand rails (1.5" Diameter)
- Digital display on conveniently located, user-friendly control panel shows exact bar height
- Memory buttons allow storage and recall of four separate bar height settings
- Microprocessor control unit ensures synchronous operation of all four columns
- Height monitoring functions ensure reliable operation and safe switch-off in case of a fall
- Slip-resistant platform
- ETL Approved

STANDARD CONFIGURATIONS:

☐ AM-TM702 - 7' Electric Platform Mount Parallel Bars 15-29"W x 84"L x 27-52"H	☐ AM-TM700 - 10' Electric Platform Mount Parallel Bars 15-29"W x 120"L x 27-52"H	☐ AM-TM703 - 15' 2PC Electric Platform Mount Parallel Bars 15-29"W x 180"L x 27-52"H
QTY:	QTY:	QTY:

OPTIONS:

☐ AM-714 - Abduction Board for Parallel Bars .75"W x 94"L x 5"H
QTY:
☐ AM-715 - Balance Beam for Parallel Bars 5"W x 80"L x 2.5"H
QTY:
☐ AM-716 - Foot Placement Ladder for Parallel Bars 24"W x 84"L x 4"H
QTY:



800.701.5122
armedicamfg.com