

# KARMA FLEXX HD

Order Form 2024

**flexx**

*Versatile design for Heavy Duty users*



### Standard Features

- 170 kg Maximum Weight Capacity
- Self Propulsion with 22" or 24" Q-Release Rear Wheels
- Transit Propulsion with 14" Fix Rear Wheel
- 6", 7" or 8" Caster Wheels
- Seat Widths: 505/ 555/ 605mm
- Seat Depths: 405/ 455/ 505mm
- Angle Adjustable Backrest: -4°/ 0°/ 4°/ 8°/ 12°/ 16°
- Backrest Heights: 430mm/ 510mm
- 80° Swing in/Swing out Legrest with Flip-up Alu. Footplates
- Push to Lock Brakes with Brake Extension

**Base Price: S\$2000**

◆ Standard feature Free of Charge

### Seat

#### Seat Width and Seat Depth

Width \ Depth	505mm(20")	555mm(22")	605mm(24")
405mm(16")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
455mm(18")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
505mm(20")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

#### Seat-to-Floor Height

Front Caster Rear Wheel	6"x1 1/4"	7"x1 1/4"	8"x1 1/4"	8"x2" (200x50)	Castor Positions
14"	<input type="checkbox"/> 18"/ L1 <input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2	<input type="checkbox"/> 20"/ L2	<input type="checkbox"/> 20"/ W2	<b>L= Long Fork</b> <b>H= for 21" seat height</b> <b>W=Wide Fork for 20" Seat Height</b> 
22"	<input type="checkbox"/> 18"/ L1 <input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2 <input type="checkbox"/> 20"/ L3	<input type="checkbox"/> 20"/ L2	<input type="checkbox"/> 20"/ W2	
24"	<input type="checkbox"/> 19"/ L3	<input type="checkbox"/> 19"/ L2 <input type="checkbox"/> 20"/ L3	<input type="checkbox"/> 20"/ L2 <input type="checkbox"/> 21"/ H2	<input type="checkbox"/> 20"/ W2	

#### Seat Upholstery

Padded Seat Upholstery

<b>Backrest</b>	
<b>Backrest</b>	
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 430mm Back Height Fixed Backrest	◆
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 510mm Back Height Fixed Backrest	○
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 430mm Back Height Foldable Backrest	\$100
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Delivered 510mm Back Height Foldable Backrest	\$100
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) Height Adjust 380mm to 460mm Fixed Backrest	\$130
<b>Backrest Upholstery</b>	
<input type="checkbox"/> Without Backrest Upholstery	○
<input type="checkbox"/> Padded Backrest Upholstery	◆
<input type="checkbox"/> Backrest Upholstery for Height Adjustable Backrest	\$100
<input type="checkbox"/> Strap Tension Adjustable Backrest (Including Detachable Padded Cover)	\$130
<b>Legrests/Footplates</b>	
<b>Legrest RH</b>	
<input type="checkbox"/> 80° RH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) RH Swing in/Swing out	\$130
<input type="checkbox"/> Amputee Legrest RH	\$100
<b>Legrest LH</b>	
<input type="checkbox"/> 80° LH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) LH Swing in/Swing out	\$130
<input type="checkbox"/> Amputee Legrest LH	\$100
<b>Footplate RH (Aluminum Footplate is recommended while 80° Swing in/Swing out Legrest is ordered)</b>	
<input type="checkbox"/> Flip-up RH Aluminum Footplate	◆
<input type="checkbox"/> Flip-up RH Tube-in-center Footplate	○
<input type="checkbox"/> Flip-up RH Angle Adjustable Footplate	\$25
<b>Footplate LH (Aluminum Footplate is recommended while 80° Swing in/Swing out Legrest is ordered)</b>	
<input type="checkbox"/> Flip-up LH Aluminum Footplate	◆
<input type="checkbox"/> Flip-up LH Tube-in-center Footplate	○
<input type="checkbox"/> Flip-up LH Angle Adjustable Footplate	\$25
<b>Armrests</b>	
<b>Armrest RH</b>	
<input type="checkbox"/> Flip-Back Ergo Armpad RH Adj. in Height	◆
<input type="checkbox"/> Flip-Back Boat Armpad RH Adj. in Height	\$15
<input type="checkbox"/> Flip-Back Hemi Armpad RH Adj. in Height	\$130
<b>Armrest LH</b>	
<input type="checkbox"/> Flip-Back Ergo Armpad LH Adj. in Height	◆
<input type="checkbox"/> Flip-Back Boat Armpad LH Adj. in Height	\$15
<input type="checkbox"/> Flip-Back Hemi Armpad LH Adj. in Height	\$160
<b>Front Castors/Forks</b>	
<b>Castors</b>	
<input type="checkbox"/> 150mm x 30mm Solid Tire (6")	○
<input type="checkbox"/> 175mm x 30mm Solid Tire (7")	◆
<input type="checkbox"/> 200mm x 30mm Solid Tire (8")	○
<input type="checkbox"/> 200mm x 50mm Solid Tire (8")	\$30
<b>Forks</b>	
<input type="checkbox"/> Stylish Fork	◆
<input type="checkbox"/> 21" Seat Height Long Fork	○
<input type="checkbox"/> 2" Wide Long Fork (While 200mm x 50mm Tire is ordered)	\$15

<b>Rear Wheels</b>	
<b>Rear Wheel</b>	
<input type="checkbox"/> 14" Solid Tire Plastic Rim	○
<input type="checkbox"/> 22" Pneumatic Tire Spoked Rim with Plastic Handrim	\$25
<input type="checkbox"/> 22" Pneumatic Tire Spoked Rim with Aluminum Handrim	\$70
<input type="checkbox"/> 24" Solid Tire Spoked Rim with Aluminum Handrim	◆
<input type="checkbox"/> 24" Pneumatic Tire Spoked Rim with Ergo Handrim	\$25
<input type="checkbox"/> 24" Pneumatic Tire Spoked Rim with Aluminum Handrim	\$70
<b>Rear Wheels Axle</b>	
<input type="checkbox"/> Fixed Axle (Only available for 14" rear wheel)	○
<input type="checkbox"/> Quick Release Axle (Karma version)	◆
<b>Brakes</b>	
<input type="checkbox"/> Push to Lock Brakes with Brake Extension	◆
<input type="checkbox"/> Attendant Brakes	\$70
<b>Safety Options</b>	
<input type="checkbox"/> Swing-away Height Adjustable Anti-tippers (Pair)	\$130
<b>Accessories</b>	
<input type="checkbox"/> Lap Belt (Safety Belt)	\$40
<input type="checkbox"/> Cane Holder	\$180
<b>Price</b>	
	<b>Total Price</b>