Stationary Lift Systems

HeliQ

Discreet and contemporary design

The Human Care HeliQ stationary overhead lift features a modern, discreet look that easily fits into any environment. Its small size ensures ease of use in either institutional or home care applications. The lift comes in three different models with three different weight capacities, 150 kg, 220 kg and 300 kg. It also allows optimal patient comfort with 2-point or 4-point hanger bars in different lengths to suit virtually any patient.

Lift intelligence

The HeliQ has the ability to store usage and performance data. With our data and service application and a PC or tablet device, it lets you find out more about the usage and the condition of your lift. Detailed information can be viewed directly in the application or exported to other media for additional analysis.

Maximum comfort and lifting efficiency

Constant lifting speed regardless of the patient's weight plus soft-start and stop makes lifting a comfortable experience. The Human Care HeliQ offers twospeed lifting, high battery capacity and convenient battery charging through the hand control. The lift can

also be equipped with power traverse (motorized movement) and in-rail charging (continuous charge). The very low weight of the Human Care HeliQ makes it easy to handle during installation and service.

Power traverse

The HeliQ lifts have an option to be connected to a power traverse, meaning motorized movement sideways. The HeliQ Power Traverse has dual speed, high water protection (IPX4), a 300 kg (661 lbs) safe-working load capacity, including the possibility of passing curtain gaps in rails.

70" /min -

118"/ min

according to ISO

10535:2006

ETL-listed

Description							Part no.	Safe wor	Safe working load	
HeliQ 150 kg incl. lift, charger, charging dock (rail roller and hanger bar, ordered separately) HeliQ 220 kg incl. lift, charger, charging dock (rail roller and hanger bar, ordered separately) HeliQ 300 kg incl. lift, charger, charging dock (rail roller and hanger bar, ordered separately)							55100-XX 55110-XX 55120-XX	150,220, 300kg (330, 484,661lbs)		7kg (15.4lbs)
HeliQ 150 kg with power traverse incl. lift, charger, charging dock, rail roller (hanger bar, ordered separately) HeliQ 220 kg with power traverse incl. lift, charger, charging dock, rail roller (hanger bar, ordered separately) HeliQ 300 kg with power traverse incl. lift, charger, charging dock, rail roller (hanger bar, ordered separately)							55130-XX 55140-XX 55150-XX			
Specify region when ordering, XX=EU(Europe), AU(Australia), GB (Great Britain), NA (North America), JP (Japan)										
Dimensions	Battery	Battery charger	Emergency lowering	Hanger bar option	IP-class	Lifting distance	Lifting speed			
(WxHxL) 235x235x140mm	2,3 Ah 2x12V	Input 100-240 V, AC 50-60 Hz	Electronic and mechanical	2-point, 4-point.	IPX4 (lift) IPX5 (hand	2080mm (81,9")	n 1.8m/m 3m/min	1		

Ordered

separately.

control)

For hanger bars and accessories, please see page 24

Output 27,8V,

0,8A DC

(9.2x9.2x5.5")

HUMAN

A stationary lift system is the optimal solution for hospitals and surgical care facilities.