




Featured Products

Improving your life everyday™





Established in 1976, K Care Healthcare Solutions is an Australian owned, family business, that manufactures, distributes and supplies over 2500 types of products to Hospital, Residential Aged Care and Community Care for both Australian and International markets.

Our team is made up of some of Australia's most knowledgeable and experienced healthcare industry trade professionals, offering excellent support right across the country. We have a proven track record for delivering equipment solutions within tight time frames to DVA, Hospitals, Occupational Therapists, Healthcare Centres and other service providers.

With values based on integrity, respect, leadership, passion and innovation, you can have confidence knowing that you are dealing with a company that has customer satisfaction at the core of everything we do.

Working closely with Healthcare Professionals, K Care Healthcare Solutions has developed trusted brands that provide an extensive range of products that meet and exceed requirements demanded in healthcare environments today. Our products have a reputation for being aesthetically pleasing in design, highly functional and long-lasting.

K Care Healthcare Solutions product range feature as a trusted brand in many Homes, Hospitals and Aged Care Facilities throughout Australia.



Contents

Bathroom and Toilet	1
Seating and Lift Chairs	29
Lifters and Slings	47
Clinical Equipment	53
Daily Living Aids	63
Mobility Range	75
Wheelchair Range	83
Beds and Mattresses	93

Note: All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information. For regular product and specifications updates please visit www.kcare.com.au

Bathroom and Toilet



35° of tilt using twin gas struts for smooth and easy tilt

Lateral thoracic supports and knee guides

Swing away footrests with padded calf supports

Conforms to AS/NZS 3973:2009 for safety and peace of mind

Tilt in Space

Rehab Shower Commode

KA11507-REHAB

KA11512-REHAB

Features

- Side access standard on 445 and 500mm seat widths
- 10-35° tilt angle
- Lateral thoracic support mounting plates
- Hip/thigh guide support mounting plates
- Height adjustable frame
- Tension adjustable backrest
- Integrated push-handle for manoeuvrability
- Height adjustable flip up arms for unrestricted side transfer
- Swingaway and removable height adjustable footrest with contoured padded calf supports
- 4x 125mm total locking castors
- Designed to be used with a K Care patented removable seats
- Available in 400, 445, 500mm seat widths
- Padded seat not included

Specifications

KA11507-REHAB

Seat Width	445mm
Seat Depth	435mm
Backrest Height	855mm
Overall Width	580mm
Overall Depth	1070mm
Floor to Seat Frame Height Range	465mm - 540mm
Max User Weight	200kg
Brand	K Care

KA11512-REHAB

Seat Width	500mm
Seat Depth	435mm
Backrest Height	855mm
Overall Width	635mm
Overall Depth	1070mm
Floor to Seat Frame Height Range	465mm - 540mm
Max User Weight	200kg
Brand	K Care

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind

Removable swing-away height adjustable footrests

3 strap PVC backrest for comfort and ease of access

Durable and long lasting with a stainless steel frame, impervious PVC straps and water resistant castors

Self Propelled Rear Wheels

Mobile Shower Commode

KA120SR

Features

- Swing-up arms for unrestricted side transfer
- Large self-propelling rear wheels are flush with seat top for unrestricted side transfer
- Three axle position holes for flexible rear wheel set-up
- Removable swing-away height adjustable footrests for improved comfort and easy access
- 2 x 125mm Locking swivel castors on front
- Range of padded seating available using K Care patented seating system
- Available in 320, 400,445, 500, 550, 600 and 650mm seat widths
- Padded seat not included

Specifications

KA120SR

Overall Width	630mm
Overall Depth	1050mm
Seat Width	445mm
Seat Depth	435mm
Min Seat Height	452mm
Max Seat Height	550mm
Backrest Height	415mm
Max User Weight	200kg
Brand	K Care

Integrated
push-handle for
manoeuvrability

Swingaway and
removable height
adjustable footrest

4x 125mm
locking
castors

Independently tested to conform to
AS/NZS 3973:2009 for safety and peace of mind

Paediatric Attendant Propelled Rehab Commode

KAT14SC-REHAB

Features

- Side access standard on 445 and 500mm seat widths
- Lateral thoracic support mounting plates
- Hip/thigh guide support mounting plates
- Height adjustable frame
- Tension adjustable backrest
- Integrated push-handle for manoeuvrability
- Height adjustable flip up arms for unrestricted side transfer
- Swingaway and removable height adjustable footrest
- 4x 125mm total locking castors
- Designed to be used with a K Care patented removable seats
- Available in 320, 400, 445, 500mm seat widths
- Padded seat not included
- Self Propelled model also available


Specifications

KAT14SC-REHAB

Overall Width	485mm
Overall Depth	750mm
Seat Width	320mm
Seat Depth	320mm
Min Floor to Seat Height	452mm
Max Floor to Seat Height	575mm
Armrest Type	Loop
Backrest Height	415mm
Footrest Type	Swing Away
Product Weight	12kg
Max User Weight	120kg
Brand	K Care



IPX6 rated powered electric actuator tilt mechanism



3x 125mm twin locking castors with 1x 125mm directional locking castor

Adjustable height swing-away foot rests with padded calf rests

Powered Tilt In Space shower commode for larger users

Electric Tilt In Space Mobile Shower Commode

KA11506

KA11506S

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Dual electric actuator tilt mechanism
- Contoured adjustable clamp on headrest
- Curved covered padded backrest
- Integrated push-handle for manoeuvrability
- Curved removable padded armrests
- Swingaway and removable height adjustable footrest with contoured padded calf supports
- 3x 125mm total locking twin castors 1x 125mm directional locking twin castor
- Designed to be used with a K Care patented removable seats
- Available in 600 and 650mm seat widths
- Padded seat not included

Specifications

KA11506

Seat width	600mm
Seat depth	470mm
Floor to top of seat frame	550mm
Backrest height	630mm
Tilt angle	10-35 degrees
Overall width	710mm
Overall depth	1140mm
Overall height	1250mm
Max user weight	300kg

KA11506S

Seat width	650mm
Seat depth	480mm
Floor to top of seat frame	550mm
Backrest height	630mm
Tilt angle	10-35 degrees
Overall width	760mm
Overall depth	1140mm
Overall height	1250mm
Max user weight	300kg
Brand	K Care

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind



Gas strut operated tilt mechanism for smooth easy tilting



Range of padded seating available using K Care patented seating system

Attendant Tilt In Space

Shower Commode

KA11507

KA11507-HA

KA11512

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Extra high PVC strap backrest for improved comfort and safety
- 10–35° tilt angle
- Integrated push-handle for manoeuvrability
- Swing-up padded arms for unrestricted side transfer
- Swingaway and removable height adjustable footrest with contoured padded calf supports
- 4x 125mm total locking castors
- Designed to be used with a K Care patented removable seats
- Available in 400, 445mm and 500mm seat widths
- Padded seat not included

Specifications

Overall Depth	1070 mm
Overall Height	1453 mm
Floor To Seat Height	540 mm
Armrest Type	Loop
Armrest Height	365 mm
Backrest Height	855 mm
Footrest Type	Swing Away
Max User Weight	200 kg
Brand	K Care

KA11507

Seat Width	445 mm
Seat Depth	430 mm
Product Weight	21 kg

KA11512

Seat Width	500 mm
Seat Depth	435 mm
Product Weight	22 kg

Electric tilt in space



Side frame access



Latest battery and controller technology



Powered electric tilt for ease of use

**NEW
PRODUCT**

Electric Tilt in Space

KA11507-E

Shower Commode

Features

- Electric tilt in a Standard Width Commode
- New Linak Lift 50 Lithium Ion Battery
 - light weight, fast charge
- Max User Weight 200kg
- Locally manufactured in Adelaide, can be customised
- Based on the incredibly popular KA11507 with the latest battery and controller technology
- Padded seat not included

Specifications

KA11507-E

Overall Width	580mm
Overall Depth	1070mm
Seat Width	445mm
Seat Depth	435mm
Floor to Seat Frame Top	540mm
Armrest Type	Loop
Backrest Height	875mm
Footrest Type	Swing Away
Max User Weight	200kg
Brand	K Care

30° of tilt



Two armrest configurations



Curved vinyl covered padded backrest



Easy operation tilt & recline
with active access

**NEW
PRODUCT**

Tilt in Space/Recline Shower Commode

KA11507-R

Features

- Max User Weight 200kg
- 30 Degrees of Tilt
- 10 – 30 Degrees of recline
- Australian Designed, Made in Adelaide
- Two armrest configurations, articulating or fixed
- Can be Customised
- Padded seat not included
- Includes all standard features from the popular KA11507

Specifications

KA11507-R

Overall Width	580mm
Overall Depth	1070mm
Seat Width	445mm
Seat Depth	435mm
Floor to Seat Frame Top	540mm
Armrest Type:	Loop or Swing-up Dynamic
Backrest Height:	875mm
Footrest Type	Swing Away
Max User Weight	200kg
Brand	K Care


Designed for strength and enhanced cleaning



Rear access



Swing-up arms for unrestricted side transfer



Range of padded seating available using K Care patented seating system

NEW
PRODUCT

Rear Access Tilt in Space

KA11507-RA

Shower Commode

Features

- Designed for strength and enhanced cleaning
- Starting point for custom commodes with rear access
- Padded seat not included
- **Not compatible with Side access**

Specifications

KA11507-RA

Overall Width	580mm
Overall Depth	1070mm
Seat Width	445mm
Seat Depth	435mm
Floor to Seat Frame Top	540mm
Armrest Type	Loop
Backrest Height	875mm
Footrest Type	Swing Away
Max User Weight	200kg
Brand	KCare


Designed for strength
and enhanced cleaning



Rear access



Swing-up arms for
unrestricted side transfer



Range of padded seating available
using K Care patented seating system

**NEW
PRODUCT**

Self Propelled Rear Access

Mobile Shower Commode

KA120SR-RA

Features

- Designed for strength & enhanced cleaning
- Starting point for custom commodes with rear access
- Padded seat not included
- **Not compatible with Side access**

Specifications

KA120SR-RA

Overall Width	625mm
Overall Depth	1050mm
Seat Width	445mm
Seat Depth	435mm
Floor to Seat Frame Top	540mm
Armrest Type	Loop
Backrest Height	335mm
Footrest Type	Swing Away
Max User Weight	200kg
Brand	KCare

Integrated
push-handle for
manoeuvrability



Slide out/retractable
weight bearing footplate



3 strap back rest for
comfort and easy access



Integrated push-handle for manoeuvrability

Weight Bearing Footplate

Self Propel Mobile Shower Commode

KA123S

KA123S50

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- 3 strap back rest for comfort and easy access
- Swing-up arms for unrestricted side transfer
- Slide out/retractable weight bearing footplate
- 2x 125mm total locking castors
- 2x 22" self-propel wheels
- Designed to be used with a K Care patented removable seats
- Available in 445 and 500mm seat widths
- Push rims available as optional extra
- Padded seat not included

Specifications

Overall Depth	930mm
Overall Height	1005mm
Seat Depth	430mm
Floor To Seat Height	550mm
Backrest Height	400mm
Footrest Type	Weight Bearing
Product Weight	15kg
Max User Weight	200kg
Brand	K Care

KA123S

Overall Width	630mm
Seat Width	445mm
Armrest Type	Loop

KA123S50

Overall Width	685mm
Seat Width	500mm
Armrest Type	Padded Loop

Curved vinyl covered padded backrest



Removable armrests with arm pads

Swingaway and removable height adjustable footrest

Designed to be used with a K Care patented removable seats

Swingaway Footrests

Transit Mobile Shower Commode

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- Curved vinyl covered padded backrest
- Integrated push-handle for manoeuvrability
- Removable armrests with arm pads
- Swingaway and removable height adjustable footrest
- 4x 125mm total locking castors
- Available in 320, 400,445, 500, 550, 600 and 650mm seat widths
- Padded seat not included

KA114SC

KA114S

KA114S50

KA114S55

KA114S60

KA114S65

Specifications


KA114S65

Overall Width	765mm
Overall Depth	930mm
Overall Height	1005mm
Seat Width	650mm
Seat Depth	480mm
Floor To Seat Height	550mm
Armrest Type	Padded Removable
Armrest Height	315mm
Backrest Height	405mm
Footrest Type	Swing Away
Product Weight	14kg
Max User Weight	300kg
Brand	K Care

Durable and long lasting with a stainless steel frame, impervious PVC straps and water resistant castors



Footrest is weight bearing and retractable



Swing-up arms for unrestricted side transfer

Range of padded seating available using K Care patented seating system

Attendant Propelled
Weight Bearing Platform

KA113S

KA113S50

Mobile Shower Commode

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- 3 strap back rest for comfort and easy access
- Integrated push-handle for manoeuvrability
- Swing-up arms for unrestricted side transfer
- Slide out/retractable weight bearing platform
- 4x 125mm total locking castors
- Designed to be used with a K Care patented removable seats
- Available in 445, and 500mm seat widths
- Padded seat not included

Specifications

KA113S

Overall Depth	830mm
Overall Height	1005mm
Overall Width	520mm
Seat Depth	430mm
Seat Height	590mm
Seat Width	445mm
Backrest Height	400mm
Max User Weight	200kg

KA113S50

Overall Depth	830mm
Overall Height	1005mm
Overall Width	630mm
Seat Depth	430mm
Seat Height	550mm
Seat Width	500mm
Backrest Height	400mm
Max User Weight	200kg
Brand	K Care

Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind



Slide out/retractable non-weight bearing platform



Swing up arms

Available in 445mm, 500mm, 550mm, 600mm and 650mm seat widths

Attendant Sliding Foot Plate

Transit Mobile Shower Commode

Features

- Independently tested to conform to AS/NZS 3973:2009 for safety and peace of mind
- 3 strap back rest for comfort and easy access
- Integrated push-handle for manoeuvrability
- Swing-up arms for unrestricted side transfer
- Slide out/retractable non-weight bearing platform
- 4x 125mm total locking castors
- Designed to be used with removable seats
- Available in 445, 500, 550, 600 and 650mm seat widths
- Padded seat not included

KA112S

KA112S50

KA112S55





KA112S60

KA112S65

Specifications

KA112S

Overall Depth	840mm
Overall Height	1005mm
Overall Width	480mm
Seat Width	445mm
Seat Depth	430mm
Seat Height	550mm
Armrest Height	285 mm
Backrest Height	400mm
Max User Weight	200kg
Brand	K Care



For use on K Care REHAB
mobile shower commodes

Lateral and Hip/Thigh Support

Features

- For use on K Care REHAB mobile shower commodes
- Available with 35mm or 85mm adjustable linkage or fixed height
- Effortless push button release mechanism
- Depth adjustable pad
- Height adjustable pads

KA1100LAT-FA1-35L


KA1100LAT-FA1-85R

KA1100LAT-FA3-35L

KA1100LAT-FA4-R

For all model specifications please visit www.kcare.com.au

Strong with larger
footplate size



Compatible with K Care mobile shower commode
swing away and elevating legrest styles

Angle Adjustable Footplates

KA11AAFP

Features

- Suitable for swingaway and elevating legrest styles
- Simple angle adjustment
- Strong with larger footplate size
- Footplate can be positioned several different ways to accommodate for a wider range of users
- Height adjustable

Specifications

KA11AAFP

Overall Width	150mm
Overall Depth	210mm
Brand	K Care



Front



Side

Top

Innovative design. Superior comfort. Introducing the 'Dar-Avoid' pressure care commode Seat

'Dar-Avoid' Pressure Care Commode Seat

KAT1CPM-0

KAT1CPM-050

KAT1CPM-055

Features

- Highest level of pressure reducing shower commode seats on the market
- Multiple layers of high quality foam are contoured, maximising surface area contact for the occupant
- Designed to offload pressure from the ischial tuberosity (IT) and coccyx areas
- The seat is customisable and can be modified to suit individual requirements
- The combination of our Dar Avoid Commode Seat and a tilt/tilt recline shower commode allows for a seated posture that is both comfortable and functional, reducing pressure on skin and soft tissue

Specifications

KAT1CPM-0

Overall Length	460mm
Overall Width	445/500/550mm
Cushion Height	80mm to 100mm
Aperture Width	200mm
Aperture Depth	260mm
Front Opening	90mm
Brand	K Care

Open, closed, left and right side opening seats available



Waterproof cover for ease of cleaning



Single density vinyl covered commode seat

Handmade Open Seat Blue Vinyl

Features

- Waterproof cover for ease of cleaning
- Anti fungal and anti bacterial
- Flame retardant and abrasion tested
- Easily maintained
- Superior tear strength
- Open, closed, left and right side opening seats available
- Standard or pressure relieving foam
- Sizes from 320 to 650mm wide. For additional information please visit www.kcare.com.au

KA11HOP-B

KA11O65-B

Specifications

KA11HOP-B

Overall Width	445mm
Overall Depth	435mm
Overall Height	55mm
Aperture Width	200mm
Brand	K Care

KA11O65-B

Overall Width	650mm
Overall Depth	480mm
Overall Height	55mm
Aperture Width	200mm
Brand	K Care

Stable stainless steel frame



Stable even at full tilt

Ocean Ergo VIP

Features

- Stable stainless steel frame
- Seat tilt from -5 to 40°
- Tool-less height adjustable seat
- Machine washable backrest
- Ergonomically designed frame

PKG-1643577

Specifications

PKG-1643577

Overall Width	565mm
Overall Depth	1000mm
Seat Width	480mm
Seat Depth	480mm
Min Seat Height	475mm
Max Seat Height	600mm
Product Weight	23.5kg
Max User Weight	150kg
Brand	Aquatec

Seat tilt from 0-35 degrees



Ergonomically contoured seat

Ocean Dual VIP Ergo

PKG-1643717

Features

- Seat tilt from -5 to 40°
- Backrest recline of 0-35°
- Tool-less height adjustable seat
- Stable even at full tilt
- Machine washable backrest
- Ergonomically contoured seat

Specifications

PKG-1643717

Overall Width	565mm
Overall Depth	1000mm
Seat Width	480mm
Seat Depth	480mm
Min Seat Height	475mm
Max Seat Height	600mm
Product Weight	28.3mm
Max User Weight	150kg
Brand	Aquatec

Large user weight rating



Removable splash guard

Padded Seat

Over Toilet Frame

Features

- Sturdy powder coated steel frame
- Large user weight rating
- Removable padded seat for comfort and easy cleaning
- Removable splash guard
- Available in three different widths to suit bariatric users (550, 600 and 650mm wide)

KA410AP55

KA410AP60

KA410AP65

Specifications

KA410AP55

Overall Width	600mm
Seat Width	550mm
Seat Depth	460mm
Floor To Seat Height	520mm
Product Weight	7kg
Max User Weight	300kg
Brand	K Care

Fully padded, upholstered commercial grade vinyl



Fixed height



Powder coated finish



Available in three different widths to suit bariatric users

Padded Back Seat and Arms

KA220ZA55

Maxi Shower Chair

Features

- Fixed height for extra stability
- Rigid, anti-corrosive zinc treated, mild steel frame
- Fully padded, upholstered commercial grade vinyl
- Powder coated finish
- Available in three different widths to suit bariatric users (550, 600 and 650mm wide)

Specifications

KA220ZA55

Overall Depth	625mm
Overall Height	940mm
Overall Width	685mm
Seat Width	550mm, 600mm, and 650mm
Seat Depth	470mm
Backrest Height	420mm
Armrest Height	170mm
Backrest Height	420mm
Max User Weight	300kg
Product Weight	8kg
Brand	K Care



Available in 450mm, 520mm and 600mm seat widths

Premium Shower Stool with PU seat

12325A52

Features

- Designed to assist those who have difficulty showering
- Height adjustable legs
- PU padded seat for extra comfort
- Powder coated white 32mm aluminium tubing
- Available in 450mm, 520mm and 600mm seat widths

Specifications

12325A52

Max Overall Height	810mm
Max Seat Height	640mm
Min Seat Height	480mm
Seat Depth	300mm
Seat Width	450mm, 520mm and 600mm
Max User Weight	250kg
Brand	R & R Healthcare Equipment

Australian standards compliant

Plastic moulded seat

Height adjustable legs

Independently certified to meet Australian Standards

Premium Shower Chair

12028A48

Features

- Allows seating in the shower
- Armrests for extra support and help positioning
- Powder coated white 32mm aluminium tubing
- Solid plastic seat
- Height adjustable legs
- Certified to meet Australian Standards

Specifications

12028A48

Max Overall Height	810mm
Max Seat Height	640mm
Min Seat Height	400mm
Seat Depth	300mm
Seat Width	480mm
Max User Weight	130kg
Brand	R & R Healthcare Equipment



Larger 28mm tubing makes it easier to grip handles when transferring



Contoured seat with inbuilt handle and water drainage slots

Anti-slip rubber tips on each leg for stability

Anti-slip rubber tips on each leg for stability

Shower Chair Heavy Duty

KA220AAHD

Features

- Designed for users preferring additional seating support while showering
- Larger 28mm tubing makes it easier to grip handles when transferring
- Adjustable height and angle to suit individual users and uneven floors
- Contoured seat with inbuilt handle and water drainage slots
- Seat clips on for easy assembly, change and cleaning
- Anti-slip rubber tips on each leg for stability
- Durable and long lasting aluminium frame
- Independently tested to AS/NZS 3973:2009
- Suitable for institutional and domestic use


Specifications

KA220AAHD

Overall Width	530mm
Overall Depth	520mm
Seat Width	440mm
Seat Depth	465mm
Min Floor To Seat Height	430mm
Max Floor To Seat Height	550mm
Width Between Arms	490mm
Product Weight	3.8kg
Max User Weight	200kg
Brand	K Care



Robust and anti corrosive aluminium frame



Moulded plastic clip-on seat for easy cleaning

Anti-slip rubber tips on each leg for stability

Durable and long lasting aluminium frame

Shower Stool Heavy Duty

KA222AAHD

Features

- Designed for users preferring additional seating support while showering
- Larger 28mm tubing makes it easier to grip handles when transferring
- Adjustable height and angle to suit individual users and uneven floors
- Moulded plastic clip-on seat for easy cleaning
- Anti-slip rubber tips on each leg for stability
- Durable and long lasting aluminium frame
- Independently tested to AS/NZS 3973:2009
- Suitable for institutional and domestic use


Specifications

KA222AAHD

Overall Width	525mm
Min Overall Depth	380mm
Max Overall Depth	410mm
Seat Width	350mm
Seat Depth	300mm
Min Floor To Seat Height	510mm
Max Floor To Seat Height	640mm
Width Between Arms	480mm
Product Weight	2.8kg
Max User Weight	200kg
Brand	K Care



Larger 28mm tubing makes it easier to grip handles when transferring



Supplied with toilet flap with removable splash guard

Anti-slip rubber tips on each leg for stability

Designed for users preferring additional seating support while toileting

Over Toilet Frame Heavy Duty

KA410AAHD

Features

- Larger 28mm tubing makes it easier to grip handles when transferring
- Adjustable height and angle to suit individual users and uneven floors
- Seat clips on for easy assembly, change and cleaning
- Supplied with toilet flap and removable splash guard
- Anti-slip rubber tips on each leg for stability
- Durable and long lasting aluminium frame
- Independently tested to AS/NZS 3973:2009
- Suitable for institutional and domestic use

Specifications

KA410AAHD

Overall Width	530mm
Min Overall Depth	570mm
Max Overall Depth	605mm
Seat Width	380mm
Seat Depth	440mm
Min Floor To Seat Height	470mm
Max Floor To Seat Height	595mm
Product Weight	2.8kg
Max User Weight	200kg
Brand	K Care

Manufactured from zinc treated mild steel the anti-corrosive frame has a powder coat finish



Frame folds up quickly and easily



This model is height adjustable



Designed for convenient storage and for travellers

Over Toilet Frame Folding

KA430Z

Features

- Frame folds up quickly and easily
- Manufactured from zinc treated mild steel the anti-corrosive frame has a powder coat finish
- Adjustable seat height and angle

Specifications

KA430Z

Overall Width	510mm
Overall Depth	600mm
Overall Height	680mm to 805mm
Seat Width	450mm
Seat Height	480mm to 605mm
Max User Weight	110kg
Brand	K Care

Ultra light and extra strong with 250kg weight capacity



Designed to assist those having difficulty using a standard toilet

Over Toilet Frame

Features

- Ultra light and extra strong with 250kg weight capacity
- Australian standards approved
- Anti-corrosive aluminium frame
- Comfortable seat with a lid
- Adjustable seat height and angle
- Moulded splash guard and bucket included
- Available in 450, 520, and 600mm widths (seat width does not change)
- Padded seat options available to provide larger seating surface

12247A45

Specifications

12247A45

Seat Width	450mm
Seat Depth	410mm
Min Seat Height	480mm
Max Seat Height	640mm
Product Weight	3kg
Max User Weight	250kg
Brand	R & R Healthcare Equipment



A range of grab rails to suit you

Grab Rails

Features

- Straight or configured models available
- Durable Stainless steel tube
- 32mm diameter tube for easy grip

30032S

45032S

60032S

75032S

90032S

120032S

013LUDL

013LUDR

Specifications

30032S to 120032S

Standard lengths available 300, 450, 600, 750, 900, 1200mm
Rails may be mounted vertically or horizontally etc. Over 1200mm length requires centre support.

013LUDL

The 013LUD is manufactured from durable stainless steel. The cane is 32mm diameter tube and locks into both the raised and lowered position. The locking mechanism is an easy to use pull pin assembly and comes in either right hand (013LUDR) or left hand configuration (013LUDL). Please note the pull button must always be fitted facing away from the toilet pan.

Brand

K Care

Seating and Lift Chairs



Lift / Tilt Function



Shiatsu Massage



Heated Seat

Zero gravity, Shiatsu massage and heated seat

NEW
PRODUCT

Donatello Plus Lift Recliner

KA554-PL

Features

- Zero gravity tilt-in-space function provides a feeling of weightlessness
- Shiatsu massage – choose from 3 pressure levels
- Heated seat – known to reduce stress, relieve pain, improve circulation, and increase comfort
- Lift and tilt function
- Individual posture supports with 5 powerful individual motors
- Wireless phone charger
- Drinks cup holder
- Airline head support


Specifications

KA554-PL

Overall Width	1010mm
Overall Depth	1060mm
Overall Height	1090mm
Seat Width	550mm
Seat Depth	530mm
Max Seat Height	495mm
Product Weight	79kg
Max User Weight	158kg
Brand	Alivio



The new lift and tilt function ensures safe entry and exit from the chair



Air Massage

Wireless Charger
& Cup Holder

Full featured lift & tilt recliner with a great range of features for comfort, support and safety



Carrington Lift Recliner

KA558

Features

- Air massage – choose from three pressure levels
- Lift and tilt function – ensure that entering and exiting the chair can be completed safely
- Individual posture supports – easy individual support adjustments with 4 independent motors
- Premium fabric – ultra-durable, comfortable, comfort, long lasting, and easy-to-clean
- Soft-touch controller – air massage, lift & tilt, and postural support functions, with USB port for extra charging capability
- Drinks cup holder
- Wireless phone charger

Specifications

KA558

Overall Depth	1060mm
Overall Width	1010mm
Overall Height	1090mm
Seat Depth	530mm
Seat Width	550mm
Seat Height	508mm
Product Weight	85kg
Max User Weight	158kg
Brand	Alivio



4 independent motors for wide range of position combinations and support

Michelangelo Lift Recliner

KA553

Features


- Individually adjustable head support
- Headrest and lumbar system
- All close and lift button
- Footrest extension
- Soft-touch controller



Specifications

KA553

Overall Depth	1080mm
Overall Width	890mm
Overall Height	1067mm
Seat Depth	546mm
Seat Width	559mm
Seat Height	470mm
Product Weight	70kg
Max User Weight	158kg
Brand	Alivio



Michelangelo petite
four motor lay flat recliner

Michelangelo Petite Lift Recliner

Features

- Individually adjustable head support
- Headrest and Lumbar System
- All Close and Lift Button
- Footrest Extension
- Soft-touch controller

KA556

Specifications

KA556

Seat to Floor	457mm
Seat Depth	457mm
Seat Width	457mm
Chair Depth	965mm
Chair Width	812mm
Chair Height	1016mm
Max User Weight	136kg
Brand	Alivio


Lateral Support Backrest



Standard Support Backrest



High Support Backrest



Donatello independent
four motor lay flat recliner

Donatello Lift Recliner

Features

- Dual motor lift lay flat recliner
- Individually adjustable head support
- Individually adjustable lumbar support
- Power backrest & power footrest
- Headrest and lumbar system
- All close and lift button
- Footrest extension
- Soft-touch controller

KA554-TB

KA554-SB

KA554-LS

Specifications

KA554-TB

KA554-SB

KA554-LS

Overall Width	890mm
Overall Depth	1016mm
Overall Height	1092mm
Seat Width	520mm
Seat Depth	520mm
Max Seat Height	508mm
Product Weight	74kg
Max User Weight	158kg
Brand	Alivio

Lateral support cushion height: 690mm

Standard support cushion height: 690mm

High support cushion height: 770mm



The new lift mechanism has a more vertical lift with a slight tilt compared to the forward tilt you find on generic chairs



Lay-flat recliner



Footrest extension

Combines modern design aesthetics with an ultra-durable fabric

Donatello Canyon Lift Recliner

KA554LSC

Features

- Dual lift function
- Lay-flat recliner
- Footrest extension
- All close and all lift buttons
- Hand controller
- Ultra durable fabrics
- Lumbar support to accommodate sitting posture
- Soft-touch controller

Specifications

KA554LSC

Overall Depth	1016mm
Overall Width	890mm
Overall Height	1092mm
Seat Depth	520mm
Seat Width	520mm
Seat Height	508mm
Reclined Length	1830mm
Product Weight	74kg
Max User Weight	158kg
Brand	Alivio


Individually adjustable lumbar support



Dual motor lift lay flat recliner



Improved lift mechanism



Leonardo independent four motor lay flat recliner

Leonardo Lift Recliner

Features

- Dual motor lift lay flat recliner
- Individually adjustable head support
- Individually adjustable lumbar support
- Power backrest & power footrest
- Headrest and lumbar system
- All close and lift button
- Footrest extension
- Soft-touch controller

KA552

Specifications

KA552

Overall Depth	1067mm
Overall Width	902mm
Overall Height	1092mm
Seat Depth	521mm
Seat Width	571mm
Seat Height	483mm
Max User Weight	158kg
Brand	Alivio



Full featured lift recliner with modern design aesthetics

Kahlo Lift Recliner

KA557

Features

- Dual lift function
- Lay-flat recliner
- Footrest extension
- All close and all lift buttons
- Hand controller
- Ultra durable fabrics
- Support to accommodate sitting posture
- Soft-touch controller

Specifications

KA557


Overall Depth	1016mm
Overall Width	890mm
Overall Height	1092mm
Seat Depth	520mm
Seat Width	520mm
Seat Height	508mm
Reclined Length	1830mm
Product Weight	74kg
Max User Weight	158kg
Brand	Alivio



Dual motor lift lay flat recliner



Footrest extension



Easy to operate handset

Picasso independent four motor lay flat chair

Picasso Lift Recliner

Features


- Dual motor lift lay flat recliner
- Individually adjustable head support
- Individually adjustable lumbar support
- Power backrest & power footrest
- Headrest and lumbar system
- All close and lift button
- Footrest extension
- Soft-touch controller

KA551

Specifications

KA551

Overall Depth	1003mm
Overall Width	813mm
Overall Height	1080mm
Seat Depth	508mm
Seat Width	521mm
Seat Height	495mm
Max User Weight	158kg
Brand	Alivio



Total surface air pressure system
for long term pressure relief

Compact Lift Chairs

Features

- Removable seat and back cushions
- Electric lift and recline mechanism
- Front and rear castors
- Wall Saver mechanism
- Back-up battery capability
- Easy to clean, repair and maintain
- VersoTela Fusion and OsGuard PU fabric
- Available in single or twin motor

AC59098

AC59087

AC59099

AC59100

AC59090

AC59101

Specifications

AC59090

Overall Length	820 mm
Overall Width	750 mm
Overall Height	1060 mm
Seat Width	500 mm
Seat Depth	510 mm
Seat Height	520 mm
Backrest Height	740 mm
Max User Weight	140 kg
Brand	Air Comfort


For all model specifications please visit www.kcare.com.au



Alternating pressure care system



Deep air cells



Total surface air pressure system
for long term pressure relief

NEW
PRODUCT

Compact Mobicare Recliner

Features

- Carilex IPS technology automatically adjusts the internal cushion pressure
- Unique combination of the Air Comfort Compact Lift Chairs and the Carilex MobiCare Plus alternating air cushion
- Electric lift and recline mechanism
- Anti-decubitus air alternating cushion to reduce sores
- This device facilitates circulation and decreases contact pressure
- Available in single or twin motor
- Mobicare cushion: Specially designed deep air cells for excellent immersion and pressure relief

AC59130

AC59131

AC59132

AC59133

AC59134

AC59135

AC59136

Specifications

AC59132

Overall Length	820 mm
Overall Width	750 mm
Overall Height	1060 mm
Seat Width	500 mm
Seat Depth	510 mm
Seat Height	520 mm
Pressure Sore Grading	3
Waterlow Scale	15+
Product Weight	65 kg
Max User Weight	140 kg
Brand	Air Comfort

For all model specifications please visit www.kcare.com.au



Total surface air pressure system
for long term pressure relief

Compact Lift Chairs Titanium

AC59125

AC59126

AC59127

Features

- Removable seat and back cushions
- Electric lift and recline mechanism
- Front and rear castors
- Wall Saver mechanism
- Back-up battery capability
- Easy to clean, repair and maintain
- Available in VersoTela Fabric
 - Theme Range - Titanium colour
- Available in Small, Medium and Large sizes
- Available in Twin Motor

Specifications

AC59126


Overall Length	820 mm
Overall Width	750 mm
Overall Height	1060 mm
Seat Width	500 mm
Seat Depth	510 mm
Seat Height	520 mm
Backrest Height	740 mm
Max User Weight	140 kg
Brand	Air Comfort

For all model specifications please visit www.kcare.com.au


Adjustable back, recline, seat, tilt and raising leg rest



Adjustable arms rest



Adjustable back wings



Market leader in pressure care seating

Standard Deluxe V2

Features

- Ease of handling for carers and maximum comfort for patients
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the user
- Adjustable back, recline, seat, tilt and raising leg rest
- Drop away/adjustable arms rest and back wings
- Carrflex extreme multi-stretch materials with sanitised treatment and enhanced durability
- Large twin castors make the chair easy to manoeuvre
- One directional lock castors to assist with straight line tracking

AC59111

Specifications

AC59111

Overall Length	1130mm
Overall Height	1260mm
Overall Width	680mm
Seat Width	470mm
Floor to Seat Height	570mm
Product Weight	45kg
Max User Weight	180kg
Brand	Air Comfort

Colours Available:



Electric adjustable back recline,
seat tilt and raising leg rest



Patented Airgonomic Seating System uses
four air cells to immerse and cushion the user

Deluxe Bed V2

AC59112

Features

- Recommended user height:
Up to 1620mm to 1880mm
- Electric adjustable back recline, seat tilt
and raising leg rest
- Drop away/adjustable arm rests and back wings
- Carrflex extreme multi-stretch covering with
sanitised treatment and enhanced durability
- Large castors make the chair easy to manoeuvre
- One directional lock castor to assist with
straight line tracking
- Optional: Power drive version now available

Specifications

AC59112

Overall Length	1130mm
Overall Width	680mm
Overall Height	1260mm
Seat Width	470mm
Seat Depth	550mm
Floor To Seat Height	570mm
Brand	Air Comfort
Product Weight	45kg
Max User Weight	180kg
Brand	Air Comfort

Colours Available:



Patented Airgonomic Seating System uses four air cells to immerse and cushion the users



Unparalleled pressure redistribution and adjustable positioning

Maxi Deluxe V2

AC59115

Features

- Recommended user height 1620 – 1880mm
- Patented Airgonomic Seating System uses four air cells to immerse and cushion the users
- Electric adjustable back recline, seat tilt and raising leg rest
- Drop away / adjustable arm rests and back wings
- Carrflex extreme multi-stretch covering with sanitised treatment and enhanced durability
- Large twin castors make the chair easy to manoeuvre
- One directional local castor to assist with straight line tracking
- Ergonomically designed handlebars


Specifications

AC59115

Overall Length	1360mm
Overall Height	1385mm
Overall Width	950mm
Seat Width	720mm
Floor to Seat Height	640mm
Product Weight	100kg
Max User Weight	300kg
Brand	Air Comfort

Colours Available:





Contoured seat to fit natural body shape and relieve pressure points

Fusion High Back Pressure Chair

Features

- Backrest is shaped with lumbar support for higher level of comfort
- Standard with depth adjustable backrest and height adjustable leg rests
- Standard with tilt handle and non-marking transit wheels
- Standard with full padded arm rests (not plastic)
- Removable head rest

18030

Specifications

18030


Seat Width	520mm
Min Seat Depth	500mm
Max Seat Depth	650mm
Backrest Height	1130mm
Min Seat Height	450mm
Max Seat Height	550mm
Max User Weight	250kg
Brand	R & R Healthcare Equipment



Contoured backrest specifically designed contoured back for lumbar support



Knocks down into a compact



Chair width is adjustable

Lightweight and sturdy, 200kg maximum user weight

HiLite Chair

KA582KDG05

KA582KDV03

Features

- Knock down design
- Put together with no tools required
- 60% smaller than shipping and storing a fully assembled chair
- Lightweight and sturdy, this chair is designed for comfort and function
- Contoured backrest to provide lumbar and back support
- Seat height and angle are adjustable
- Wide variety of models and accessories available

Specifications

Greystone KA582KDG05

Seat Width	490mm to 550mm
Seat Depth	480mm to 500mm
Seat Height	430mm to 580mm
Backrest Height	700mm
Overall Width	640mm to 700mm
Max User Weight	200kg
Brand	K Care

Champagne KA582KDV03

Seat Width	490mm to 550mm
Seat Depth	480mm to 500mm
Seat Height	430mm to 580mm
Backrest Height	700mm
Overall Width	640mm to 700mm
Max User Weight	200kg
Brand	K Care



Anti fungal



Fire resistant



Wear resistant



Anti bacterial

Colours Available:





Suitable for home and hospital environments

Kingston Chair

Features

- Multi purpose height adjustable chair
- Large padded backrest, seat and armrests
- Suitable for home and hospital environments
- Available in blue or fawn vinyl
- Available in steel or aluminium frames
- Available in mobile model with dropside arms

- 18005AF
- 18005AF52
- 18005AF60
- 18005SF
- 18005AB
- 18005SB
- 18005SMF

Specifications

18005AF

Seat Width	450 – 600mm
Seat Depth	440mm
Min Seat Height	440mm
Max Seat Height	560mm
Backrest Height	460mm
Max User Weight	200kg
Brand	R & R Healthcare Equipment

Colours Available:





Suitable for home and hospital environments

**NEW
PRODUCT**

Kingston Mobile Chair

18005SMG

Features

- All New Frame designed in Australia
- Adjustable Backrest Depth
- Adjustable Height
- Light Weight
- Maximum user weight to 160 kg
- Independently tested to Australian Standards
- Dual Locking Castors
- Armrests lower level with seat for easy side transfers
- Integrated Push Handle
- Perfect for dining table
- Not to be used as a wheelchair

Specifications

18005SMG

Weight	13kg
Overall Height	915 - 1015mm
Overall Width	600mm
Overall Length	630mm
Seat Width	450mm
Seat Height	460 - 560mm
Seat Depth	440 - 490mm
Max User Weight	160kg
Castors	4 dual castors, swivel locking 75mm / 3"
Materials	Powder coat steel frame, stainless steel arm rests, PU foam seat/backrest textiles, plastic accessories



KA586KDDV03



KA586KDDG05

Suitable for home and hospital environments

Desk Utility Chair

Features

- Width adjustable through Knock Down design
- Unique desk/table friendly frame design allows user to get closer to desks/tables
- Shorter padded armrests
- Contoured backrest for increased lumbar support
- Anti fungal and anti bacterial fabric
- Fire resistant and wear resistance tested
- Seat height adjustability

KA586KDDV03

KA586KDDG05

Specifications

KA586KDDV03

Min Width	600 mm
Max Width	660 mm
Overall Depth	510 mm
Min Seat Width	460 mm
Max Seat Width	510 mm
Seat Depth	460 mm
Min Seat Height	430 mm
Max Seat Height	580 mm
Backrest Height	440 mm
Max User Weight	200 kg
Brand	K Care

Colours Available:





Contoured backrest specifically designed for lumbar support and comfort



Knocks down into a compact size

Chair width is adjustable

Seat width and depth adjustability as standard

Utility Chair Standard

KA586KDG05

KA586KDV03

Features

- Adjustable width and height
- Padded arms and back
- Low back version has curved backrest providing extra postural support
- Rubber stoppers for stability
- Suitable for use in any area
- Fire resistant and wear resistance tested
- Anti fungal and anti bacterial fabric

Specifications

Greystone KA586KDG05

Seat Width	460mm to 510mm
Seat Depth	460mm to 520mm
Overall Width	600mm to 660mm
Overall Depth	510mm
Backrest Height	440mm
Seat Height	430mm to 580mm
Max User Weight	200kg
Brand	K Care

Champagne KA586KDV03

Seat Width	460mm to 510mm
Seat Depth	460mm to 520mm
Overall Width	600mm to 660mm
Overall Depth	510mm
Backrest Height	440mm
Seat Height	430mm to 580mm
Max User Weight	200kg
Brand	K Care



Anti fungal



Fire resistant



Wear resistant



Anti bacterial

Colours Available:



New knock down design



Knocks down into a compact size



Height adjustability with total range of 150mm



Angled, height adjustable

Padded Leg Rest

Features

- Height adjustable leg rest
- Angled for added comfort
- Designed to support the legs while sitting
- New knock down design
- Put together with no tools required

KA58LV03

KA58LG05

Specifications

Champagne KA58LV03

Min Overall Height	420mm
Max Overall Height	570mm
Overall Width	500mm
Overall Length	400mm
Padded Thickness	40mm
Product Weight	4.2kg
Max User Weight	200kg
Brand	K Care

Greystone KA58LG05

Min Overall Height	420mm
Max Overall Height	570mm
Overall Width	500mm
Overall Length	400mm
Padded Thickness	40mm
Product Weight	4.2kg
Max User Weight	200kg
Brand	K Care



Anti fungal



Fire resistant



Wear resistant



Anti bacterial

Colours Available:



Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Lifters and Slings



With the touch of a button, the Multi Lift is capable of lifting patients in either a standing or sitting position



Powered lifting arm and leg gating control



The Multi Lift has many features and is compatible with a range of slings

Multi Lift

KH401

Features

- With the touch of a button, the Multi Lift is capable of lifting patients in either a standing or sitting position
- Powered lifting arm and leg gating control
- Multi-function handset makes for easier, faster and safer patient lifting and transport
- Intelligent electronics within the lifter record usage data, battery charge and the condition of the battery
- Electronic overload feature prevents the lifter operating if it detects a weight above its capacity
- A maintenance capacity required indicator prompts users to service the Multi Lift to ensure ongoing safety and extend the operating life

Specifications

KH401

Overall Width	700mm
Overall Length	1310mm
Overall Height	1220mm
Min Lift Height	430mm
Max Lift Height	1700mm
Turning Radius	650mm
Product Weight	58kg
Max User Weight	180kg
Brand	Kerry

Electric lifting arm and leg adjustment handset



Ideal for use in aged, community and hospital environments

Patient Lifter **Easy Lift**

KH405F

Features

- Independently tested to conform to ISO 10535 for safety and peace of mind
- Linak Actuation
- Electric lifting arm and leg adjustment handset
- Plug in rechargeable battery system
- Emergency stop button

Specifications

KH405F

Overall Width	700mm
Overall Length	1300mm
Overall Height	1315mm
Product Weight	53kg
Max User Weight	225kg
Brand	Kerry

Wide lateral push handles
for manoeuvrability



Independently tested to conform to
ISO 10535 for safety and peace of mind

Patient Lifter

Quick Stand

KH400G


Features

- Linak Actuation
- Electric lifting arm and leg adjustment using handset
- Plug in rechargeable battery system
- Hand operated emergency lowering
- Emergency stop button
- Wide lateral push handles for manoeuvrability
- Adjustable height knee pad
- Fold up foot plate

Specifications

KH400G

Overall Length	1050mm
Overall Height	910mm
Overall Width	670mm
Min Lift Height	820mm
Max Lift Height	1620mm
Product Weight	45kg
Max User Weight	200kg
Brand	Kerry



Easily installed between boom and yoke

Digital Lift Scale for Patient Lifters

KH40EWS

Features

- Digital lift scale for patient lifters
- Weigh scale attachment for K Care patient lifters
- Easy to read LCD display
- Easily fitted to the Kerry lifter range
- Maximum weight measuring capacity of 400kg (subject to weight capacity of patient lifter)


Specifications

KH40EWS

Maximum weight measuring capacity of 400kg

Brand

Kerry



Padding offers excellent patient comfort and support

Sling Options

Features

- Ideal for transferring and bathing patients
- Excellent back and leg support
- Quick and easy to fit
- Machine washable
- Max user weight 250kg
- Various models and sizes available

KH406NE2

KH400SSL

KH408NE3

KH407N2

For all model specifications please visit www.kcare.com.au

Types of slings


- General purpose; mesh or poly
- Standing
- Sitting
- Full Body

Clinical Equipment





Self-locking stainless steel drop side rails & Stainless steel push handles at either end



3 x 175mm diameter total lock castors and 1 x 175mm diameter directional lock castor



Pneumatic assisted 10° tilt top for drainage

Practical and supportive

Electric Shower Trolley

KH600E

Features

- Heavy duty padded vinyl liner and pillow (included)
- Self-locking stainless steel drop side rails & Stainless steel push handles at either end
- Pneumatic assisted 10° tilt top for drainage
- 3 x 175mm diameter total lock castors and 1 x 175mm diameter directional lock castor
- Electric height adjustment with plug in rechargeable battery system (Waterproof)
- Also available with hydraulic height adjustment and in short, long and wide liner platforms

Specifications

KH600E

Overall Length	1850mm
Overall Height (Floor to Platform)	635mm – 935mm
Overall Height (Lowest Setting)	835mm
Liner Height	180mm
Overall Width	650mm
Liner Length	1650mm
Product Weight	80kg
Liner Width	550mm
Max. User Weight	150kg
Brand	Kerry


Precision-welded steel frame
with high-impact panels



Defibrillator Shelf



Cardiac Board



Easily manoeuvrable
with excellent stability

Bravo Series

Emergency Cart

KH810S

Specifications

KH810S

Overall Length	790mm
Overall Width	610mm
Overall Height	1000mm
Number of Drawers	5
Drawers 75mm	2
Drawers 160mm	2
Drawers 250mm	1
Brand	Bravo

Features

- Precision-welded steel frame with high-impact panels
- Seal system with central locking drawers
- Expansive work area with slide-out utility surface
- Flexible divider system

Included

- Cardiac Board
- Extended Utility Surface
- Defibrillator Shelf
- IV Pole
- Oxygen Bottle Holder
- Self Closing, Soft Close Drawers
- Waste Bin with Lid
- 4x 125mm Castors, 2 with Brakes


Easily manoeuvrable
with excellent stability



Digital Lock Systems



Extended Utility Surface



Expansive work area with
slide-out utility surface

Bravo Series

Anaesthesia Cart

Features

- Precision-welded steel frame with high-impact panels
- Easily manoeuvrable with excellent stability
- Digital lock system
- Flexible divider system

Included

- Tilt Bins
- Extended Utility Surface
- Self Closing, Soft Close Drawers

KH820D

Specifications

KH820D

Overall Length	790mm
Overall Width	610mm
Overall Height	1000mm
Number of Drawers	5
Drawers 75mm	2
Drawers 160mm	2
Drawers 250mm	1
Brand	Bravo


- Waste Bin with Lid
- 4 x 125mm Castors, 2 with Brakes



Precision-welded steel frame with high-impact panels



Digital Lock Systems



Extended Utility Surface

Improved workflow and performance

Bravo Series

Medication Carts

KH830D

KH840D

Features

- Precision-welded steel frame with high-impact panels
- Easily manoeuvrable with excellent stability
- Digital lock system
- Expansive work area with slide-out utility surface
- Flexible divider system

Included

- Extended Utility Surface
- Self Closing, Soft Close Drawers
- Waste Bin with Lid
- 4 x 125mm Castors, 2 with Brakes

Specifications

KH830D


Overall Length	790mm
Overall Width	610mm
Overall Height	900mm
Number of Drawers	5
Drawers 75mm	3
Drawers 160mm	1
Drawers 250mm	1
Brand	Bravo

KH840D

Overall Length	790mm
Overall Width	610mm
Overall Height	1000mm
Number of Drawers	5
Drawers 75mm	2
Drawers 160mm	2
Drawers 250mm	1
Brand	Bravo



5 Drawers for medication,
instrument and equipment storage



Ergonomically
positioned push handle

Large 125mm diameter
castors for easy
manoeuvrability

Customised options available

Medication Storage Cart

KH220D

Features

- A.B.S Plastic removable recessed plain top for ease of cleaning & replacement
- 4 Small drawers 125mm high, 1 large draw 225mm high
- Universal locking drawers
- Corner buffers to prevent damage to the cart or surrounding equipment & structures
- 2 Locking castors at the front, 2 swivel at the rear
- Compact and stable design
- Colour options to match facility decor, location positioning or identification


Specifications

KH220D

Overall Width	700mm
Overall Depth	500mm
Overall Height	1000mm
Internal Drawer Width	550mm
Internal Drawer Depth	455mm
Castor Diameter	125mm
Product Weight	48kg

Colours Available





Strong and dependable

Stainless Steel

Instrument & Dressing Trolley

Features

- Stainless steel side frames and shelves
- Fitted with heavy duty, grey rubber, easy to push and durable 100mm diameter castors
- Integrated push-handle for manoeuvrability
- Guard rails on 3 sides of top and shelf
- Rotating buffers to protect surrounding items
- Available with sliding drawer to hold dressings and instruments (KH141SD)

Most of the products in our Kerry range can be customised to suit the specific needs of our customers.

KH141S

KH141SD

Specifications

KH141S

Overall Width	500mm
Overall Length	500mm
Overall Height	1050mm
Product Weight	10kg
Max Load	25kg Per Shelf
Brand	Kerry

KH141SD

Overall Width	500mm
Overall Length	500mm
Overall Height	1050mm
Product Weight	10kg
Max Load	25kg Per Shelf
Brand	Kerry

Robust construction

75mm Twin Wheel Swivel Castors

Fits most linen bags

Strong and hygienic

Stainless Steel

Double Soiled Linen Trolley

KH541SL

Features

- Stainless steel only
- Robust construction
- Foot operated lid - Optional lid colours available
- Robust stainless steel sheet metal base
- Angled vertical frame for improved stability
- Fits most linen bags (optional extra)
- 75mm Twin Wheel Swivel Castors - Optional 100mm diameter swivel or locking castors
- Single and triple versions available with and without lids



Specifications

KH541SL

Overall Width	900mm
Overall Depth	460mm
Overall Height	885mm
Max Load	60kg
Product Weight	10kg
Brand	Kerry



Recessed 610mm recessed stainless steel 6 star base gives low centre of gravity and excellent stability



Stainless steel 4 prong IV hook to securely hold bags

Integrated power board and cable

Free flowing stability

Stainless Steel
Mobile Intravenous Stand

Features

- Stainless steel upright
- Height adjustable from 1525mm to 1985mm
- Damper to prevent upright dropping rapidly
- Offset mount pole to ensure syringe pump mounts over the centre of gravity
- Accepts most volumetric infusion, syringe and PCA pumps

Most of the products in our Kerry range can be customised to suit the specific needs of our customers.

KH313

Specifications

KH313

Overall Width	610mm
Overall Depth	610mm
Min Overall Height	525mm
Max Overall Height	1985mm
Product Weight	12kg
Brand	Kerry

Materials

Upright	Stainless Steel
Base	Stainless Steel
Castor housing	Stainless Steel
Castors	Rubber tyre castors for smooth movement
Product Weight	12kg

Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Daily Living Aids



Its non-weight bearing 'tilt and pour' action directs hot water to the fixed position of the cup, without removing the kettle from its cradle



Provides effortless pouring action with its ergonomic handle

Uccello Kettle

Features

- The Uccello Kettle is specifically designed using unique powerpour technology to provide effortless pouring action with its ergonomic handle
- Its non-weight bearing 'tilt and pour' action directs hot water to the fixed position of the cup, without removing the kettle from its cradle

PTKD01BW



PTKD01RW



PTKD01WW



PTKD01BB





The one-touch automatic jar opener will open stubborn lids at the touch of a button

One Touch **Jar Opener**

KC12

Features

- Ideal for people with weak or infirm hands
- Simply fit One Touch over the lid and let it do the hard work
- Suitable even for new, factory sealed jars
- Requires 2 x AA batteries



Automatic can opener with one-touch operation

One Touch

KC26

Can Opener

Features

- The innovative hands-free One Touch Automatic Can Opener opens most cans between 53mm and 153mm diameter
- Perfect for seniors, or anyone with weak hands
- The opener 'walks' around the can and switches itself off automatically, leaving an easy-to-remove lid with no sharp edges
- Requires 2 x AA batteries



Easy to use, store and clean, the one-touch bottle opener is perfect for people with weak or infirm hands

One Touch

KC85

Bottle Opener

Features

- The One Touch Twizoff Automatic Bottle Opener is suitable for either metal or plastic screw-top bottle caps up to 55mm in diameter
- Its rubber jaws adjust automatically to the lid size and will break seals and open caps with a smooth, effortless motion
- Requires 2 x AA batteries



Provides great stability, control and balance

Webequ Knives

120001W

11001

10001

Features

- Developed to secure the best, most effective and least demanding cutting and slicing movement
- Features a soft grip
- Ergonomic design minimizes strain on the user's wrist
- The shape and non-stick feature add to the elegant design of these knives



Minimising the risk of entrapment
when a bed stick is required

Patient Assist Rail

9077

Features

- Designed to assist patients when getting in and out of bed
- Slides under the mattress for easy application
- Suitable for use in any size beds
- Powder coated tubular steel frame
- Assists user to change positions in bed
- Multiple grip points
- Enclosed return to minimise injury
- 45mm Gap between shield and bar
- Aust Patent 2018100644

Simple handle design
for extra grip



Provides great stability, control and balance

Comfi Grip Sponge

HA075010CCG

Features


- Long handle allows users to reach hard to get places
- Simple handle design for extra grip
- Sponges do not absorb too much water
- Light weight

Specifications

HA075010CCG

Overall Length	650mm
Product Weight in Grams	150g

Magnetic tip and
dressing post



Designed to avoid
finger and wrist strain



Non-slip grip for firmer hold.

Classic Pro

Reacher Folding

Features

- Non-slip grip for firmer hold
- Ergonomic comfy grip handle hand grip
- Designed to avoid finger and wrist strain
- Magnetic tip and dressing post
- Available in short, standard and long

HA6421CCG

HA6232CCG

Specifications

HA6421CCG

Overall Length 530mm

HA6232CCG

Overall Length 850mm

Large grip handle



Stainless steel



1x each spoon, knife and fork



Perfect for those with
arthritis or limited hand grip

Knife, Fork and Spoon Package

KADL160ABC

Features

- Stainless steel utensil with large grip handle
- Perfect for those with arthritis or limited hand grip
- Dishwasher safe
- Hygienically sealed and easy to clean
- 1x each spoon, knife and fork

Specifications

KADL160ABC

Handle Length	110mm
Overall Length	210mm
Overall Width	40mm
Product Weight	105g

Long Handle Comb

KADL110

Features

- Made from durable plastic
- Rounded tooth design helps untangle knots and gentle on sensitive scalps
- Ergonomic rubberised handle for extra grip
- Light weight for extended use

Specifications

KADL110

Overall Length	380mm
Overall Height	100mm
Product Weight	200g



Long Handle Brush

KADL111

Features

- Made from durable plastic
- Padded bristles accommodate for sensitive scalps
- Ergonomic rubberised handle for extra grip
- Light weight for extended use

Specifications

KADL111

Overall Length	360mm
Overall Height	60mm
Product Weight	200g



Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Mobility Range





The Deluxe Knee Walker is intended to maximise and provide additional mobility to individuals during recovery from lower leg injury or surgery

Deluxe Knee Walker

Features

- Independently tested to conform to ISO 11199-2:2005 R2016 for safety and peace of mind
- Adjustable handle height
- Adjustable knee pad height
- Folds down to compact size to optimise
- Transport in vehicles
- Removable carry bag
- Easy to operate brake handle

KA227


Specifications

KA227


Overall Length	820mm
Overall Height	740-990mm
Overall Width	480mm
Knee Pad Height	510-610mm
Max User Weight	160kg



Ultimate safety and strength to support the user



Extra wide polyurethane seat



Soft Touch locking brakes

Ergonomic design to facilitate spacious support for the user

Premium Walker Medium Red

KA367MR

Features

- Adjustable to suit all users of various size
- Bag and basket are standard
- Lightweight aluminium frame
- Easily foldable frame
- Available in various sizes and colors

Specifications

KA367MR

Min Handle Height	830mm
Max Handle Height	95mm
Overall Length	690mm
Overall Width	620mm
Seat Width	410mm
Seat Depth	350mm
Floor To Seat Height	550mm
Castor Diameter	200mm
Product Weight	8kg
Max User Weight	200kg
Brand	K Care

Suitable for
indoor use only



Extra wide
polyurethane seat



Soft Touch
locking brakes



Ultimate safety and strength
to support the user

Petite Walker

KA367P

Features

- Adjustable to suit all users of various size
- Bag and basket are standard
- Ergonomic design to facilitate spacious support for the user
- Lightweight aluminium frame
- Easily foldable frame

Specifications

KA367P


Min Handle Height	755mm
Max Handle Height	895mm
Overall Width	600mm
Overall Height	350mm
Seat Depth	350mm
Floor To Seat Height	450mm
Diameter	180mm
Castor Diameter	150mm
Product Weight	6.3kg
Max User Weight	136kg
Brand	K Care



Lightweight
aluminium frame



Bag and basket
are standard



Soft Touch
locking brakes

Approved walking frame
by Bowls Australia

BOWLS Walker

KA367BW

Features

- Extra wide polyurethane seat
- Ultimate safety and strength to support the user
- Ergonomic design to facilitate spacious support for the user
- Easily foldable frame


Specifications

KA367BW

Min Handle Height	780mm
Max Handle Height	920mm
Overall Width	620mm
Seat Width	410mm
Seat Depth	350mm
Floor To Seat Height	520mm
Diameter	200mm
Product Weight	6kg
Max User Weight	200kg
Brand	K Care



T-Handle



T Handle



Swan Neck



Regular Foot



Quad Foot*

Stability when you need it

Cane Range

Features

- The ergonomically-designed Contour Comfort Grip is shaped to comfortably fit your hand
- The durable soft grip material, along with finger grooves on the underside, ensures maximum grip and traction
- The Contour Comfort Grip has a sturdy wrist-strap to instill added confidence and security
- The canes are intended for indoor and outdoor use

KA305BLK	T Handle Folding Cane Anodised Black
KA305BLU	T Handle Folding Cane Anodised Blue
KA305BRN	T Handle Folding Cane Anodised Brown
KA305BRZ	T Handle Folding Cane Anodised Bronze
KA322BLK	T Handle Cane Anodised Black
KA322BLU	T Handle Cane Anodised Blue
KA322BRN	T Handle Cane Anodised Brown
KA322BRZ	T Handle Cane Anodised Bronze
KA322FI4	T Handle Cane Pattern Light Pink Flower
KA322GRN	T Handle Cane Anodised Green
KA326BLK	Swan Neck Cane Anodised Black
KA326BLU	Swan Neck Cane Anodised Blue
KA326BRN	Swan Neck Cane Anodised Brown
KA326BRZ	Swan Neck Cane Anodised Bronze
KA326GRN	Swan Neck Cane Anodised Green
KA332BLK	Swan Neck Cane Anodised Black with X Base
KA19MQ	Quad Foot 19mm Walking Stick

For all model specifications please visit www.kcare.com.au

* Quad foot does not come standard with Swan Neck

Wheelchair Range





Flip back armrest design allows the user closer access to tables, desks & transfer



Compact fold for easy transportation



Swing away footrests

Easily manoeuvrable
with excellent stability

Karma Soma Agile Transit

KMAGLBLK-18T

Features


- Durable aluminium wheelchair
- Tension adjustable backrest and seat
- Flip back armrest design allows the user closer access to tables and desks and transfer
- Detachable swingaway footrests
- Wheelchair frame is compatible with 14 and 24" rear wheel sizes
- AEGIS microbe shield inhibits bacterial and fungal growths, eliminates odours and is easy to clean .
- Crash test approved in accordance with ISO 7176-19
- Wide range of options and accessories available

Specifications

KMAGLBLK-18T

Max User Weight	125kg
Overall Length	1095mm
Overall Width	650mm
Overall Height	915mm
Seat Width	460mm
Seat Depth	460mm
Floor To Seat Height	480mm
Folded Width	280mm
Folded Depth	830mm
Armrest Height	220mm
Backrest Height	510mm
Product Weight	13kg
Brand	Karma

Durable aluminium wheelchair



Detachable swingaway footrests



Tension adjustable backrest and seat

The new standard for transferring wheelchairs

Karma Soma Agile Self Propel

KMAGLBLK-18SP

Features

- Flip back armrest design allows the user closer access to tables and desks and transfer
- Wheelchair frame is compatible with 14 and 24" rear wheel sizes
- AEGIS microbe shield inhibits bacterial and fungal growths, eliminates odours and is easy to clean
- Crash test approved in accordance with ISO 7176-19
- Wide range of options and accessories available
- Available in both Transit & Self Propelled models

Specifications

KMAGLBLK-18SP

Max User Weight	125kg
Overall Length	1120mm
Overall Width	655mm
Overall Height	940mm
Seat Width	460mm
Seat Depth	460mm
Floor To Seat Height	505mm
Folded Width	280mm
Folded Depth	1070mm
Armrest Height	220mm
Backrest Height	510mm
Product Weight	14kg
Brand	Karma

Ellipse-shaped frame design provides style and durability



Multiple self-propel wheel settings



Detachable swingaway footrests



An elliptical-frame design that's durable, stylish, and robust

Karma Flexx Self Propel

KM8022-1818

Features

- Designed with key adjustments to create a suitable sitting environment for the user
- Height adjustable flip back padded armrests
- AEGIS microbe shield inhibits bacterial and fungal growths, eliminates odours and is easy to clean
- Patented ergonomic hand rim provide extra grip and easy propulsion
- Colour coded features for intuitive adjustment
- Crash test approved in accordance with ISO 7176-19
- Available in both Transit & Self Propelled models

Specifications

KM8022-1818

Max User Weight	130kg
Overall Length	1040mm
Overall Width	690mm
Overall Height	1060mm
Seat Width	505mm
Seat Depth	505mm
Floor To Seat Height	510mm
Folded Width	340mm
Folded Depth	1060mm
Min Armrest Height	230mm
Max Armrest Height	280mm
Backrest Height	430mm
Product Weight	14.7kg
Brand	Karma



The 2-in-1 Tilt-in-Space and reclining wheelchair

Karma VIP 2 Tilt Self-Propel

KMVIP2-1618

KMVIP2-1818

Features

- Allows users with partial trunk control to traverse easily and in complete comfort
- Height adjustable and detachable padded armrests
- Detachable swingaway footrests
- Curved multi adjustable padded headrest
- AEGIS microbe shield for hygiene & easy cleaning
- Colour coded features for intuitive adjustment
- Modularized design allows easy upgrades
- Crash test approved in accordance with ISO 7176-19
- Available in 16" / 18"

Specifications

KMVIP2-	1618	1818
Max User Weight	115kg	115kg
Overall Length	1100mm	1130mm
Overall Width	635mm	640mm
Overall Height	1080mm	1100mm
Seat Width	400mm	400mm
Seat Depth	410mm	460mm
Floor to Seat	490mm	490mm
Folded Width	380mm	420mm
Folded Depth	790mm	840mm
Armrest Height	270mm	270mm
Backrest Height	490mm	490mm
Product Weight	17.4kg	17.4kg
Brand	Karma	Karma



S-Ergo patented seating system for added comfort

Ultra lightweight and foldable for easy transport

Karma Ergo Lite Deluxe Transit

KM2512F14-16

KM2512F14-18

Features

- Ultra lightweight, folding, attendant propelled or self-propelled style wheelchair
- Push and pull brake system for 14" rear wheel
- Detachable swingaway footrests
- Optional fixed or quick release 20" self-propel wheel
- Armrests designed to allow closer seating to desks
- Top backrest section & frame folds for ease of transportation
- Double cross braces provide strength & durability

Specifications

KM2512F14-	16	18
Max User Weight	100kg	100kg
Overall Length	935mm	935mm
Overall Width	570mm	615mm
Overall Height	900mm	900mm
Seat Width	410mm	450mm
Seat Depth	420mm	420mm
Floor To Seat	480mm	480mm
Folded Width	280mm	295mm
Folded Depth	635mm	635mm
Armrest Height	215mm	215mm
Backrest Height	420mm	420mm
Product Weight	8.7kg	8.8kg
Brand	Karma	Karma



Easily manoeuvrable
with excellent stability

Karma S-Ergo 125 Self-Propel

- KMS125-16
- KMS125-18
- KMS125-20

Features

- S-Ergo patented seating for added comfort
- Detachable swingaway footrests
- Optional fixed or quick release 20" self-propel wheel
- Flip back armrest design allows the user closer access to tables and desks
- Top backrest section and frame folds for ease of transportation
- Double cross braces for strength & durability
- Crash test approved in accordance with ISO 7176-19
- Wide range of options & accessories available
- AEGIS microbe shield for hygiene & easy cleaning

Specifications

KMS125-	16	18	20
Max User Weight	115kg	115kg	115kg
Overall Length	990mm	960mm	1080mm
Overall Width	570mm	620mm	690mm
Overall Height	910mm	900mm	910mm
Seat Width	400mm	450mm	510mm
Seat Depth	430mm	430mm	430mm
Floor To Seat	500mm	510mm	500mm
Armrest Height	210mm	210mm	200mm
Backrest Height	410mm	400mm	400mm
Product Weight	14.2kg	14.5kg	14.8kg
Brand	Karma	Karma	Karma

Lightweight & compact design



Removable battery



Height adjustable armrests



Compact design allows more room through standard doorways

Karma KP25.2 Power Wheelchair

KMP25.2-16

KMP25.2-18

KMP25.2-20

Features


- Easily manoeuvrable with excellent stability
- Ideal for medium-short travel distances
- Lightweight & compact design
- Features a toolless foldable frame
- Removable battery pack easy to access
- Narrower width for greater turning radius
- Height adjustable armrests
- Flip up, swing away & detachable footrests provide clearance for lateral transfers
- Crash test approved for transportation in accordance with ISO7176-19

Specifications

KMP25.2-	16	18	20
Max User Weight	135kg	135kg	135kg
Overall Width	590mm	640mm	690mm
Min Height	960mm	960mm	960mm
Max Height	1035mm	1035mm	1035mm
Seat Width	410mm	460mm	510mm
Seat Depth	410mm	460mm	460mm
Turning Radius	780mm	955mm	1130mm
Folded Width	380mm	395mm	395mm
Product Weight	67kg	67.6kg	67.6kg
Brand	Karma	Karma	Karma



Lightweight, foldable and durable design



Options and accessories available



Detachable swingaway footrests

Tilt in space and anti-sliding reclining system features anti-shear

Karma VIP 2 Transit Wheelchair 18"

KMVIP2T-18

Features

- Allows users with partial trunk control to traverse easily and in complete comfort
- Lightweight, foldable and durable design
- Height adjustable and detachable padded armrests
- Detachable swingaway footrests
- Curved multi adjustable padded headrest
- AEGIS microbe shield inhibits bacterial and fungal growths, eliminates odours and is easy to clean
- Tension adjustable backrest
- Colour coded features for intuitive adjustment
- Modularized design allows easy upgrades such as Jay and Matrix backrest systems
- Crash test approved in accordance with ISO 7176-19
- Wide range of options and accessories available

Specifications

KMVIP2T	18
Max User Weight	115kg
Overall Length	1130mm
Overall Width	680mm
Overall Height	1100mm
Seat Width	400mm
Seat Depth	460mm
Floor To Seat	490mm
Folded Width	420mm
Folded Depth	840mm
Armrest Height	270mm
Backrest Height	490mm
Product Weight	17.7kg
Brand	Karma

Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.


Beds and Mattresses



Designed for Australian standard mattress sizes



Multiple One Touch Buttons for simple operation



Feature-packed homecare bed with multiple positioning options

NEW PRODUCT

Jannali King Single

KA3200

Features

- Hi - Lo profiling bed
- Patented Autoregression feature
- Wireless Bluetooth remote (Wired option coming)
- Support frame never touches the floor – protection for carers, pets, children
- USB Sockets as standard
- Under bed light with motion detection
- Backup Battery Standard
- Multiple One Touch Buttons for simple operation (Zero Gravity, TV, Flat, CPR)
- Control from smart phone via Bluetooth & App
- Long Single, King Single
- Trendelenburg

Specifications

KA3200

Bed Size	King Single
Safe Working Load	200kg*
Min Height	220 mm
Max Height	745 mm
Length	2020 mm
Width	1060 mm
Trendelenburg	14°
Backrest	0° - 73°
Backrest to Legrest angle	180° - 90°
Minimum frame to floor	80mm
Brand	Alivio

*This is not Max User Weight

Independently tested to conform to AS/NZS3200.2.38 for safety and peace of mind



Light weight removable head and foot boards



Light weight removable head and foot boards

Vienna Bed, Central Lock Castors, Arctic White

MC100C5AW

Features


- Independently tested to conform to AS/NZS3200.2.38 for safety and peace of mind
- Light weight removable head and foot boards
- Four section profiling bed for added comfort and extra adjustability for the user
- Integrated 200mm bed extension
- Rattle free plastic sockets
- Clip on/off breathable mattress panels
- Knee break with plastic lift handles converts to calf platform

Specifications

MC100C5AW

Max User Weight	225kg
Overall Length	2180mm
Overall Width	960mm
Product Weight	125kg
Brand	My Comfort

Low interface pressure for improved blood circulation and wound healing



Auto-set with 5 manual firmness settings



Clinically proven foam is breathable, medical grade and temperature resistant

Dual Flex Hybrid System Mattress Replacement & Control Unit

HYB-DF-S

HYB-DF-CU

Features

- Clinically proven foam is breathable, medical grade and temperature resistant
- Zonal area design to provide appropriate care to sacrum, head, shoulder blades & heel
- Low interface pressure for improved blood circulation and wound healing
- 2 way stretch washable coverlet to prevent shearing
- Dynamic mode 2:1 alternating therapy, 10 minute cycles
- Ultra quiet pump emitting less than 28db
- Available in single and king single sizes

Specifications

HYB-DF-M

Overall Length	2000mm
Overall Width	900mm
Thickness	180mm
Product Weight	8kg
Max User Weight	250kg
Brand	Carilex

HYB-DF-CU

Pressure Sore Grading	4
Waterlow Scale	20+
Product Weight	2.5kg
Max User Weight	250kg
Brand	Carilex

Extra soft 100% TPU material



Dual Plus Pump



CPR dial deflates the entire mattress in seconds



Alternating 8" mattress replacement system designed to suit individual user needs

Dual Plus Mattress and Pump

APR-DU-PL

Features

- Alternating 8" mattress replacement system designed to suit individual user needs
- Extra soft 100% TPU material
- Durable and washable coverlet
- Available in single and king single, Euro Single and Euro King Single
- Dynamic & Constant Low Pressure therapy modes to meet the individual's needs
- Green bridging cells provides constant pressure relief and localised healing
- Square air cell design with no gap between cells for maximum support & pressure redistribution

Specifications

APR-DU-PL

Overall Length	2000mm
Overall Width	900mm
Thickness	210mm
Product Weight	11.3kg
Max User Weight	250kg
Brand	Carilex


No fuss cable and cord management



Centrius Pump



CPR pull tab deflates the entire mattress in seconds



Alternating 7" mattress replacement system designed to suit individual user needs

Centrius Pump, Entrix 7" Mattress

CEN-7-XL

Features

- Alternating 7" alternating mattress replacement system designed to suit individual user needs
- Low profile alternating mattress envelopes and off loads patient body mass
- Easy connection between power unit and mattress
- No fuss cable and cord management
- 2 way stretch washable cover
- CPR pull tab deflates the entire mattress in seconds
- Extra 2" zip-on foam protection for user support
- Available in single, king single & Euro king single

Specifications

CEN-7-XL

Overall Length	2000mm
Overall Width	1070mm
Thickness	180mm
Product Weight	4.3kg
Max User Weight	150kg
Brand	Carilex

Reduced pressure on heels & feet



Help alleviate pressure on the heels



Pressure Care Plus Sloped Heel Mattress

- MC11MPPS
- MC11MPPKS
- MC11MPPEs
- MC11MPPEKS

Features

- Australian designed and made
- Static pressure relieving foam mattress
- Integrated sloped heel area helps offload pressure
- Carrflex Extreme & Defence MT520 cover
- Antimicrobial, fire retardant and hypoallergenic
- Welded 360° zipper with waterfall flap
- Latex free and anti static
- Good average pressure below 32mmHg in supine and lateral positions

Specifications

Overall Length	2000mm
Overall Height	150mm
Max User Weight	250kg
Brand	Carilex

MC11MPPS

Overall Width	900mm
Product Weight	9kg

MC11MPPKS

Overall Width	1050mm
Product Weight	10.7kg

MC11MPPEs

Overall Width	850mm
Product Weight	8.5kg

MC11MPPEKS

Overall Width	1000mm
Product Weight	10kg

Specially designed deep air cells for excellent immersion and pressure relief



Carichair Pump



Extra Support at the Sacrum Area



Alternating air overlay for recliners and chairs with battery backup

Carichair

APC-CARI-O

Features

- Extra Support at the Sacrum Area
- Seamless Transition from Bed to Chair
- Specially designed deep air cells for excellent immersion and pressure relief
- Continuous, uninterrupted care whenever, wherever
- Battery Backup up to 8 hours
- Dynamic and Constant Low Pressure Relief to meet the individual needs
- Available in cushions sizes

Specifications

APC-CARI-O

Safe Working Load	Up to 150kg
Available Size / Dimensions	180cm x 50cm x 14cm
Travel / Weight	2.2kg
Compatible Pumps	CariChair
Brand	Carilex

Bariatric APM, SWL 450kg



Battery back up, 2 hours



Envelop Mattress Replacement System

NEW
PRODUCT

Envelop Mattress & Pump

ENV-BAR-M

Features

- Bariatric & Standard Modes up to SWL 450kg
- S.M.A.R.T. Technology automatically determines patient load on the mattress
- S31 tested, with 79% immersion, one of the best immersion effect in the industry
- Dynamic and Constant Low Pressure therapy modes to meet the individual's needs
- Green bridging cells provides constant pressure relief and localised healing
- Square air cell design with no gap between cells for maximum support and pressure redistribution
- Bariatric extension, bolsters extend the mattress from 90cm to 105cm and 120cm

Specifications

ENV-BAR-M

Overall Length	2000mm
Overall Width	900mm
Thickness	280mm
Product Weight	11.2kg
Pump Weight	6kg
Pump Dimensions	375 x 320 x 185mm (W,L,H)
Max User Weight	150kg
Brand	Carilex

Notes

[illegible]

Improving your life everyday™

At Your Service

Sales, administration
and general enquiries

National 1300 783 783

Your Customer Service Team
customerservice@kcare.com.au

kcare.com.au

Last Published
date: **Apr 2023**

All information is considered to be true and correct at the date of publication. Changes in circumstances after the time of publication may impact on the accuracy of the information.

